

A P P L I C A T I O N

for

UNITED STATES LETTERS PATENT

on

METHODS OF DIAGNOSING LIVER FIBROSIS

by

Steven L. Rose

Esther H. Oh

and

Michael J. Walsh

CERTIFICATE OF MAILING BY "EXPRESS MAIL"

"EXPRESS MAIL" MAILING LABEL NUMBER: EL 857042971 US

DATE OF DEPOSIT: February 28, 2002

I HEREBY CERTIFY THAT THIS PAPER OR FEE IS BEING
DEPOSITED WITH THE UNITED STATES POSTAL SERVICE
"EXPRESS MAIL POST OFFICE TO ADDRESSEE" SERVICE UNDER
37 C.F.R. 1.10 ON THE DATE INDICATED ABOVE, AND IS
ADDRESSED TO: BOX PATENT APPLICATION,
COMMISSIONER FOR PATENTS WASHINGTON, D.C. 20231.

Karly Tambara

(TYPED OR PRINTED NAME OR PERSON MAILING PAPER OR FEE)

K. Tambara

(SIGNATURE OF PERSON MAILING PAPER OR FEE)

Attorneys

CAMPBELL & FLORES LLP
4370 La Jolla Village Drive, 7th Floor
San Diego, California 92122
USPTO CUSTOMER NO. 23601

Sheets of Drawings: 2

Docket No.: P-PM 4978

20087183 022802

BACKGROUND OF THE INVENTION

The present invention relates generally to the
5 fields of hepatology and fibrosis and, more specifically,
to a panel of serological markers which together are
diagnostic of liver fibrosis.

Progressive fibrosis of the liver, kidney,
lungs and other organs frequently results in organ
failure that leads to organ transplantation or death,
affecting millions in the United States and worldwide.
Hepatic fibrosis, for example, is the leading
non-malignant gastrointestinal cause of death in the
United States, and the progression of fibrosis is the
single most important determinant of morbidity and
mortality in patients with chronic liver disease.
Furthermore, the process of fibrosis is common to liver
diseases of many etiologies, including chronic viral
hepatitis B and C, autoimmune liver disease such as
autoimmune hepatitis, alcoholic liver disease, fatty
liver disease; primary biliary cirrhosis; and
drug-induced liver disease. The fibrosis seen in these
disorders results from chronic insults to the liver such
as viral infection, alcohol or drugs.

Hepatitis C, for example, is one of the leading causes of chronic liver disease in the United States,

where an estimated 3.9 million people are chronically infected with hepatitis C virus (HCV) and approximately 30,000 new cases of acute HCV occur each year (Alter, Semin. Liver Dis. 15:5-14 (1995)). The prevalence of
5 hepatitis C is estimated to be 1.8% in the United States, with up to 10,000 deaths per year likely resulting from chronic hepatitis C infection (Alter, *supra*, 1995).

While hepatic fibrosis is a reversible process resulting in the accumulation of extracellular matrix,
10 liver cirrhosis is an irreversible process characterized by thick bands of matrix which completely encircle the parenchyma to form nodules. Untreated, fibrosis of the liver leads to cirrhosis and eventually end-stage liver disease or cancer. Cirrhosis of the liver is a common
15 condition that frequently goes undetected. For example, in a large sample of the general Danish population, the prevalence of liver cirrhosis was 4.5%, of which one-third were undiagnosed at the time of death (Graudal, J. Intern. Med. 230:165-171 (1991)).

20 Timely and accurate diagnosis of liver fibrosis is important to effective medical treatment. As an example, patients with hepatitis C and cirrhosis are less likely to respond to treatment with α -interferon compared to patients with less advanced disease (Davis, Hepatology
25 26(Supp. 1):122-127S). Similarly, treatments for chronic HCV infection can be contra-indicated in patients with histologically advanced and decompensated disease (NIH Consensus Development Conference Panel Statement, Hepatology 26 (Suppl. 1):25-105S (1997)). The importance
30 of early diagnosis is further emphasized by the serious early complications such as variceal rupture that are

20087183, 0000002

associated with cirrhosis; these complications can be prevented by early detection of cirrhosis (Calés and Pasqual, Gastroenterol. Clin. Biol. 12:245-254 (1988)).

Diagnosis of the presence or severity of
 5 fibrotic liver disease is difficult, with liver biopsy currently the most reliable method available. Unfortunately, liver biopsy has several limitations: pain in about 30% of patients; the risk of severe complications such as hemorrhage or infection; a death
 10 rate of 3 in 10,000; and the cost of hospitalization (Nord, Gastrointest. Endosc. 28:102-104 (1982); Cadranet et al., Hepatology 32:47-481 (2000); and Poynard et al., Can. J. Gastroenterol. 14:543-548 (2000)). Furthermore, slowly progressive diseases such as hepatitis C require
 15 repeated biopsies for continual assessment of disease progression, thus compounding the risks and costs of the procedure. Finally, biopsy can fail to detect disease because of the heterogeneous distribution of pathological changes in the liver; it is not surprising, then, that
 20 false negatives are seen in a significant percentage of cases biopsied (Nord, *supra*, 1982).

For years there has been a search for biochemical or serological markers which reflect fibrotic processes in liver disease and which can serve as a
 25 surrogate for liver biopsy. However, the performance of any single marker has not been good enough to substitute for the biopsy procedure in detecting or staging fibrosis. Thus, there is a need for a non-invasive method of diagnosing the presence or severity of liver
 30 fibrosis. The present invention satisfies this need by providing a convenient and reliable method for detection

20087188.023002

of liver fibrosis that is suitable for serial testing. Related advantages are provided as well.

SUMMARY OF THE INVENTION

20087188-028302

The present invention provides a method of
5 diagnosing the presence or severity of liver fibrosis in
an individual by detecting α 2-macroglobulin (α 2-MG) in a
sample from the individual; detecting hyaluronic acid
(HA) in a sample from the individual; detecting tissue
inhibitor of metalloproteinases-1 (TIMP-1) in a sample
10 from the individual; and diagnosing the presence or
severity of liver fibrosis in the individual based on the
presence or level of α 2-MG, HA and TIMP-1. A method of
the invention can be useful, for example, for
differentiating no or mild (F0-F1) liver fibrosis from
15 moderate to severe (F2-F4) liver fibrosis.

The methods of the invention for diagnosing the
presence or severity of liver fibrosis can be useful in a
variety of patient populations including, but not limited
to, those with viral hepatitis, autoimmune liver disease
20 such as autoimmune hepatitis, alcoholic liver disease,
fatty liver disease and drug-induced liver disease. In
one embodiment, a method of the invention is used to
diagnose the presence or severity of liver fibrosis in an
individual infected with hepatitis C virus.

25 A variety of means can be useful for detecting
 α 2-MG, HA and TIMP-1 in a sample. In one embodiment, the
invention is practiced by determining the level of α 2-MG
protein in a sample from the individual to be diagnosed
using, for example, one or more α 2-MG-specific binding

agents such as anti- α 2-MG antibodies. In another embodiment, a method of the invention is practiced by determining the level of α 2-MG activity in a sample from the individual.

5 A variety of means also can be used in a method of the invention to detect hyaluronic acid in a sample. In one embodiment, the invention is practiced by determining the level of HA in a sample, for example, using one or more HA-specific binding agents such as
10 HA-binding proteins (HABPs) or anti-HA antibodies.

Similarly, a variety of means can be used in a method of the invention to detect TIMP-1 in a sample. In one embodiment, the invention is practiced by determining the level of TIMP-1 protein in a sample from the
15 individual to be diagnosed. The level of TIMP-1 protein can be determined, for example, using one or more TIMP-1-specific binding agents such as anti-TIMP-1 antibodies. In another embodiment, the invention is practiced by assaying for TIMP-1 activity in a sample
20 from the individual to be diagnosed.

The invention provides, for example, a method of diagnosing the presence or severity of liver fibrosis in an individual by determining the level of α 2-MG protein in a sample from the individual; determining the
25 level of HA in a sample from the individual; and determining the level of TIMP-1 protein in a sample from the individual; and diagnosing the presence or severity of liver fibrosis in the individual based on the levels of α 2-MG protein, HA and TIMP-1 protein. If desired, the

20087483.0238002

level of α 2-MG protein, HA and TIMP-1 protein each can be determined using an enzyme-linked assay.

A variety of samples can be useful in practicing the methods of the invention including, for example, blood, serum, plasma, urine, saliva and liver tissue. In one embodiment, a single sample is obtained from the individual to be diagnosed. Such a sample can be, for example, a serum sample. Such a sample also can be, for example, a tissue sample, for example, a liver biopsy sample.

The present invention further provides a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual. The method includes the steps of contacting an appropriate dilution of a sample from the individual with anti- α 2-MG antibody under conditions suitable to form a first complex of α 2-MG and anti- α 2-MG antibody; washing the first complex to remove unbound molecules; determining the amount of α 2-MG-containing first complex; contacting an appropriate dilution of a sample from the individual with a HA-binding protein under conditions suitable to form a second complex of HA and HA-binding protein; washing the second complex to remove unbound molecules; determining the amount of HA-containing second complex; contacting an appropriate dilution of a sample from the individual with anti-TIMP-1 antibody under conditions suitable to form a third complex of TIMP-1 and anti-TIMP-1 antibody; washing the third complex to remove unbound molecules; determining the amount of TIMP-1-containing third complex; and differentiating no or mild liver fibrosis from moderate to severe liver

20087189, 00000000

fibrosis in the individual based on the amounts of α 2-MG, HA and TIMP-1-containing complexes.

The methods of the invention can be practiced by detecting the three markers α 2-MG, HA and TIMP-1, 5 without detecting additional serological markers, or can be combined with a detection method for one or more additional markers. Thus, in one embodiment, the invention is practiced by detecting α 2-MG, HA and TIMP-1 and also detecting at least one of the following markers 10 of fibrosis: N-terminal procollagen III propeptide (PIIINP), laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- β 1, IL-10, apoA1, apoA2 or apoB. In a further embodiment, the presence or severity of liver 15 fibrosis is diagnosed by detecting α 2-MG, HA, TIMP-1 and YKL-40 in a sample from an individual.

The present invention also provides a method of monitoring the efficacy of anti-fibrotic therapy in a patient by detecting α 2-macroglobulin in a sample from a 20 patient administered an anti-fibrotic therapy; detecting hyaluronic acid (HA) in a sample from the patient; detecting tissue inhibitor of metalloproteinases-1 (TIMP-1) in a sample from the patient; and determining the presence or severity of liver fibrosis in the patient 25 based on the presence or level of α 2-MG, HA and TIMP-1, thereby monitoring the efficacy of anti-fibrotic therapy. Such a method can further include, if desired, comparing the presence or severity of liver fibrosis determined in step (d) to the presence or severity of liver fibrosis in 30 the patient at an earlier time. The methods of the invention can be used to monitor, for example, the

10087188.022802

progression or regression of fibrosis over time in a patient treated with one or more anti-fibrotic therapies, or to compare, for example, the efficacies of two or more anti-fibrotic therapies.

5 In one embodiment, at most three markers of fibrosis are detected. In another embodiment, the method includes the step of detecting in a sample from the patient at least one marker selected from the group consisting of: PIIINP, laminin, tenascin, collagen
10 type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- β 1, IL-10, apoA1, apoA2, and apoB.

A variety of means can be useful for detecting α 2-MG, HA and TIMP-1 in a method of the invention. Step
15 (a) can be practiced, for example, by determining the level of α 2-MG protein in the sample. In one embodiment, the level of α 2-MG protein is determined using one or more anti- α 2-MG antibodies. Step (b) can be practiced, for example, by determining the level of HA in the
20 sample. In one embodiment, the level of HA is determined using one or more HA-binding proteins. Step (c) can be practiced, for example, by determining the level of TIMP-1 protein in said sample. In one embodiment, the level of TIMP-1 protein is determined using one or more
25 anti-TIMP-1 antibodies.

Further provided herein is a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual by determining an α 2-MG level in a sample from the individual;
30 determining a HA level in a sample from the individual;

10087188 00000000

determining a TIMP-1 level in a sample from the individual; and diagnosing the individual as having no or mild liver fibrosis when the α 2-MG level is below an α 2-MG cut-off value X1, the HA level is below a HA cut-off value Y1 or the TIMP-1 level is below a TIMP-1 cut-off value Z1; diagnosing the individual as having moderate to severe liver fibrosis when the α 2-MG level is above an α 2-MG cut-off value X2, the HA level is above a HA cut-off value Y2 and the TIMP-1 level is above a TIMP-1 cut-off value Z2; and diagnosing remaining individuals as having an indeterminate status.

The methods of the invention based on dual cut-off values for the levels of the α 2-MG, HA and TIMP-1 markers can be useful in differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in a variety of patient populations. The methods of the invention can be useful, for example, in diagnosing an individual having a liver disease such as viral hepatitis, autoimmune liver disease such as autoimmune hepatitis, alcoholic liver disease, fatty liver disease or drug-induced liver disease. In one embodiment, the methods of the invention are used to differentiate no or mild liver fibrosis from moderate to severe liver fibrosis in an individual infected with hepatitis C virus. Samples useful in the methods of the invention include, but are not limited to, blood, serum, plasma, urine, saliva and liver tissue. In one embodiment, a method of the invention is practiced by determining the α 2-MG level, HA level and TIMP-1 level in one or more serum samples from the individual to be diagnosed.

Thus, the present invention provides, for example, a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual in which the differentiation is based on an X1 cut-off value between 1.8 and 2.2 mg/ml; a Y1 cut-off value between 31 and 39 ng/ml; a Z1 cut-off value between 900 and 1100 ng/ml; an X2 cut-off value between 1.8 and 2.2 mg/ml; a Y2 cut-off value between 54 and 66 ng/ml; and a Z2 cut-off value between 1415 and 1735 ng/ml. In a particular embodiment, the differentiation is based on an X1 cut-off value of 2.0 mg/ml; a Y1 cut-off value of 35 ng/ml; a Z1 cut-off value of 1000 ng/ml; an X2 cut-off value of 2.0 mg/ml; a Y2 cut-off value of 60 ng/ml; and a Z2 cut-off value of 1575 ng/ml. In another embodiment, the differentiation is based on an X1 cut-off value of 2.0 mg/ml; a Y1 cut-off value of 37 ng/ml; a Z1 cut-off value of 1100 ng/ml; an X2 cut-off value of 2.0 mg/ml; a Y2 cut-off value of 60 ng/ml; and a Z2 cut-off value of 1575 ng/ml. In a further embodiment, X1, Y1, Z1, X2, Y2 and Z2 are selected such that, in a population having up to 30% liver fibrosis prevalence, at least 65% of individuals in the population are diagnosed as having no/mild fibrosis or moderate/severe fibrosis with an accuracy of at least 90%. In another embodiment, X1, Y1, Z1, X2, Y2 and Z2 are selected such that, in a population having up to 30% liver fibrosis prevalence, at least 65% of individuals in said population are diagnosed as having no/mild fibrosis or moderate/severe fibrosis with a positive predictive value of at least 90% and a negative predictive value of at least 90%. In yet a further embodiment, X1, Y1, Z1, X2, Y2 and Z2 are selected such that, in a population having up to 10% fibrosis prevalence, at least 70% of individuals in the population

are diagnosed as having no/mild fibrosis or moderate/severe fibrosis with an accuracy of at least 90%.

The present invention also provides a method of
5 diagnosing the presence or severity of liver fibrosis in
an individual by comparing a level of a first fibrotic
marker X in the individual to a cut-off value X1 to
determine whether the individual is positive for the
first fibrotic marker X; comparing a level of a second
10 fibrotic marker Y in the individual to a cut-off value Y1
to determine whether the individual is positive for the
second fibrotic marker Y; and diagnosing the presence or
severity of liver fibrosis in the individual based on
positivity or negativity for X and Y, where, in a
15 population with up to 40% fibrosis prevalence, at least
65% of individuals in the population are diagnosed with
an accuracy of at least 90%.

A method of the invention can include, if
desired, comparing a level of a third fibrotic marker Z
20 in the individual to a cut-off value Z1 to determine
whether the individual is positive for the third fibrotic
marker Z and diagnosing the presence or severity of liver
fibrosis in the individual based on positivity or
negativity for X, Y and Z. In one embodiment, the first
25 fibrotic marker is α 2-MG, the second fibrotic marker is
HA, and the third fibrotic marker is TIMP-1.

In another embodiment, the levels of at least
three fibrotic markers are compared, and, in a further
embodiment, the levels of exactly three fibrotic markers
30 are compared. In additional embodiments, the levels of

at least four or at least five fibrotic markers are compared. A method of the invention can be useful, for example, to differentiate no or mild liver fibrosis from moderate to severe liver fibrosis.

5 In a specific embodiment, a method of the invention serves to diagnose at least 65% of individuals in a population with up to 30% fibrosis prevalence with an accuracy of at least 93%. In a further embodiment, a method of the invention serves to diagnose at least 70%
10 of individuals in a population with up to 20% fibrosis prevalence with an accuracy of at least 94%. In yet a further embodiment, a method of the invention serves to diagnose at least 70% of individuals in a population with up to 10% fibrosis prevalence with an accuracy of at
15 least 96%.

The present invention further provides a method of diagnosing the presence or severity of liver fibrosis in an individual by comparing a level of a first fibrotic marker X in the individual to a cut-off value X1 to
20 determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second fibrotic marker Y in the individual to a cut-off value Y1 to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or
25 severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where the cut-off values X1 and Y1 are optimized individually to give a desired performance characteristic.

If desired, a method of the invention can
30 include the steps of comparing a level of a third

10087188 022802

fibrotic marker Z in the individual to a cut-off value Z1 to determine whether the individual is positive for the third fibrotic marker Z and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z, where the cut-off values X1, Y1 and Z1 are optimized individually to give a desired performance characteristic. In one embodiment, levels of α 2-MG, HA and TIMP-1 are compared. In another embodiment, the cut-off values are optimized using design of experiments (DOE) analysis. In further embodiments, the levels of exactly three, at least three, at least four, or at least five fibrotic markers are compared. A method of the invention can be useful, for example, in differentiating no or mild liver fibrosis from moderate to severe liver fibrosis.

Further provided by the invention is a method of diagnosing the presence or severity of liver fibrosis in an individual by comparing a level of a first fibrotic marker X in the individual to two cut-off values X1 and X2 to determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second fibrotic marker Y in the individual to two cut-off values Y1 and Y2 to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where the cut-off values X1, Y1, X2 and Y2 are optimized individually to give a desired performance characteristic. A method of the invention can further include the steps of comparing a level of a third fibrotic marker Z in the individual to two cut-off values Z1 and Z1 to determine whether the individual is positive

for the third fibrotic marker Z; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z, where the cut-off values X1, Y1, Z1, X2, Y2 and Z2 are
 5 optimized individually to give a desired performance characteristic. Cut-off values can be conveniently optimized, for example, using DOE analysis.

BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1 shows the nucleic acid sequence (SEQ
 10 ID NO: 1) and corresponding amino acid sequence (SEQ ID NO: 2) for mature human α 2-macroglobulin available from Genbank accession M36501.

Figure 2 shows the nucleic acid sequence (SEQ
 ID NO: 3) and corresponding amino acid sequence (SEQ ID
 15 NO: 4) for human tissue inhibitor of metalloproteinases-1 (TIMP-1) available from Genbank accession NM_003254.

DETAILED DESCRIPTION OF THE INVENTION

As disclosed herein, the serum levels of a number of biochemical markers were analyzed in a patient
 20 population with confirmed hepatitis C and having a known Metavir stage (fibrosis score) of F0 to F4, where F0 represents very low or no fibrosis; F1, F2 and F3 represent intermediate fibrosis stages; and F4 represents severe fibrosis (Knodell et al., Hepatology 1:431-435
 25 (1981)). See Tables 2 and 3. Using Design of Experiments (DOE) analysis for simultaneous variation of multiple cut-off values, a four-marker panel made up of hyaluronic acid (HA), PIIINP, collagen type IV and

α2-macroglobulin (α2-MG) was identified which was capable of differentiating F0-F1 (no or mild) fibrosis from F2-F4 (moderate to severe) fibrosis with an accuracy of about 77% in a patient population with a fibrosis prevalence of 60%.

As further disclosed herein in Example I, two three-marker panels, α2-MG/HA/TIMP-1 and α2-MG/HA/YKL-40, also performed well in differentiating F0-F1 fibrosis from F2-F4 fibrosis when cut-offs were optimized using DOE analysis. In particular, the α2-MG/HA/TIMP-1 and α2-MG/HA/YKL-40 panels each performed better than the four-marker panel and were capable of differentiating F0-F1 from F2-F4 fibrosis with about 80% accuracy in the study population. As can be seen in Table 6, line 15, for example, the α2-MG/HA/TIMP-1 panel performed with a sensitivity of 83.48% and a specificity of 75.95% in the study population having 60% fibrosis prevalence. These results demonstrate that the α2-MG/HA/TIMP-1 three-marker panel can be useful for differentiating no or mild fibrosis from moderate to severe fibrosis.

Based on these findings, the present invention provides a method of diagnosing the presence or severity of liver fibrosis in an individual by detecting α2-MG in a sample from an individual; detecting HA in a sample from the individual; detecting TIMP-1 in a sample from the individual; and diagnosing the presence or severity of liver fibrosis in the individual based on the presence or levels of α2-MG, HA and TIMP-1. A method of the invention can be useful, for example, for differentiating no or mild (F0-F1) liver fibrosis from moderate to severe (F2-F4) liver fibrosis.

Liver and other fibrotic disorders

The methods of the invention can be useful for diagnosing the presence or severity of liver fibrosis in a variety of individuals including those at risk for, or
5 having one or more symptoms of, a liver disorder characterized by fibrosis. The methods of the invention can be used to diagnose liver fibrosis in an individual having, for example, viral hepatitis such as hepatitis A, B or C virus or a human immunodeficiency
10 virus (HIV) such as HIV-1; chronic persistent hepatitis or chronic active hepatitis; autoimmune liver disease such as autoimmune hepatitis; alcoholic liver disease; fatty liver disease; primary biliary cirrhosis; primary sclerosing cholangitis, biliary atresia; liver disease
15 resulting from medical treatment (drug-induced liver disease); or a congenital liver disease. The methods of the invention can be extremely useful, for example, in alleviating concerns of potential liver damage due to methotrexate treatment. Periodic monitoring of liver
20 fibrosis in individuals treated with methotrexate or other drugs associated with risk of liver damage can be conveniently performed using the non-invasive methods of the invention, without the risks associated with liver biopsy.

25 In one embodiment, the methods of the invention are useful for differentiating individuals having a Metavir score of F0 or F1 from individuals having a Metavir score of F2, F3 or F4. Metavir scoring is a well accepted system for grading liver biopsy specimens and is
30 described in Knodell, *supra*, 1981. F0 is equivalent to the absence of fibrosis; F1 signifies portal fibrosis

10087188, 0333002

without septa. F2 signifies portal fibrosis with a few septa. F3 signifies numerous septa without cirrhosis. F4 signifies cirrhosis.

It is understood that the methods of the invention are useful for diagnosing the presence or severity of fibrosis associated with a variety of fibrotic disorders, including but not limited to liver fibrosis, pulmonary fibrosis, kidney fibrosis, prostate fibrosis and breast fibrosis. The methods of the invention can be applied, without limitation, to diagnosing the presence or severity of pulmonary fibrosis, for example, idiopathic pulmonary fibrosis or emphysema; kidney fibrosis; bladder fibrosis; periureteric fibrosis or retroperitoneal fibrosis; endomyocardial fibrosis, aortic aneurysm disease; rheumatoid diseases such as rheumatoid arthritis or systemic lupus erythematosus; or another fibrotic disorder such as Alzheimer's disease. It is understood that a α 2-MG/HA/TIMP-1, α 2-MG/HA/YKL-40 or α 2-MG/HA/TIMP-1/YKL-40 panel or other combination of markers disclosed herein as useful for diagnosing the presence or severity of liver fibrosis also can be used to diagnose the presence or severity of fibrosis in another disorder.

It is understood that the diagnostic methods of the invention are applicable to a variety of individuals including individuals with chronic or active disease, individuals with one or more symptoms of fibrotic disease, asymptomatic or healthy individuals and individuals at risk for one or more fibrotic diseases. It further is clear to the skilled person that the

methods of the invention can be useful, for example, to corroborate an initial diagnosis of disease or to gauge the progression of fibrosis in an individual with a previous definitive diagnosis of fibrotic disease. The
5 methods of the invention can be used to monitor the status of fibrotic disease over a period of time and further can be used, if desired, to monitor the efficacy of therapeutic treatment. If desired, the results obtained from a sample from an individual undergoing
10 therapy can be compared, for example, to the individual's baseline results prior to treatment, to results earlier during treatment, or to a historic or reference value.

Samples

A variety of samples can be useful in
15 practicing the methods of the invention including, for example, blood, serum, plasma, urine, saliva and liver tissue. In one embodiment, a single sample is obtained from the individual to be diagnosed. Such a sample can be, for example, a serum sample.

20 As used herein, the term "sample" means a biological specimen that contains one or more fibrotic markers such as α 2-MG, HA or TIMP-1. A sample can be, for example, a fluid sample such as whole blood, plasma, saliva, urine, synovial fluid or other bodily fluid, or a
25 tissue sample such as a lung, liver, kidney, prostate or breast tissue sample. One skilled in the art understands that fluid samples can be diluted, if desired, prior to analysis.

One skilled in the art understands that a single sample can be obtained from the individual to be diagnosed and can be subdivided prior to detecting $\alpha 2$ -MG-, HA- and TIMP-1. One skilled in the art also
5 understands that, if desired, two or more samples can be obtained from the individual to be diagnosed and that the samples can be of the same or a different type. In one embodiment, $\alpha 2$ -MG-, HA- and TIMP-1 each are detected in serum samples. In another embodiment, a single serum
10 sample is obtained from an individual and subdivided prior to detecting $\alpha 2$ -MG-, HA- and TIMP-1.

$\alpha 2$ -macroglobulin

The methods of the invention rely, in part, on detecting $\alpha 2$ -macroglobulin in a sample. $\alpha 2$ -MG is a
15 conserved, highly abundant component of plasma that functions as a broad spectrum protease-binding protein to clear active proteases from tissue fluids. Unlike active site protease inhibitors, members of the $\alpha 2$ -macroglobulin family do not inactivate the catalytic activity of their
20 protease substrates but act by physical entrapment of the target protease within the folds of the $\alpha 2$ -MG family member. $\alpha 2$ -MG is itself cleaved by target proteases; reorganization of the $\alpha 2$ -MG molecule results in sequestering of the target protease within an internal
25 pocket of the $\alpha 2$ -MG molecule (Starkey et al., Biochem. J. 131:823-831 (1973)). While an $\alpha 2$ -MG entrapped protease is sterically prevented from interacting with macromolecular substrates such as proteins, it remains active against low molecular mass
30 substrates, such as amide and ester compounds, able to diffuse into the $\alpha 2$ -MG cage to access the enzymatic site.

10057188.022802

Thus, α 2-MG activity is characterized, in part, by the ability to inhibit proteolytic activity but not amidolytic activity of a protease substrate. α 2-MG also is characterized by the ability to shield entrapped
5 proteases from antibodies and high molecular mass active site inhibitors. For example, trypsin bound by α 2-MG is protected from inhibition by soybean trypsin inhibitor (STI).

In contrast to the restricted specificity of
10 active-site protease inhibitors, α 2-MG acts on a broad spectrum of proteases with diverse substrate specificity and catalytic activity. Such target proteases include trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin and papain. Substrate diversity is
15 determined, in part, by the α 2-MG "bait" region, a highly flexible and solvent-exposed sequence of 30-40 residues that contains at least one site sensitive to cleavage by each of the major classes of proteolytic enzyme.

As used herein, the term " α 2-macroglobulin" is
20 synonymous with " α 2-MG" and means a protein with significant structural homology to human α 2-MG (SEQ ID NO: 2) and having broad spectrum protease inhibitory activity. α 2-MG contains a unique thiol ester bond that is inactivated by small primary amines such as
25 methylamine. Thus, α 2-MG activity can be characterized, in part, by methylamine-sensitive protease inhibitory activity. α 2-MG can be distinguished, if desired, from other members of the α 2-macroglobulin family such as related protease-binding proteins and C3, C4 and C5 of
30 the complement system (Sottrup-Jensen, " α 2-Macroglobulin and Related Thiol Ester Plasma Proteins," in Putnam

10087188.023802

(Ed.), The Plasma Proteins: Structure, Function and Genetic Control Second edition, Orlando: Academic Press (1987), pages 191-291. It is understood that an assay for detecting α 2-MG can be specific for α 2-MG or can
5 additionally detect one or more other members of the α 2-macroglobulin family.

The methods of the invention rely, in part, on detecting α 2-macroglobulin in a sample. As used herein, the phrase "detecting α 2-MG" means any quantitative or
10 qualitative assay for determining the presence of α 2-MG. As used herein, the phrase "determining the level of α 2-MG" means any direct or indirect quantitative assay for α 2-MG.

Similarly, detecting any specified fibrotic
15 marker in a sample means determining whether the marker is present in the sample, said fibrotic marker having a positive or negative correlation with liver fibrosis or with another fibrotic disorder such as are described herein above. It is understood that detection can refer
20 to non-quantitative analysis, for example, the presence or absence of a particular trait, variable or biochemical or serological substance.

Diagnosis is based on analyzing the sample for the presence or level of the fibrotic marker or other
25 characteristic and comparing it to a reference value, where the reference value serves to assist in differentiating those with a fibrotic disorder from other individuals. Where the fibrotic marker is a biochemical or serological marker, determining a "level" in a sample
30 means quantifying the fibrotic marker by determining, for

10087488, 023802

5

10

25

- the sample α 2-MG concentration. α 2-MG also can be detected, for example, by laser immunonephelometry using a Behring Nephelometer Analyzer (Fink et al., J. Clin. Chem. Clin. Biol. Chem. 27:261-276 (1989)) and rabbit
- 5 anti-human α 2-MG antiserum as described in Naveau et al., Dig. Diseases Sci. 39:2426-2432 (1994), or using the nephelometry assay commercially available from Beckman Coulter (Brea, CA; kit #449430). Furthermore, monoclonal and polyclonal anti- α 2-MG antibodies useful in
 - 10 immunoassays can be readily obtained from a variety of sources. As examples, affinity purified goat anti-human α 2-MG and peroxidase-labeled goat anti-human α 2-MG antibodies suitable for immunoassays such as ELISA assays and western blotting are available from Cedarlane
 - 15 Laboratories Limited (Ontario, Canada; CL20010AP and CL20010APHP) and Affinity Biologicals Incorporated (Ontario, Canada; GAA2M-AP and GAA2M-APHRP). Levels of α 2-MG protein also can be determined by quantifying the amount of purified α 2-MG protein.
 - 20 Purification of α 2-macroglobulin can be achieved, for example, by HPLC, alone or in combination with mass spectrophotometry, or as described, for example, in Hall and Roberts, Biochem. J. 171:27-38 (1978) or Imber and Pezzo, J. Biol. Chem. 256:8134-8139 (1981)).
 - 25 Quantitation can be determined by well known methods including Bradford assays, Coomassie blue staining and assays for radiolabeled protein.

A variety of assays for α 2-MG activity also can be useful for detecting α 2-MG or determining a level of

- 30 α 2-MG in a sample according to a method of the invention. α 2-MG can be detected or a level of α 2-MG can be determined indirectly, for example, as a function of

inhibition of target protease activity, without a corresponding inhibition of amidolytic activity. As discussed above, α 2-MG-bound proteases retain the ability to hydrolyze amide and ester bonds of small substrates, even while high molecular mass substrates such as proteins cannot be hydrolyzed (see, for example, Armstrong et al., Develop. Compar. Immunol. 23:375-390 (1999)). As an example, α 2-MG can be detected or the level of α 2-MG can be determined by assaying for inhibition of trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin, or papain activity without inhibition of amidolytic activity. Convenient substrates to be analyzed include ^{14}C -labeled casein and ^{125}I -fibrin.

The characteristic of broad protease substrate specificity distinguishes α 2-MG from inhibitors of protease active sites. Based on this characteristic, α 2-MG can be detected or the level of α 2-MG can be determined by assaying for inhibition of the activity of two or more proteases with different active site specificities. α 2-MG can be detected or the level of α 2-MG in a sample can be determined, for example, by analyzing the reduction in protease activity of two or more target proteases such as two or more of the following proteases: trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin and papain. Labeled protease substrates such as ^{14}C -casein or ^{125}I -fibrin can be useful in such methods (Armstrong et al., *supra*, 1999).

α 2-MG also can be detected or the level of α 2-MG determined based on the ability of α 2-MG to shield a bound protease from an antibody or a high molecular

weight inhibitor. A target protease such as trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin, or papain can be added to a plasma sample. Following removal of unbound protease, for example, by

- 5 immunoprecipitation with anti-protease antibody, the amount of protease bound by α 2-MG can be determined using a low molecular mass amide or ester substrate. The amount of hydrolyzed low molecular mass substrate is an indicator of the amount of protected, α 2-MG-bound,
- 10 protease and, therefore, of the concentration of α 2-MG. Similarly, a sample can be reacted first with a protease such as trypsin and subsequently with excess protease inhibitor such as soybean trypsin inhibitor before assaying residual trypsin activity with a low molecular
- 15 mass substrate, such as the amide BApNA (N^{α} -benzoyl-DL-arginine p-nitroanilide (Ganrot, Clin. Chem. Acta 14:493-501 (1966); Armstrong et al., J. Exp. Zool. 236:1-9 (1985)). Trypsin not sequestered by α 2-MG is inactivated by the trypsin inhibitor, with only α 2-MG-protected
- 20 trypsin remaining capable of substrate hydrolysis. Thus, a positive reaction in a soybean trypsin inhibitor assay detects α 2-MG and is a quantitative measure of the amount of α 2-MG (Armstrong et al., *supra*, 1999). One skilled in the art understands that the presence of low molecular
- 25 mass protease inhibitors capable of inactivating α 2-MG-bound enzyme can affect the results obtained with such an assay. It is further understood that these and other routine assays for α 2-MG activity, as well as α 2-MG RNA or protein levels, can be useful for detecting α 2-MG
- 30 or determining a level of α 2-MG in a method of the invention.

200807183-022802

Hyaluronic acid

The methods of the invention further rely, in part, on detecting hyaluronic acid or determining a level of hyaluronic acid in a sample. Hyaluronic acid, also known as hyaluronate or hyaluronan, is a high molecular weight polysaccharide with an unbranched backbone made up of alternating glucuronic acid and $\beta(1,3)$ -N-acetylglucosamine moieties linked by β -1,4 linkages. Hyaluronic acid can have a length of a few to more than 1,000 dimeric units, with each dimeric unit having a molecular weight of about 450 D. Hyaluronic acid, which is produced principally by fibroblasts and other specialized connective tissue cells, plays a structural role in the connective tissue matrix. Furthermore, hyaluronic acid is widely distributed throughout the body and can be found as a free molecule in, for example, plasma, synovial fluid and urine. In plasma, hyaluronic acid has a relatively short half-life.

Serum HA levels can be elevated in liver diseases including cirrhosis (Bramley et al., J. Hepatol. 13:8-13 (1991); Ueno et al., Gastroenterol. 105:475-481 (1993); Oberti et al., Gastroenterol. 113:1609-1616 (1997); and McHutchison et al., J. Gastroenterol. Hepatol. 15:945-951 (2000)). Serum HA levels also can be elevated during synovial inflammation and cartilage destruction seen in rheumatoid arthritis; these levels have been found to correlate with disease activity and degree of synovial involvement (Konttinen et al., Clin. Chimica Acta 193:39-48 (1990); Poole et al., Arthritis Rheum. 37:1030-1038 (1994); Goldberg et al., Arthritis Rheum. 34: 799-807 (1991); and Emlem et al., J.

20087189, 022302

Rheum. 23:974-978 (1996)). Elevated serum levels of HA also can be present, for example, in patients with osteoarthritis (OA), progressive systemic sclerosis (PSS) and systemic lupus erythematosus (SLE).

5 As used herein, the term "hyaluronic acid" is synonymous with "HA" and means a polymer of two or more dimeric units of alternating glucuronic acid and $\beta(1,3)$ -N-acetylglucosamine moieties linked by β -1,4 linkages. As used herein, the phrase "detecting HA" means any quantitative or qualitative assay for
10 determining the presence of HA, and the phrase "determining the level of HA" means any direct or indirect quantitative assay for HA. In view of the above, it is understood that the phrase "detecting HA" encompasses "determining the level of HA."

HA can be detected or a level of HA can be determined using one of a variety of well known assays based on HA-binding proteins or anti-HA antibodies, or by quantitation of purified HA. HA-binding proteins, for
20 example, can be useful in detecting HA; a radiometric assay for HA based on ^{125}I -labelled HA-binding protein is available from Pharmacia (Guehot et al., Clin. Chem. 42:558-563 (1996). Other commercial assays based on HA-binding proteins are available, for example, from
25 Corgenix (Westminster, CO; kit 029001). In addition, HA can be detected or a level of HA can be determined using hyaluronectin as described in Maingonnat and Delpech, Ann. Clin. Biochem. 28:305-306 (1991), or using the kit available from Nalgenunc International (Rochester, NY;
30 Delpech and Bertrand, Anal. Biochem. 149:555-565 (1985)). Assays for detecting HA or determining a level of HA

10087188 022802

- include a variety of competitive and non-competitive binding assays, for example, competitive binding assays using ^{125}I -labeled HA binding protein; competitive binding assays based on alkaline phosphatase
- 5 labeled-hyaluronectin (HN); and non-competitive binding assays based on peroxidase-labeled proteoglycan or peroxidase-labeled HA-binding protein, among others (Lindquist et al., Clin. Chem. 38:127-132 (1992)). See, also, Delpech and Bertrand, *supra*, 1985; Engstrom-Laurent
- 10 et al., Scand. J. Clin. Lab. Invest. 45:497-504 (1985); Brandt et al., Acta Otolaryn. 442 (Suppl.):31-35 (1987); Goldberg, Anal. Biochem. 174:448-458 (1988); Chichibu et al., Clin. Chim. Acta 181:317-324 (1989); Li et al., Conn. Tissue Res. 19:243-254 (1989); Poole et al., Arth.
- 15 Rheum. 33:790-799 (1990); Poole et al., J. Biol. Chem. 260:6020-6025 (1985); and Laurent and Tengblad, Anal. Biochem. 109:386-394 (1980)). Assays for detecting HA or determining a level of HA in a sample can be performed using a variety of immunoassay formats,
- 20 including radioimmunoassays and enzyme-linked immunoassays. Anti-HA antiserum useful in immunoassays can be, for example, affinity purified sheep anti-HA antiserum available from Biotrend (Cologne, Germany; #5029-9990).
- 25 A level of HA also can be determined by purifying HA from a sample, and quantifying the amount of purified polysaccharide. High performance liquid chromatography can be used alone or in conjunction with mass spectrophotometry. As an example, HPLC can be used
- 30 to determine HA levels after digestion of samples containing an internal standard with hyaluronidase, separation by a reversed phase octadecylsilyl column and

20087183 122802

elution with 0.01 M tetrabutylammonium
phosphate-acetonitrile (83:17, v/v) at pH 7.35 (Payan et
al., J. Chromatogr. 566:9-18 (1991)).

HA levels have been shown to correlate with
5 hyaluronidase levels (Bray et al., Am. Rev. Respir.
Dis. 3:284-288 (1991)). Thus, HA can be detected or a
level of HA can be determined indirectly by assaying for
hyaluronidase activity. Assays for hyaluronidase
activity are known in the art, as described in Bray et
10 al., *supra*, 1991. One skilled in the art understands
that these and other routine assays for determining
hyaluonidase or HA levels are encompassed by the phrases
"detecting HA" and "determining the level of HA" and can
be useful in diagnosing the presence or severity of liver
15 fibrosis according to a method of the invention.

TIMP-1

The methods of the invention also are based on
detecting TIMP-1 in a sample and, in particular
embodiments, on determining a level of TIMP-1 in a
20 sample. Tissue inhibitors of metalloproteinases (TIMPs)
regulate the activity of the matrix metalloproteinases
(MMPs), which are an important group of ECM-degradative
enzymes that include gelatinase A (MMP-2) and gelatinase
B (MMP-9). In normal liver, matrix components such as
25 collagens, fibronectin, laminin, tenascin, undulin and
entactin are constantly remodeled by matrix degrading
enzymes to control deposition of extracellular matrix.
Elevation of TIMP levels results in inhibition of MMP
activity and favors the accumulation of extracellular
30 matrix. The TIMPs, which include TIMP-1, TIMP-2, TIMP-3

20082188 "022002

and TIMP-4, interact with the matrix metalloproteinases with a 1:1 stoichiometry and inhibit metalloprotease activity through reversible non-covalent binding. TIMP-1, TIMP-2 and TIMP-3 have similar MMP-inhibitory activities, inhibiting the proteolytic activity of collagenase, gelatinase, stromelysin, proteoglycanase and metalloelastases although their localization and regulation differ (Cawston et al., "Protein Inhibitors of Metalloproteinases" in Barrett and Salvesen (Eds), Proteinase Inhibitors Amsterdam Elsevier pages 589-610 (1986)).

Human TIMP-1 is a 184 amino acid sialoglycoprotein with a molecular weight of 28.5 kDa (Murphy et al., Biochem. J. 195:167-170 (1981); Dockerty et al., Nature 318:66-69 (1985); and Bodden et al., J. Biol. Chem. 269:18943-18952 (1994)). TIMP-1 inhibits all active metalloproteinases, for example, interstitial collagenase MMP-1 as well as stromelysin and gelatinase B (MMP-9). The nucleic acid sequence (SEQ ID NO: 3) and corresponding amino acid sequence (SEQ ID NO: 4) of human TIMP-1 are shown in Figure 2.

As used herein, the term "tissue inhibitor of metalloproteinase-1" is synonymous with "TIMP-1" and means a protein with significant structural homology to human TIMP-1 (SEQ ID NO: 4) that inhibits the proteolytic activity of metalloproteinases with a specificity similar to human TIMP-1. The presence of human TIMP-1 can be conveniently detected by the presence of epitopes reactive with a known specific anti-TIMP-1 antibody such as 7-6Cl or 7-23G9.

As used herein, the phrase "detecting TIMP-1" means any quantitative or qualitative assay for determining the presence of TIMP-1, and the phrase "determining the level of TIMP-1" means any direct or indirect quantitative assay for TIMP-1. In view of the above, it is understood that the phrase "detecting TIMP-1" encompasses "determining the level of TIMP-1."

Assays for detecting TIMP-1 and for determining a level of TIMP-1 include well known assays for TIMP-1 RNA, protein and enzymatic activity. Methods of determining TIMP-1 RNA levels by Northern analysis or RT-PCR are well known in the art (Yoshiji et al., Int. J. Cancer 69:131-134 (1996); Janowska-Wieczorek et al., Exp. Hematol. 28:1274-1285 (2000); and Groft et al., Br. J. Cancer 85:55-63 (2001)) as described further below. TIMP-1 protein can be detected or the level of TIMP-1 protein can be conveniently determined, for example, by radioimmunoassay as described in Brophy et al., Biochem. Biophys. Res. Comm. 167:898-903 (1990) or by two-antibody sandwich assay as described in Murawaki et al., Clinica Chimica Acta 218:47-58 (1993). Plasma concentrations of TIMP-1 protein can be assayed by ELISA with a kit commercially available from Amersham Pharmacia (see, also Example III). Levels of TIMP-1 protein also can be determined by quantifying the amount of purified TIMP-1 protein. Purification of TIMP-1 can be achieved, for example, by HPLC, alone or in combination with mass spectrophotometry, or as described, for example, in Murphy et al., Biochem. J. 195:167-170 (1981), or Stricklin and Welgus, J. Biol. Chem. 258:12252-12258 (1983). TIMP-1 also can be detected or a level of TIMP-1 determined by assaying for inhibition of the activity of

2092220" 8872807

one or more metalloproteases, for example, using reverse gelatin zymography as described in Kossakowska et al., Amer. J. Pathology 153:1895-1902 (1998). Assays for TIMP-1 RNA, protein or activity are described further
5 hereinbelow, and one skilled in the art understands that these and other routine assays for detecting TIMP-1 are encompassed by the methods of the invention.

Rule-in/Rule-out analysis

As disclosed herein, two sets of cut-off values
10 can be used to increase the accuracy of an assay based on the α 2-MG/HA/TIMP-1 three-marker panel. As set forth in Example II, a first set of cut-off values for α 2-MG, HA and TIMP-1 were selected based on optimization for sensitivity in order to first rule out fibrosis, followed
15 by analysis of the "positive" population using a second set of cut-off values optimized for specificity to determine the presence of significant fibrosis. Table 7 shows the results of the dual optimization strategy on the 194 HCV patient study population. The primary
20 cut-offs were set at 2.0 mg/ml, 35 ng/ml and 1000 ng/ml for α 2-MG, HA and TIMP-1, respectively, to achieve a high sensitivity in the primary analysis. Any samples with all three of α 2-MG, HA and TIMP-1 levels above the primary cut-off values were tentatively indicated to be
25 positive for F2-F4 fibrosis and were further evaluated using a second set of cut-off values of 2.0 mg/ml, 60 ng/ml and 1575 ng/ml for α 2-MG, HA and TIMP-1, respectively, which were obtained by optimizing for specificity.

Using the second set of cut-off values optimized for high specificity, 54 of the 122 patients initially designated as positive for F2-F4 fibrosis were confirmed positive, only one of which was a false
5 positive. In sum, of the 194 HCV patients in the study population, 72 were classified as negative (having F0-F1 fibrosis) and 54 were classified as positive (having F2-F4 fibrosis), while 68 samples had indeterminate results and were not classified. When the indeterminate samples
10 were excluded, the α 2-MG/HA/TIMP-1 panel performed with a positive predictive value of about 98% and a negative predictive value of about 79%. Furthermore, in a more typical patient population having 30% fibrosis prevalence, the same panel performs with positive and
15 negative predictive values of close to 93%. These results indicate that the use of primary and secondary cut-off levels, whereby sensitivity is initially optimized followed by optimization for specificity, can increase the overall accuracy of a three-marker test,
20 resulting in a panel test with about 93% accuracy for non-indeterminate samples, which make up about 70% of the samples tested.

Thus, the present invention provides a method of differentiating no or mild liver fibrosis from
25 moderate to severe liver fibrosis in an individual by determining an α 2-MG level in a sample from the individual; determining a HA level in a sample from the individual; determining a TIMP-1 level in a sample from the individual; and diagnosing the individual as having
30 no or mild liver fibrosis when the α 2-MG level is below an α 2-MG cut-off value X1, the HA level is below a HA cut-off value Y1 or the TIMP-1 level is below a TIMP-1

20250320 13:02:30

cut-off value Z1; diagnosing the individual as having moderate to severe liver fibrosis when the α 2-MG level is above an α 2-MG cut-off value X2, the HA level is above a HA cut-off value Y2 and the TIMP-1 level is above a
5 TIMP-1 cut-off value Z2; and diagnosing remaining individuals as having an indeterminate status.

The methods of the invention based on dual cut-off values for the levels of the α 2-MG, HA and TIMP-1 markers can be useful in differentiating no or mild liver
10 fibrosis from moderate to severe liver fibrosis in a variety of patient populations. Such methods can be useful, for example, in diagnosing an individual having a liver disease such as viral hepatitis, autoimmune liver disease such as autoimmune hepatitis, alcoholic liver
15 disease, fatty liver disease or drug-induced liver disease. In one embodiment, a method of the invention is used to differentiate no or mild liver fibrosis from moderate to severe liver fibrosis in an individual infected with hepatitis C virus. Samples useful in a
20 method of the invention based on dual cut-off values include, but are not limited to, blood, serum, plasma, urine, saliva and liver tissue. In one embodiment, a method of the invention is practiced by determining the α 2-MG level, HA level and TIMP-1 level in one or more
25 serum samples.

In a further embodiment, the present invention provides a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual, where the differentiation is based on an X1
30 cut-off value between 1.8 and 2.2 mg/ml; a Y1 cut-off value between 31 and 39 ng/ml; a Z1 cut-off value between

10087188 022802

900 and 1100 ng/ml; an X2 cut-off value between 1.8 and 2.2 mg/ml; a Y2 cut-off value between 54 and 66 ng/ml; and a Z2 cut-off value between 1415 and 1735 ng/ml. In another embodiment, the differentiation is based on an X1
5 cut-off value of 2.0 mg/ml; a Y1 cut-off value of 35 ng/ml; a Z1 cut-off value of 1000 ng/ml; an X2 cut-off value of 2.0 mg/ml; a Y2 cut-off value of 60 ng/ml; and a Z2 cut-off value of 1575 ng/ml. In yet another embodiment, the differentiation is based on an X1 cut-off
10 value of 2.0 mg/ml; a Y1 cut-off value of 37 ng/ml; a Z1 cut-off value of 1100 ng/ml; an X2 cut-off value of 2.0 mg/ml; a Y2 cut-off value of 60 ng/ml; and a Z2 cut-off value of 1575 ng/ml. In a further embodiment, X1, Y1, Z1, X2, Y2 and Z2 are selected such that, in a population
15 having up to 30% liver fibrosis prevalence, at least 65% of individuals in the population are diagnosed as having no or mild fibrosis or moderate to severe fibrosis with an accuracy of at least 90%. In yet a further embodiment, X1, Y1, Z1, X2, Y2 and Z2 are selected such
20 that, in a population having up to 10% liver fibrosis prevalence, at least 70% of individuals in the population are diagnosed as having no or mild fibrosis or moderate to severe fibrosis with an accuracy of at least 90%.

As set forth above, the methods of the
25 invention are highly accurate for determining the presence or severity of fibrosis in a subgroup of the entire patient population assayed. For example, as shown in Table 7, the methods of the invention perform with more than 93% accuracy in determining the F0-F1 or F2-F4
30 fibrosis status in about 70% of a patient population having a liver fibrosis prevalence of 30%. The remaining 30% of the patient population are indicated to have an

indeterminate status. As used herein, the term "indeterminate status" means that the individual cannot be confidently diagnosed with sufficient predictive value.

5 As used herein, the term "X1" or "X2" refers to an α 2-MG cut-off value, against which an experimental α 2-MG sample level is compared. Similarly, as used herein, the term "Y1" or "Y2" refers to an HA cut-off value, against which an experimental HA level is
10 compared. The term "Z1" or "Z2," as used herein, refers to a TIMP-1 cut-off value against which an experimental TIMP-1 level is compared. X1, Y1 and Z1 cut-offs are combined to determine the presence or severity of fibrosis in a sample. Similarly, X2, Y2 and Z1 cut-off
15 values are combined to determine the presence or severity of fibrosis in a sample. A sample having an α 2-MG level less than X1, an HA level less than Y1, or a TIMP-1 level less than Z1 is classified as having F0-F1 fibrosis. A sample having an α 2-MG level above X1, an HA level above
20 Y1, and a TIMP-1 level above Z1 is possibly positive for F2-F4 fibrosis and warrants further analysis. Furthermore, a sample having an α 2-MG level above X2, an HA level above Y2, and a TIMP-1 level above Z2 is classified as having F2-F4 fibrosis. A sample having an
25 α 2-MG level above X1, an HA level above Y1, and a TIMP-1 level above Z1 but one or more levels below X2, Y2 or Z2 is classified as having an "indeterminate status." It is understood that X2 generally is equal to or greater than X1; Y2 generally is equal to or greater than Y1; and Z2
30 generally is equal to or greater than Z1.

10087188-072802

One skilled in the art can select α 2-MG, HA and TIMP-1 cut-offs X1, Y1, Z1, X2, Y2 and Z2 to achieve one or more clinically useful parameters, such as a desired sensitivity or specificity, or a desired negative predictive value, positive predictive value or accuracy for a patient population having a particular fibrosis prevalence. Factorial Design Optimization, also known as Design of Experiments, methodology can be used, for example, to select the appropriate cut-off values. As disclosed herein in Example II, optimization software (DOE Keep It Simple Statistically from Air Academy Associates (Colorado Springs, CO) was used in a central composite design experiment to simultaneously vary the three cut-offs X1, Y1 and Z1, and then to simultaneously vary the three cut-offs X2, Y2 and Z2. In particular, the α 2-MG cut-off was varied from 2.0 to 5.0 mg/ml; the HA cut-off was varied from 25-75 ng/ml; and the TIMP-1 cut-off was varied from 1000-1700 ng/ml. By comparing the test results determined for the 194 patients in the database (see Table 4) with the assigned X1, Y1 and Z1 cut-offs, each of the 194 samples were determined to be a true positive, true negative, false positive or false negative, and the clinical parameters of sensitivity, specificity, negative predictive value, positive predictive value and accuracy were determined for the study patient population. Although determination of the α 2-MG, HA and TIMP-1 cut-off values is illustrated herein using the DOE KISS program, one skilled in the art understands that other computer programs for identifying cooperative interactions among multiple variables and for performing simultaneous equation calculations also can be used. For example, ECHIP optimization software, available from ECHIP, Incorporated (Hockessin, DE), or

Statgraphics optimization software, available from STSC, Incorporated (Rockville, MD), also can be useful in determining α 2-MG, HA and TIMP-1 cut-off values useful in the methods of the invention.

5 The clinical parameters of sensitivity, specificity, negative predictive value, positive predictive value and accuracy are calculated using true positives, false positives, true negatives and false negatives. A "true positive" sample is a sample positive
10 for the indicated stage of fibrosis according to clinical biopsy, which is also diagnosed positive according to a method of the invention. A "false positive" sample is a sample negative for the indicated stage of fibrosis by biopsy, which is diagnosed positive according to a method
15 of the invention. Similarly, a "false negative" is a sample positive for the indicated stage of fibrosis by biopsy, which is diagnosed negative according to a method of the invention. A "true negative" is a sample negative for the indicated stage of fibrosis by biopsy, and also
20 negative for fibrosis according to a method of the invention. See, for example, Motulsky (Ed.), Intuitive Biostatistics New York: Oxford University Press (1995).

As used herein, the term "sensitivity" means the probability that a diagnostic method of the invention
25 gives a positive result when the sample is positive, for example, fibrotic with a Metavir score of F2-F4. Sensitivity is calculated as the number of true positive results divided by the sum of the true positives and false negatives. Sensitivity essentially is a measure of
30 how well a method correctly identifies those with fibrotic disease. In a method of the invention, the X1,

Y1, Z1, X2, Y2 and Z2 values can be selected such that the sensitivity of diagnosing an individual is at least about 70%, and can be, for example, at least 75%, 80%, 85%, 90% or 95% in at least 60% of the patient population
5 assayed, or in at least 65%, 70%, 75% or 80% of the patient population assayed.

As used herein, the term "specificity" means the probability that a diagnostic method of the invention gives a negative result when the sample is not positive,
10 for example, not of Metavir fibrosis stage F2-F4. Specificity is calculated as the number of true negative results divided by the sum of the true negatives and false positives. Specificity essentially is a measure of how well a method excludes those who do not have
15 fibrosis. In a method of the invention, the cut-off values X1, Y1, Z1, X2, Y2 and Z2 can be selected such that, when the sensitivity is at least about 70%, the specificity of diagnosing an individual is in the range of 70-100%, for example, at least 75%, 80%, 85%, 90% or
20 95% in at least 60% of the patient population assayed, or in at least 65%, 70%, 75% or 80% of the patient population assayed. As illustrated in Example II, a specificity of greater than 98% and a sensitivity of about 77% were achieved in the non-indeterminate patient
25 population, which was about 70% of the patient population having a fibrosis prevalence of 30%.

The term "negative predictive value," as used herein, is synonymous with "NPV" and means the probability that an individual diagnosed as not having
30 fibrosis actually does not have the disease. Negative predictive value can be calculated as the number of true

10087188.023802

negatives divided by the sum of the true negatives and false negatives. Negative predictive value is determined by the characteristics of the diagnostic method as well as the prevalence of fibrosis in the population analyzed.

5 In a method of the invention, the α 2-MG, HA and TIMP-1 cut-off values can be selected such that the negative predictive value in a population having a liver fibrosis prevalence of up to 10% is in the range of 75-99% and can be, for example, at least 80%, at least 85%, at
10 least 90%, or at least 95%, in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed. The α 2-MG, HA and TIMP-1 cut-off values also can be selected such that the negative predictive value in a population
15 having a liver fibrosis prevalence of up to 20% is in the range of 75-99% and can be, for example, at least 80%, at least 85%, at least 90%, or at least 95%, in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population
20 assayed. In addition, α 2-MG, HA and TIMP-1 cut-off values can be selected such that the negative predictive value in a population having a liver fibrosis prevalence of up to 30% is in the range of 75-99% and can be, for example, at least 80%, at least 85%, at least 90%, or at
25 least 95%, in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed.

The term "positive predictive value," as used herein, is synonymous with "PPV" and means the
30 probability that an individual diagnosed as having fibrosis actually has the condition. Positive predictive value can be calculated as the number of true positives

5 In a method of the invention, the α 2-MG, HA and TIMP-1 cut-off values can be selected such that, in a patient population having up to 10% liver fibrosis prevalence, the positive predictive value of the method is at least about 75%, and can be at least 80%, at least 85%, at
10 least 90% or at least 95% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed. The α 2-MG, HA and TIMP-1 cut-off values also can be selected such that, in a patient population having up to 20% liver
15 fibrosis prevalence, the positive predictive value of the method is at least about 75%, and can be at least 80%, at least 85%, at least 90% or at least 95% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population
20 assayed. Similarly, the α 2-MG, HA and TIMP-1 cut-off values can be selected such that, in a patient population having up to 30% liver fibrosis prevalence, the positive predictive value of the method is at least about 75%, and can be at least 80%, at least 85%, at least 90% or at
25 least 95% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed.

Predictive values, including negative and positive predictive values, are influenced by the prevalence of the disease in the population analyzed. In the methods of the invention, the cut-off values X_1 , Y_1 , Z_1 , X_2 , Y_2 and Z_2 can be selected to produce a desired

clinical parameter for a clinical population with a particular liver fibrosis prevalence. For example, cut-off values can be selected for a liver fibrosis prevalence of up to 10%, 12%, 15%, 18%, 20%, 25% or 30% which can be seen, for example, in a hepatologist's office. Cut-off values also can be selected for a liver fibrosis prevalence of up to 1%, 2%, 3%, 4%, 5%, 6%, 7% or 8%, which can be representative of the fibrosis prevalence seen in a general practitioner's office.

As used herein, the term "accuracy" means the overall agreement between the diagnostic method and the disease state. Accuracy is calculated as the sum of the true positives and true negatives divided by the total number of sample results and is affected by the prevalence of fibrosis in the population analyzed. The α 2-MG, HA and TIMP-1 cut-off values can be selected such that the accuracy of a method of the invention in a patient population having a liver fibrosis prevalence of up to 10% is at least about 80% and can be, for example, at least 85%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, or 99% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed. The α 2-MG, HA and TIMP-1 cut-off values also can be selected such that the accuracy of a method of the invention in a patient population having a liver fibrosis prevalence of up to 20% is at least about 80% and can be, for example, at least 85%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, or 99% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed. Similarly, the α 2-MG, HA and TIMP-1 cut-off values can be selected such that the

accuracy of a method of the invention in a patient population having a liver fibrosis prevalence of up to 30% is at least about 80% and can be, for example, at least 85%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%,
5 or 99% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed.

Methods not limited to specific markers

The present invention also provides a method of
10 diagnosing the presence or severity of liver fibrosis in an individual by comparing a level of a first fibrotic marker X in the individual to a cut-off value X1 to determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second
15 fibrotic marker Y in the individual to a cut-off value Y1 to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where, in a
20 population with up to 40% fibrosis prevalence, at least 65% of individuals in the population are diagnosed with an accuracy of at least 90%.

A method of the invention can include, if desired, comparing a level of a third fibrotic marker Z
25 in the individual to a cut-off value Z1 to determine whether the individual is positive for the third fibrotic marker Z and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z. In one embodiment, the first

10087188.02802

fibrotic marker is α 2-MG, the second fibrotic marker is HA, and the third fibrotic marker is TIMP-1.

In another embodiment, the levels of at least three fibrotic markers are compared, and, in a further
5 embodiment, the levels of exactly three fibrotic markers are compared to their respective cut-off values. In additional embodiments, the levels of at least four or at least five fibrotic markers are compared. A method of
10 the invention can be useful, for example, to differentiate no or mild liver fibrosis from moderate to severe liver fibrosis.

In a specific embodiment, a method of the invention serves to diagnose at least 65% of individuals in a population with up to 30% fibrosis prevalence with
15 an accuracy of at least 93%. In a further embodiment, a method of the invention serves to diagnose at least 70% of individuals in a population with up to 20% fibrosis prevalence with an accuracy of at least 94%. In yet a further embodiment, a method of the invention serves to
20 diagnose at least 70% of individuals in a population with up to 10% fibrosis prevalence with an accuracy of at least 96%.

The methods of the invention provide unparalleled performance in diagnosing the presence or
25 severity of liver fibrosis. While not all patients are provided with a diagnosis, the majority are diagnosed with extremely good accuracy. As an example, in a patient population with about 40% fibrosis prevalence, almost 70% of the population are diagnosed with more than
30 91% accuracy and with a positive predictive value of more

20087188.022822

than 96% and a negative predictive value of more than 89%. This excellent performance contrasts with alternative methods such as the method of Poynard et al., Lancet 357:1069 (2001). Using the method of Poynard et al. based on analysis of the six markers α 2-MG, α 2-globulin, total bilirubin, γ -globulin, apoA1 and GGT, only about 50% of a population having about 40% fibrosis prevalence are diagnosed, and only with an accuracy of about 89% (see Table 8). Thus, the methods of the invention provide an improvement, in that a significantly greater percentage of a patient population (about 70% as compared to about 50%) are diagnosed, and with an accuracy of more than 91% as compared to an accuracy of around 89% (see Table 8). Due to the novel performance characteristics of a method of the invention, biopsy is typically unnecessary in at least 65% of a patient population, and the patients diagnosed can have confidence in a diagnosis that is more than 90% accurate.

Like other methods of the invention, a method of the invention based on comparison of at least two fibrotic markers can be used to diagnose the presence or severity of liver fibrosis in an individual having or suspected of having any liver disorder, including viral hepatitis, autoimmune liver disease such as autoimmune hepatitis, alcoholic liver disease, fatty liver disease or drug-induced liver disease, or any of the other liver diseases described herein above. Similarly a method of the invention based on comparison of at least two fibrotic markers can be used to diagnose the presence or severity of fibrotic disorders including pulmonary fibrosis, kidney fibrosis, prostate fibrosis, breast

fibrosis or a rheumatoid disease, or another fibrotic disorder described herein or known in the art.

A method of the invention relies on comparison of the level of a fibrotic marker to a predetermined cut-off value. For markers that positively correlate with fibrosis, positivity is indicated by a level that is greater than the predetermined cut-off value. For markers that negatively correlate with fibrosis, positivity is indicated by a level that is less than the predetermined cut-off value. Cut-off values useful in the methods of the invention can be determined as described herein, for example, using design of experiments (DOE) analysis.

As for the other diagnostic methods of the invention, these methods can be practiced using a variety of fibrotic markers known in the art or described herein. Such fibrotic markers include, without limitation, α 2-MG, HA, TIMP-1, PIIINP, laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- β 1, IL-10, apoA1, apoA2 or ApoB. Additional serological, biochemical, clinical and echographic fibrotic markers are described herein above or are known in the art and can be included in any combination in a method of the invention. Furthermore, it is understood that comparison of the first and second fibrotic markers and any additional fibrotic markers can be performed simultaneously or in any order and using any combination of assay formats.

As described above, the "level" of a fibrotic marker can be a relative or absolute amount of, for

example, RNA, protein or activity and can be a direct or indirect measurement of the fibrotic marker. In addition, the value of the level can be obtained from a secondary source, such as a physician or diagnostic laboratory or can be determined using any convenient sample and assay, including but not limited to those described herein above. Methods useful in determining the level of a fibrotic marker in order to perform the comparisons included in the methods of the invention encompass, for example, hybridization methods such as RT-PCR and RNA blot analysis, immunoassays including enzyme-linked immunosorbent assays (ELISAs) and radioimmunoassays (RIAs), sandwich immunoassays, quantitative western blotting and other standard assays for determining protein levels, and, where applicable, assays for the activity of the fibrotic marker. Such assays are routine in the art and described herein above.

The present invention further provides a method of diagnosing the presence or severity of liver fibrosis in an individual by comparing a level of a first fibrotic marker X in the individual to a cut-off value X1 to determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second fibrotic marker Y in the individual to a cut-off value Y1 to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where the cut-off values X1 and Y1 are optimized individually to give a desired performance characteristic.

If desired, a method of the invention can include the steps of comparing a level of a third fibrotic marker Z in the individual to a cut-off value Z1 to determine whether the individual is positive for the third fibrotic marker Z and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z, where the cut-off values X1, Y1 and Z1 are optimized individually to give a desired performance characteristic. In one embodiment, levels of α 2-MG, HA and TIMP-1 are compared. In other embodiments, the levels of exactly three, at least three, at least four, or at least five fibrotic markers are compared. A method of the invention can be useful, for example, in differentiating no or mild liver fibrosis from moderate to severe liver fibrosis. Cut-off values can be optimized as described herein, for example, using DOE analysis.

Further provided by the invention is a method of diagnosing the presence or severity of liver fibrosis in an individual by comparing a level of a first fibrotic marker X in the individual to two cut-off values X1 and X2 to determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second fibrotic marker Y in the individual to two cut-off values Y1 and Y2 to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where the cut-off values X1, Y1, X2 and Y2 are optimized individually to give a desired performance characteristic. Such performance characteristics include

particular sensitivities, specificities, PPVs, NPVs and accuracies, as described herein above.

A method of the invention can further include the steps of comparing a level of a third fibrotic marker Z in the individual to two cut-off values Z1 and Z1 to determine whether the individual is positive for the third fibrotic marker Z; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z, where the cut-off values X1, Y1, Z1, X2, Y2 and Z2 are optimized individually to give a desired performance characteristic. In a method of the invention, cut-off values can be conveniently optimized, for example, using DOE analysis.

15 Methodology

A variety of means can be useful for detecting α 2-MG, HA and TIMP-1 and for determining a level of α 2-MG, HA and TIMP in a sample. In one embodiment, the invention is practiced by determining the level of α 2-MG protein in a sample from the individual to be diagnosed using, for example, one or more α 2-MG-specific binding agents such as anti- α 2-MG antibodies. In another embodiment, a method of the invention is practiced by assaying for α 2-MG activity in a sample from the individual.

A variety of means also can be used in a method of the invention to detect HA or determine a level of HA in a sample. In one embodiment, the invention is practiced by determining the level of HA in a sample

using one or more HA-specific binding agents such as HA-binding proteins or anti-HA antibodies.

Similarly, a variety of means can be used in a method of the invention to detect TIMP-1 or determine a level of TIMP-1 in a sample. In one embodiment, the invention is practiced by determining the level of TIMP-1 protein in a sample from the individual to be diagnosed. The level of TIMP-1 protein can be determined, for example, using one or more TIMP-1-specific binding agents such as anti-TIMP-1 antibodies. In another embodiment, the invention is practiced by assaying for TIMP-1 activity in a sample from the individual to be diagnosed.

In a particular embodiment, the invention provides a method of diagnosing the presence or severity of liver fibrosis in an individual by determining the level of α 2-MG protein in a sample from the individual; determining the level of HA in a sample from the individual; and determining the level of TIMP-1 protein in a sample from the individual; and diagnosing the presence or severity of liver fibrosis in the individual based on the levels of α 2-MG protein, HA and TIMP-1 protein. If desired, the level of α 2-MG protein, HA and TIMP-1 protein each can be determined using an enzyme-linked assay.

In a further embodiment, the present invention provides a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual by contacting an appropriate dilution of a sample from the individual with anti- α 2-MG antibody under conditions suitable to form a first complex of α 2-MG and

anti- α 2-MG antibody; washing the first complex to remove unbound molecules; determining the amount of α 2-MG-containing first complex; contacting an appropriate dilution of a sample from the individual with a HA-binding protein under conditions suitable to form a second complex of HA and HA-binding protein; washing the second complex to remove unbound molecules; determining the amount of HA-containing second complex; contacting an appropriate dilution of a sample from the individual with anti-TIMP-1 antibody under conditions suitable to form a third complex of TIMP-1 and anti-TIMP-1 antibody; washing the third complex to remove unbound molecules; determining the amount of TIMP-1-containing third complex; and differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in the individual based on the amounts of α 2-MG, HA and TIMP-1-containing complexes.

It is understood that detecting α 2-MG, HA and TIMP-1, or detecting α 2-MG, HA and YKL-40, as discussed further below, can be accomplished by assaying for the amount of protein or polysaccharide directly, or, in the case of α 2-MG and TIMP-1, can be determined by assaying for RNA levels or enzyme activity of a protease regulated by α 2-MG or TIMP-1. Similarly, where one or more additional fibrotic markers is detected in a method of the invention, the marker can be assayed directly, or a precursor such as RNA, or a breakdown or proteolytic product, or an activity correlated with levels of the marker can be assayed. It is understood that determining a level of α 2-MG, HA, TIMP-1 and YKL-40, or a level of any additional marker of fibrosis, can be performed using absolute values, for example, for RNA or protein levels

or enzyme activity, or can be determined as relative values in comparison to one or more reference values.

It further is understood that each of the three fibrotic marker assays (α 2-MG/HA/TIMP-1 or
5 α 2-MG/HA/YKL-40), as well as any additional assays, is performed independently of the others, in any order, and that any combination of assay formats is encompassed by the invention. As an example, a level of α 2-MG and HA can be determined by assaying for the concentration of
10 α 2-MG and HA while a level of TIMP-1 is determined by assaying for TIMP-1 enzyme activity. As another example, a level of α 2-MG can be determined using a radioimmunoassay, while levels of HA and TIMP-1 are determined using enzyme-linked assays. One skilled in
15 the art understands that detection of the three fibrotic markers (α 2-MG/HA/TIMP-1 or α 2-MG/HA/YKL-40) and detection of any additional markers can be performed simultaneously or in any order. Furthermore, a single sample such as a serum sample can be obtained from an
20 individual and subdivided into three portions for detecting α 2-MG, HA and TIMP-1 or α 2-MG, HA and TIMP-1, or the markers can be detected using different samples, which can be of the same or a different type and can be undiluted or diluted to the same or different extents.
25 Where two or more samples are used, the samples are usually obtained from the individual within a relatively short time frame, for example, several days to several weeks.

2025 RELEASE UNDER E.O. 14176

RNA methods

Hybridization methods can be used to detect α 2-MG or TIMP-1 mRNA or determine the level of α 2-MG or TIMP-1 mRNA or the mRNA of another fibrotic marker useful in the invention such as YKL-40. Numerous methods are well known in the art for determining mRNA levels by specific or selective hybridization with a complementary nucleic acid probe. Such methods include solution hybridization procedures as well as solid-phase hybridization procedures in which the probe or sample is immobilized on a solid support. Specific examples of useful methods include amplification methods such as target and signal amplification methods and include PCR (polymerase chain reaction) and reverse-transcriptase-PCR (RT-PCR); transcription mediated amplification (Gen-Probe Incorporated; San Diego, CA); branched chain DNA (bDNA) amplification (Bayer Diagnostics; Emeryville, CA); strand displacement amplification (SDA; Becton Dickinson; Franklin Lakes, NJ); and ligase chain reaction (LCR) amplification (Abbott Laboratories; Abbott Park, IL). Additional methods useful in the invention include RNase protection; Northern analysis or other RNA blot, dot blot or membrane-based technology; dip stick; pin; and two-dimensional array immobilized onto a chip. Conditions are well known in the art for quantitative determination of mRNA levels using both solution and solid phase hybridization procedures as described, for example, in Ausubel et al., Current Protocols in Molecular Biology (Supplement 47), John Wiley & Sons, New York (1999).

The polymerase chain reaction (PCR) RT-PCR can be useful in the methods of the invention. PCR or RT-PCR can be performed with isolated RNA or crude or partially fractionated samples, for example, cells pelleted from a whole blood sample. PCR methods are well known in the art as described, for example, in Dieffenbach and Dveksler, PCR Primer: A Laboratory Manual, Cold Spring Harbor Press, Plainview, New York (1995). Multisample formats such as two-dimensional arrays offer the advantage of analyzing numerous different samples in a single assay. Solid-phase dip stick-based methods also can be useful in the invention and offer the advantage of being able to rapidly analyze a fluid sample and obtain an immediate result.

Probes for detecting $\alpha 2$ -MG and TIMP-1 mRNA or for determining $\alpha 2$ -MG and TIMP-1 mRNA levels are well known in the art. One skilled in the art can use, for example, a probe corresponding to some or all of the human $\alpha 2$ -MG nucleic acid sequence shown in Figure 1 (SEQ ID NO: 1) or some or all of the human TIMP-1 nucleic acid sequence shown in Figure 3, respectively. Appropriate conditions for various assay formats for detecting $\alpha 2$ -MG and TIMP-1 mRNA or for determining $\alpha 2$ -MG and TIMP-1 mRNA levels are well known in the art or can be established using routine methods. As an example, conditions and probes for Northern analysis of $\alpha 2$ -MG RNA in human samples are described, for example, in Ortego et al., *supra*, 1997. As another example, conditions and probes for RNA slot blot hybridization to determine $\alpha 2$ -MG RNA expression in human samples are described in Simon et al., *supra*, 1996. Similarly, Northern analysis of TIMP-1 RNA in human samples can be performed as described, for

example, in Yoshiji et al., *supra*, 1996; RT-PCR assays for TIMP-1 in human samples also are well known in the art as described, for example, in Janowska-Wieczorek et al., *supra*, 2000, and Groft et al., *supra*, 2001. The
 5 skilled person understands that these and other assays can be useful for detecting α 2-MG, TIMP-1 or YKL-40 RNA or for determining α 2-MG, TIMP-1 or YKL-40 RNA levels or the levels of other fibrotic markers useful in the methods of the invention.

10 Immunoassays

A variety of immunoassay formats, including competitive and non-competitive immunoassay formats, antigen capture assays and two-antibody sandwich assays also are useful the methods of the invention (Self and
 15 Cook, Curr. Opin. Biotechnol. 7:60-65 (1996)). In one embodiment, a method of the invention relies on one or more antigen capture assays. In an antigen capture assay, antibody is bound to a solid phase, and sample is added such that α 2-MG, HA, TIMP-1, YKL-40 or another
 20 fibrotic marker antigen is bound by the antibody. After unbound proteins are removed by washing, the amount of bound antigen can be quantitated, if desired, using, for example, a radioassay (Harlow and Lane, Antibodies A Laboratory Manual Cold Spring Harbor Laboratory: New
 25 York, 1988)). One skilled in the art understands that immunoassays useful in the invention are performed under conditions of antibody excess, or as antigen competitions, to quantitate the amount of antigen and, thus, determine a level of α 2-MG, HA, TIMP-1 or YKL-40.

30 Enzyme-linked immunosorbent assays (ELISAs) can be useful in the methods of the invention. An enzyme

such as horseradish peroxidase (HRP), alkaline phosphatase (AP), β -galactosidase or urease can be linked, for example, to an anti- α 2-MG, anti-HA, anti-TIMP-1 or anti-YKL-40 antibody or to a secondary
5 antibody for use in a method of the invention. A horseradish-peroxidase detection system can be used, for example, with the chromogenic substrate tetramethylbenzidine (TMB), which yields a soluble product in the presence of hydrogen peroxide that is
10 detectable at 450 nm. Other convenient enzyme-linked systems include, for example, the alkaline phosphatase detection system, which can be used with the chromogenic substrate *p*-nitrophenyl phosphate to yield a soluble product readily detectable at 405 nm. Similarly, a
15 β -galactosidase detection system can be used with the chromogenic substrate *o*-nitrophenyl- β -D-galactopyranoside (ONPG) to yield a soluble product detectable at 410 nm, or a urease detection system can be used with a substrate such as urea-bromocresol purple (Sigma Immunochemicals,
20 St. Louis, MO). Useful enzyme-linked primary and secondary antibodies can be obtained from a number of commercial sources such as Jackson Immuno-Research (West Grove, PA) as described further below.

Chemiluminescent detection also can be useful
25 for detecting α 2-MG, HA, TIMP-1 or YKL-40 or for determining a level of α 2-MG, HA, TIMP-1 or YKL-40 or another fibrotic marker according to a method of the invention. Chemiluminescent secondary antibodies can be obtained commercially from various sources such as
30 Amersham.

200807188.023802

Fluorescent detection also can be useful for detecting α 2-MG, HA, TIMP-1 or YKL-40 or for determining a level of α 2-MG, HA, TIMP-1 or YKL-40 or another fibrotic marker in a method of the invention. Useful
 5 fluorochromes include, without limitation, DAPI, fluorescein, Hoechst 33258, R-phycoerythrin, B-phycoerythrin, rhodamine, Texas red and lissamine. Fluorescein or rhodamine labeled α 2-MG-, HA-, TIMP-1- or YKL-40-specific binding agents such as
 10 anti- α 2-MG, anti-HA, anti-TIMP-1, or anti-YKL-40 antibodies, or fluorescein- or rhodamine-labeled secondary antibodies can be useful in the invention. Useful fluorescent antibodies can be obtained commercially, for example, from Tago Immunologicals
 15 (Burlingame, CA) as described further below.

Radioimmunoassays (RIAs) also can be useful in the methods of the invention. Such assays are well known in the art. For example, Brophy et al., Biochem. Biophys. Res. Comm. 167:898-903 (1990)), describes a
 20 radioimmunoassay for detection of TIMP-1, and Pharmacia makes a radiometric assay for quantitation of HA using an 125 I-labelled HA-binding protein (Guechot et al., Clin. Chem. 42:558-563 (1996)). Radioimmunoassays can be performed, for example, with 125 I-labeled primary or
 25 secondary antibody (Harlow and Lane, *supra*, 1988).

A signal from a detectable reagent can be analyzed, for example, using a spectrophotometer to detect color from a chromogenic substrate; a radiation counter to detect radiation, such as a gamma counter for
 30 detection of 125 I; or a fluorometer to detect fluorescence in the presence of light of a certain wavelength. Where

an enzyme-linked assay is used, quantitative analysis of the amount of α 2-MG, HA, TIMP-1 or YKL-40 or another fibrotic marker can be performed using a spectrophotometer such as an EMAX Microplate Reader

5 (Molecular Devices; Menlo Park, CA) in accordance with the manufacturer's instructions. It is understood that the assays of the invention can be automated or performed robotically, if desired, and that the signal from multiple samples can be detected simultaneously.

10 The methods of the invention also encompass the use of capillary electrophoresis based immunoassays (CEIA), which can be automated, if desired. Immunoassays also can be used in conjunction with laser-induced fluorescence as described, for example, in Schmalzing and

15 Nashabeh, Electrophoresis 18:2184-93 (1997), and Bao, J. Chromatogr. B. Biomed. Sci. 699:463-80 (1997). Liposome immunoassays, such as flow-injection liposome immunoassays and liposome immunosensors, also can be used to detect α 2-MG, HA, TIMP-1 or YKL-40 or to determine a

20 level of α 2-MG, HA, TIMP-1 or YKL-40 or another fibrotic marker according to a method of the invention (Rongen et al., J. Immunol. Methods 204:105-133 (1997)).

Sandwich enzyme immunoassays also can be useful in the methods of the invention. In a two-antibody

25 sandwich assay, a first antibody is bound to a solid support, and the antigen is allowed to bind to the first antibody. The amount of α 2-MG, HA, TIMP-1, YKL-40 or another fibrotic marker antigen is quantitated by measuring the amount of a second antibody that binds the

30 fibrotic marker.

20087188 00000000

As an example, a two-antibody sandwich immunoassay can be useful to determine a level of TIMP-1 as described in Murawaki et al., *supra*, 1993. Briefly, serum (25 μ l) is diluted 41-fold with 10 mM sodium phosphate buffer, pH 7.0 (1.0 ml). The diluted sample (20 μ l) is mixed with 0.3 ml of 10 mM sodium phosphate buffer, pH 7.0, containing 50 ng/ml monoclonal antibody (Fab of clone 7-6C1) labeled with horseradish peroxidase, 1% bovine serum albumin, 0.1% Tween 20, 0.1 M NaCl and 0.005% thimerosal. A 0.1 ml aliquot of the mixed solution is transferred to each microplate well previously coated with a second monoclonal antibody (clone 7-23G9) having a different epitope specificity, and the plate incubated for 30 minutes at room temperature without shaking. The plate is washed three times with 0.3 ml 10 mM sodium phosphate buffer, pH 7.0, containing 0.1% Tween 20 and 0.1 M NaCl. Peroxidase activity bound to the plate is assayed by a 15 minute incubation at room temperature with 0.1 ml 0.15 M citric acid sodium phosphate buffer, pH 4.9, containing 0.5 mg/ml o-phenylenediamine and 0.02% H_2O_2 . After stopping the reaction by addition of 0.1 ml 2 N H_2SO_4 , the absorbance at 492 nm is measured in a microplate reader using a standard of human serum TIMP-1. Linearity between the amount of TIMP-1 and absorbance at 492 nm is demonstrated by graphing with logarithmic scales and yields an assay range of about 1.5 to 300 μ g/well.

Quantitative western blotting also can be used to detect α 2-MG, HA, TIMP-1 or YKL-40 or to determine a level of α 2-MG, HA, TIMP-1 or YKL-40 or a level of another fibrotic marker antigen in a method of the invention. Western blots can be quantitated by well

known methods such as scanning densitometry. As an example, protein samples are electrophoresed on 10% SDS-PAGE Laemmli gels. Primary murine monoclonal antibodies, for example, against human α 2-MG, HA, TIMP-1 or YKL-40 are reacted with the blot, and antibody binding confirmed to be linear using a preliminary slot blot experiment. Goat anti-mouse horseradish peroxidase-coupled antibodies (BioRad) are used as the secondary antibody, and signal detection performed using chemiluminescence, for example, with the Renaissance chemiluminescence kit (New England Nuclear; Boston, MA) according to the manufacturer's instructions. Autoradiographs of the blots are analyzed using a scanning densitometer (Molecular Dynamics; Sunnyvale, CA) and normalized to a positive control. Values are reported, for example, as a ratio between the actual value to the positive control (densitometric index). Such methods are well known in the art as described, for example, in Parra et al., J. Vasc. Surg. 28:669-675 (1998).

Sources for antibodies

As described herein above, immunoassays including but not limited to enzyme-linked immunosorbent assays, radioimmunoassays and quantitative western analysis, can be useful in the diagnostic methods of the invention. Such assays rely on one or more antibodies, for example, anti- α 2-MG, anti-HA, anti-TIMP-1 or anti-YKL-40 antibodies. As used herein, the term "antibody" is used in its broadest sense to include polyclonal and monoclonal antibodies, as well as polypeptide fragments of antibodies that retain binding activity for α 2-MG, HA,

TIMP-1, YKL-40 or the relevant fibrotic marker antigen of at least about $1 \times 10^5 \text{ M}^{-1}$. One skilled in the art understands that antibody fragments such as anti- $\alpha 2$ -MG, anti-HA, anti-TIMP-1 and anti-YKL-40 antibody fragments and including Fab, F(ab')_2 and Fv fragments can retain binding activity for the relevant fibrotic marker antigen and, thus, are included within the definition of the term antibody as used herein. Methods of preparing monoclonal and polyclonal antibodies are routine in the art, as described, for example, in Harlow and Lane, *supra*, 1988.

The term antibody, as used herein, also encompasses non-naturally occurring antibodies and fragments containing, at a minimum, one V_H and one V_L domain, such as chimeric antibodies, humanized antibodies and single chain Fv fragments (scFv) that specifically bind $\alpha 2$ -MG, HA, TIMP-1, YKL-40 or the relevant fibrotic marker antigen. Such non-naturally occurring antibodies can be constructed using solid phase peptide synthesis, produced recombinantly or obtained, for example, by screening combinatorial libraries consisting of variable heavy chains and variable light chains as described by Borrebaeck (Ed.), Antibody Engineering (Second edition) New York: Oxford University Press (1995).

A variety of useful anti- $\alpha 2$ -MG, anti-HA, anti-TIMP-1 and anti-YKL-40 monoclonal and polyclonal antibodies are well known in the art and, in many cases, are commercially available. For example, a nephelometry assay for $\alpha 2$ -macroglobulin is available from Beckman Coulter (kit #449430), and affinity purified goat anti-human $\alpha 2$ -MG and peroxidase-labeled goat anti-human $\alpha 2$ -MG antibodies suitable for ELISA and western blotting

are available, for example, from Cedarlane Laboratories Limited (CL20010AP and CL20010APHP) and Affinity Biologicals Incorporated (GAA2M-AP and GAA2M-APHRP). Similarly, affinity purified sheep anti-HA antiserum can
5 be obtained from Biotrend (#5029-9990).

Anti-human TIMP-1 antibodies also are readily available from a variety of commercial sources. For example, the anti-human TIMP-1 monoclonal antibody 147-6D11 is suitable for ELISA or western blotting
10 analysis and can be obtained from Medicorp, Inc. (Montreal, Canada), and the anti-human TIMP-1 monoclonal antibody MAB970 is available from R&D Systems, Inc., for use, for example, in western blotting or sandwich ELISA assays. MAB970 can be combined, for example, with
15 biotinylated anti-human TIMP-1 antibody (BAF970) from R&D Systems, Inc., for detection of TIMP-1 by sandwich ELISA. In addition, rabbit anti-human TIMP-1 polyclonal antiserum and mouse anti-human monoclonal antibodies suitable, for example, for western blotting with enhanced
20 chemiluminescence detection can be obtained from Research Diagnostics Inc. (RDI-TIMPlabr and RDI-TIMP1-C1).

Assays for activity

As discussed above, assays based on the activity of a fibrotic marker also can be useful for
25 detecting α 2-MG, HA or TIMP-1 or for determining a level of α 2-MG, HA or TIMP-1 or another fibrotic marker and, therefore, are useful in the methods of the invention. As an example, a variety of assays for α 2-MG activity can be useful for detecting α 2-MG or determining a level of
30 α 2-MG in a sample in a method of the invention. Because

α 2-MG-bound proteases display inhibited proteolytic activity but retain the ability to hydrolyze amide and ester bonds of small substrates, α 2-MG can be detected, or a level determined, by assaying for inhibition of
 5 trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin, or papain activity or the activity of another target protease without inhibition of amidolytic activity. Substrates such as labeled casein or labeled fibrin can be useful for assaying for inhibition of
 10 target protease activity. Furthermore, based on its broad protease substrate specificity, a level of α 2-MG can be determined by assaying for inhibition of the activity of two or more target proteases using, for example, ^{14}C -casein and ^{125}I -fibrin (Armstrong et al.,
 15 *supra*, 1999). α 2-MG also can be detected or a level of α 2-MG determined based on the ability of α 2-MG to shield a bound protease from an antibody or a high molecular weight inhibitor. Following reaction of a sample with, for example, trypsin and then trypsin inhibitor, residual
 20 trypsin activity is assayed with a low molecular mass substrate such as the amide BApNA (Ganrot, *supra*, 1966; Armstrong et al., *supra*, 1985). Trypsin activity following treatment with trypsin inhibitor is indicative of α 2-MG. These and other well known assays for α 2-MG
 25 activity can be useful in the methods of the invention.

Similarly, assays for TIMP-1 activity are well known in the art. In particular, one assays for the ability to inhibit protease activity of a matrix metalloproteinase, for example, using reverse gelatin
 30 zymography. Reverse gelatin zymography is performed by including a gelatinase such as gelatinase A in a gel mix with the gelatin substrate. Conditioned media, such as

conditioned media from baby hamster kidney cells can be used as a convenient source of gelatinase. Plasma samples are electrophoresed, and the resulting pattern analyzed, for example, with scanning digitization using a Hewlett Packard scanner. TIMP-1 activity is observed as a reduction of gelatin degradation. See, for example, Kossakowska et al., *supra*, 1998. The skilled person recognizes that these and other routine assays for TIMP-1 activity can be useful in the methods of the invention.

10 Additional markers

It is clear that the methods of the invention can be practiced, if desired, by detecting the three markers α 2-MG, HA and TIMP-1 without assaying for any additional markers or evaluating any other clinical or echographic characteristics. In addition, these three assays can be used as a panel in combination with one or more additional fibrotic marker assays or evaluation of one or more clinical or echographic variables. In specific embodiments, the invention provides a method of diagnosing the presence or severity of liver fibrosis in an individual by detecting α 2-MG, HA and TIMP-1 in a sample and also detecting at least one of the following markers: PIIINP, laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- β 1, IL-10, apoA1, apoA2 or ApoB. In one embodiment, a method of the invention for diagnosing the presence or severity of liver fibrosis includes the steps of detecting α 2-MG, HA, TIMP-1 and YKL-40 in a sample. In a further embodiment, a method of

the invention is limited to detecting α 2-MG, HA, TIMP-1 and YKL-40, and no additional fibrotic markers are detected.

In view of the above, it is clear that assays
 5 for one or more additional biochemical or serological markers of fibrosis or evaluation of one or more clinical or echographic variables associated with fibrosis can be combined with detection of α 2-MG, HA, and TIMP-1 to diagnose the presence or severity of liver fibrosis.
 10 Examples of additional biochemical and serological markers include, yet are not limited to, PIIINP, laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- β 1, IL-10, apoA1, apoA2 and apoB. Additional biochemical and
 15 serological markers useful in the invention include, without limitation, fibronectin, pseudocholinesterase, manganese superoxide dismutase, N-acetyl- β -glucosaminidase (β -NAG), glutathione peroxidase, connective tissue growth factor (CTGF); platelet derived
 20 growth factor (PDGF), PDGF receptor, inducible nitric oxide synthetase, nitrotyrosine, bilirubin, ferritin and α -fetoprotein, γ -glutamyl transpeptidase (GGT), aspartate aminotransferase (AST), alanine aminotransferase (ALT), AST/ALT ratio, albumin, γ -globulins, β γ -block,
 25 prothrombin index, Child-Pugh score, PGA index (prothrombin time, GGT concentration and apoA1 concentration), PGAA index (PGA score with α 2-macroglobulin level), hemoglobin, mean corpuscular volume, lymphocyte count, cholesterol, urea, creatinine,
 30 sodium and platelet count.

A clinical or echographic variable also can be a fibrotic "marker" useful in the methods of the invention. Thus, analysis of one or more clinical or echographic variables can be combined with detection of $\alpha 2$ -MG, HA and TIMP-1 to diagnose the presence or severity of liver fibrosis, or another fibrotic disorder as described hereinabove. As examples, such a clinical variable can be patient age or gender or the presence of palmar erythema, Dupuytren's contracture, finger clubbing, spider nevi, firm liver, splenomegaly or collateral circulation. Echographic variables useful in a method of the invention include, for example, liver length (right kidney), irregular liver surface, liver heterogeneity, spleen length, ascites or collateral circulation. See, for example, Oberti et al., Gastroenterol. 113:1609-1616 (1997). It is understood that the analysis of these and other well known clinical or echographic variables can be useful in a method of the invention. Furthermore, a method of the invention encompasses determination of the clinical or echographic variable, for example, liver palpation, or can rely on one or more historic, or previously determined clinical or echographic variables.

Assays for detection of biochemical or serological markers useful in the invention are well known in the art and in many cases commercially available. Such assays include, but are not limited to, amplification based methods such as RT-PCR and other methods for quantitative analysis of RNA levels; immunoassays such as radioimmunoassays, enzyme-linked assays, two-antibody sandwich assays and quantitative western analysis; and assays for biological activity such

as enzyme activity. Assays for PIIINP, laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- β 1, IL-10, apoA1, apoA2 and apoB are commercially available
5 from various sources as summarized in Table 1.

10067199 000000

TABLE 1
COMMERCIAL SOURCES FOR FIBROTIC MARKER ASSAYS

	Marker	Company	Assay	Catalog number
	PIIINP	Orion Diagnostica (Espoo, Finland)	RIA	05903
5	laminin	Chemicon Intl. (Temecula, CA)	ELISA	ECM310
	tenascin	Chemicon Intl. (Temecula, CA)	ELISA	ECM320
	collagen IV	Iatron Laboratories (Tokyo, Japan)	RIA	KCAD1
	YKL-40	Metra Biosystems (Mountain View, CA)	ELISA	8020
	MMP-3	Amersham Pharmacia (Piscataway, NJ)	ELISA	RPN 2613
10	MMP-2	Amersham Pharmacia (Piscataway, NJ)	ELISA	RPN 2617
	MMP-9/TIMP-1 complex	SBA Sciences (Turku, Finland)	ELISA	MP2215
	sFas ligand	Bender MedSystems Diagnostics (Vienna, Austria)	ELISA	BMS260/2
	TGF- β 1	R&D Systems (Minneapolis, MN)	ELISA	DB100
15	IL-10	R&D Systems (Minneapolis, MN)	ELISA	HS100B
	apoA1	AlerChek, Inc. (Portland, ME)	ELISA	A70101
	apoA2	AlerChek, Inc. (Portland, ME)	ELISA	A70102
	apoB	Sigma Diagnostics (St. Louis, MO)	IT*	357-A
* designates immunoturbidimetric				

20037488-02380

Assays for additional biochemical or serological markers that can be combined with detection of α 2-MG, HA and TIMP-1 in a method of the invention also are well known in the art. Fibronectin, for example, can be conveniently assayed by turbidimetric assay available from Roche Diagnostics (Mannheim, Germany). Pseudocholinesterase (PCHE) can be assayed using standard methodology available from Boehringer. Levels of N-acetyl- β -glucosaminidase (β -NAG) can be determined by assaying for enzymatic activity using a kit available from Cortecs diagnostics. Manganese superoxide dismutase (Mn-SOD) levels can be conveniently determined by ELISA using a kit available, for example, from Bender MedSystem. Glutathione peroxidase levels can be determined by assaying for enzymatic activity using, for example, a kit available from Randox Laboratories Ltd (Oceanside, CA).

Total or direct bilirubin, GGT, AST and ALT levels can be determined using an autoanalyser such as Hitachi 917 Automate (Mannheim, Germany) with Roche Diagnostics reagents. Albumin levels can be determined, for example, by the bromocresol green method as described in Doumas et al., Clin. Chim Acta 31:87-96 (1971), and ferritin and α -fetoprotein levels can be conveniently determined using, for example, an immunoassay available from Boehringer. In addition, levels of α_1 globulin, α_2 globulin, β globulin and γ -globulin can be determined, for example, by serum protein electrophoresis in an automatic system (Hydrasys and Hyrys, Sebia; Issy-Les-Moulineaux, France). Methods of determining prothrombin activity also are well known in the art and include the clotting method available from Organon

Technika (West Orange, NJ). PGA index can be determined as described in Poynard et al., Gastroenterol. 100:1397-1402 (1991), and PGAA index also can be determined by well known methods as described in Naveau et al., Dig.
 5 Dis. Sci. 39:2426-2432 (1994)).

Platelet counts, lymphocyte counts, mean corpuscular volume and related variables can be determined by a variety of methodologies using, for example, a Bayer-Technicon H2 analyser (Bayer-Technicon
 10 Instruments; Tarrytown, NY). Cholesterol levels can be determined by standard methodologies available, for example, from Boehringer. Thus, it is clear to the skilled person that a variety of methodologies, including but not limited to the above, are well known in the art
 15 and can be useful in the diagnostic methods of the invention.

α 2-MG/HA/YKL-40 panel

The present invention also provides a method of diagnosing the presence or severity of liver fibrosis in
 20 an individual by detecting α 2-MG in a sample; detecting HA in a sample; detecting YKL-40 in a sample; and diagnosing the presence or severity of liver fibrosis in the individual based on the presence or level of α 2-MG, HA and YKL-40. A method of the invention can be useful,
 25 for example, or differentiating no or mild (F0-F1) liver fibrosis from moderate to severe (F2-F4) liver fibrosis.

In one embodiment, the invention provides a method of diagnosing the presence or severity of liver fibrosis in an individual by determining the level of

20087183 033300

α 2-MG protein in a sample from the individual;
 determining the level of HA in a sample from the
 individual; determining the level of YKL-40 protein in a
 sample from the individual; and diagnosing the presence
 5 or severity of liver fibrosis in the individual based on
 the levels of α 2-MG protein, HA and YKL-40 protein. If
 desired, the levels of α 2-MG protein, HA and YKL-40
 protein each can be determined using an enzyme-linked
 assay.

10 Thus, in particular embodiments, the present
 invention provides diagnostic methods which rely, in
 part, on determining a level of the fibrotic marker
 YKL-40 in a sample. YKL-40, also known as human
 cartilage glycoprotein 39 (HC gp-39), is named for a
 15 molecular weight of 40 kDa and the amino-terminal
 sequence of the protein, tyrosine-lysine-leucine (YKL).
 This glycoprotein, a mammalian member of the chitinase
 family (18-glycosylhydrolases), is a lectin that binds
 heparin and chitin and is produced by chondrocytes,
 20 synovial cells, activated macrophages, neutrophils and
 MG-63 osteosarcoma cells (Hakala et al., J. Biol. Chem.
 268:25803-15810 (1993); Nyirkos and Golds, Biochem.
J. 268:265-268 (1990); Renkema et al., Eur. J. Biochem.
 251:504-509 (1998); Volck et al., Proc. Assoc. Am.
 25 Physicians 110:351-360 (1998); and Johansen et al., J.
Bone Miner. Res. 7:501-511 (1992)). The pattern of
 YKL-40 expression in normal and diseased tissue indicates
 that this glycoprotein can function in extracellular
 matrix remodelling or tissue inflammation (Nyirkos and
 30 Golds, *supra*, 1990; Renkema et al., *supra*, 1998; and
 Verheijden et al., Arthritis Rheum. 40:1115-1125 (1997)).
 Furthermore, YKL-40 mRNA is expressed in liver, and

initial studies have shown that YKL-40 expression is elevated in patients with chronic liver disease and that increased serum YKL-40 can be associated with fibrosis and fibrogenesis (Johansen et al., Scand. J.

5 Gastroenterol. 32:582-590 (1997); and Johansen, J. Hepatol. 32:911-920 (2000)).

Methods of determining a level of YKL-40 in samples such as serum and synovial fluid are well known in the art. For example, a radioimmunoassay for YKL-40
10 based on a rabbit antibody raised against YKL-40 is described in Johansen et al., Br. J. Rheumatology 32:949-955 (1993). In addition, a sandwich enzyme immunoassay in a microliter stripwell format is commercially available from Metra Biosystems. In the Metra Biosystems
15 assay, the Fab fragment of a biotin-conjugated monoclonal anti-YKL-40 antibody binds to streptavidin on the strip and captures YKL-40 in a sample. Alkaline phosphatase-conjugated polyclonal anti-YKL-40 antiserum binds the captured YKL-40 antigen, and alkaline phosphatase
20 activity is detected with p-nitrophenyl phosphate substrate as an indication of YKL-40 concentration. It is understood that the methods of the invention can be practiced with these or other routine assays for detecting or determining a level of YKL-40 RNA or
25 protein.

The following examples are intended to illustrate but not limit the present invention.

MARKER PANELS FOR NON-INVASIVE DIAGNOSIS
OF LIVER FIBROSIS

This example demonstrates that several
5 serological markers can be combined together as a panel
which is useful in differentiating fibrosis stages F2, F3
and F4 from stages F0 and F1 in HCV-infected patients.

Serum samples from 194 HCV patients positive for hepatitis C virus by RNA and immuno-analysis and having elevated alanine aminotransferase (ALT) levels were chosen at random from an existing serum library. Each of the patients had had a liver biopsy as part of their care. Patient samples were chosen to allow comparison of other routine blood markers and physical examination results attendant to routine medical care, including HCV viral load.

The inclusion criteria for the study were that the patient 1) have a confirmed hepatitis C infection at the time of liver biopsy and serum draw; 2) have
20 undergone a liver biopsy as part of their medical care independent of the study; and 3) have previously given informed consent. Patients who did not give informed consent or who were incarcerated were excluded from the study.

25 Fibrosis scores (Metavir stage) for the 194
patients were established by histopathologic examination
of a needle biopsy specimen prior to therapy according to
the criteria set forth in The French Metavir Cooperative
Study Group, Hepatology, 20:15-20 (1994). All Metavir

fibrosis scores were established by the same pathologist. For all analyses, Metavir scores of F0 and F1 were grouped together as "no/mild" fibrosis, while scores of F2, F3 and F4 were grouped together as "moderate/severe" 5 fibrosis. The fibrosis prevalence in the 194 patient group was determined to be 60% was based on the proportion of F2-F4 scores in the group as shown in Table 2.

10

TABLE 2		
COMPOSITION OF THE 194 HCV PATIENT STUDY		
POPULATION BY FIBROSIS STAGE		
Fibrosis stage	Number	Total F0-F1 or F2-F4
F0	38	F0-F1 = 78
F1	40	
F2	40	F2-F4 = 116
F3	39	
F4	37	
Total	194	Prevalence = 59.8%

15

As shown in the table above, the panel of HCV 20 patient samples included 37 samples with very high fibrosis stage (F4); 39 samples from patients with very low or zero fibrosis stage (F0); and 158 samples from patients with fibrosis stage F1, F2 or F3.

Serum samples were assayed for the presence of 25 several putative fibrosis markers, including laminin, YKL-40, HA, TIMP-1, PIIINP, type IV collagen and α 2-MG. Assays were performed using commercial kits according to manufacturers' instructions (see Table 3). The results obtained for the 194 samples analyzed for laminin,

YKL-40, HA, TIMP-1, PIIINP, collagen type IV and α 2-MG are shown in Table 4.

5

TABLE 3 COMMERCIALLY AVAILABLE KITS FOR DETECTION OF FIBROSIS MARKERS			
Marker	Manufacturer	Assay type	Catalogue number
Laminin	Chemicon Intl. (Temecula, CA)	ELISA	ECM310
YKL-40	Metra Biosystems (Mountain View, CA)	ELISA	8020
HA	Corgenix (Westminster, CO)	ELISA	029001
TIMP-1	Amersham Pharmacia (Piscataway, NJ)	ELISA	RPN 2611
PIIINP	Orion Diagnostica (Espoo, Finland)	RIA	05903
collagen IV	Iatron Laboratories (Tokyo, Japan)	RIA	KCAD1
α 2-MG	Beckman Coulter	Nephelometry	449430

10

20087188 022802
2008220 SET 28001

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100010	B-A	175.244	81.608	15.730	1308.802	2.288	1.737	3.03
100038	P-B	151.888	67.220	9.288	917.104	2.049	2.617	2.01
100044	C-B	187.811	60.757	44.127	1610.690	3.883	3.408	5.03
100059	T-B	232.082	51.002	22.583	1077.343	2.297	1.901	2.29
10 100069	N-C	285.269	131.726	73.851	2381.222	8.034	2.954	4.05
100077	H-C	268.685	47.709	18.066	1122.818	2.260	3.159	1.75
100090	B-C	263.426	26.370	47.339	1380.182	3.526	2.561	3.30
100127	D-D	279.580	166.113	105.505	1180.879	3.343	2.804	3.54
100167	G-F	274.533	482.708	341.132	2523.637	9.745	5.110	3.98
15 100175	B-G	266.903	95.808	27.721	1178.105	4.345	3.911	2.67
100178	M-G	211.613	159.040	25.669	1176.718	2.357	3.795	1.93
100182	T-G	246.686	55.391	8.889	1308.815	2.924	2.468	2.86
100198	A-G	226.372	48.441	13.901	1126.962	2.595	0.819	3.03
100209	J-G	288.524	83.925	5.051	1081.470	5.173	0.801	2.73
20 100229	S-H	253.561	110.020	46.568	1391.433	4.905	4.410	3.12

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100238	T-J	229.781	38.076	29.516	1190.567	2.626	3.141	2.68
100245	D-K	279.768	270.250	171.481	2310.561	5.876	4.713	4.01
100247	C-K	244.559	131.482	10.219	1405.454	2.297	1.089	2.50
100250	J-K	262.136	101.729	54.821	1155.963	4.192	2.655	4.04
5 100252	J-K	260.998	61.366	57.275	1560.856	3.498	7.040	3.24
100253	E-K	292.189	173.917	168.768	1652.033	9.252	7.336	3.59
100254	W-K	288.551	477.560	102.775	1580.756	5.300	4.188	4.01
100271	M-L	278.201	89.900	69.651	1140.761	4.092	4.066	2.94
100276	V-L	257.309	176.112	12.196	1088.369	1.985	4.147	3.57
10 100284	M-L	224.339	130.263	31.822	1104.885	3.653	2.542	3.14
100290	D-L	199.542	541.552	50.429	1550.943	4.399	6.368	3.40
100294	TRL	281.501	217.328	200.436	2340.630	11.006	7.016	4.20
100301	P-M	285.543	430.475	29.772	1884.061	4.453	3.524	4.26
100313	J-M	301.751	188.062	45.539	1852.125	3.610	4.773	4.22
15 100323	M-M	223.002	626.140	144.334	2382.232	8.873	8.042	4.09

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100334	K-R	173.320	63.317	38.516	1290.231	4.327	4.753	2.46
100339	S-M	184.085	36.125	12.049	1268.153	1.865	2.211	2.39
100340	K-M	206.582	57.496	22.015	974.475	2.315	0.819	2.80
100341	T-M	257.580	76.834	91.748	1492.882	3.976	4.793	3.00
5 100343	D-M	334.202	140.629	49.322	1098.199	4.092	2.412	3.12
100357	P-M	419.291	27.368	14.971	784.932	4.225	0.801	2.01
100374	K-N	300.366	28.231	36.608	1697.678	5.478	6.690	3.16
100379	T-O	233.496	75.711	26.906	1437.939	2.086	2.437	3.70
100382	C-P	206.796	44.461	6.034	989.007	2.214	4.409	2.14
10 100397	R-R	223.006	66.474	13.912	981.736	4.091	3.814	3.53
100410	R-S	224.775	36.605	37.499	1152.258	4.592	6.466	3.43
100438	D-Q	228.008	149.349	75.452	1682.636	7.734	7.039	3.89
100451	H-R	248.528	526.840	226.386	1961.029	10.505	7.113	3.36
100453	A-R	225.862	65.956	33.169	1421.632	2.786	3.465	3.79
15 100454	O-R	220.481	56.892	32.531	1125.960	3.003	7.014	3.68

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100456	S-S	241.591	46.274	30.745	1337.355	5.182	2.411	2.81
100466	S-S	210.562	48.605	34.443	1482.475	4.286	4.754	3.27
100470	D-S	229.912	162.039	55.053	1684.159	6.806	3.640	3.10
100485	C-T	229.811	113.523	38.389	1247.547	3.901	3.291	4.11
5 100486	M-T	265.326	281.257	706.557	2716.589	15.362	11.975	2.82
100505	L-W	229.363	33.843	23.305	1149.367	3.397	2.889	3.29
100519	R-W	204.646	68.632	10.431	845.571	3.852	6.815	2.52
100547	S-G	223.959	75.711	8.257	1000.623	4.286	3.090	2.92
100638	J-P	265.819	264.250	68.361	2095.698	15.945	7.418	4.67
10 100640	M-V	170.293	43.770	17.728	1200.584	4.755	5.561	3.50
100006	L-A	135.628	75.349	79.430	1354.782	6.612	3.477	3.14
100009	R-A	157.239	72.429	21.947	932.635	2.080	3.050	1.96
100011	A-B	136.197	251.237	149.932	2004.294	7.600	4.853	3.57
100016	E-A	161.133	272.434	186.536	1900.341	9.341	9.071	3.08
15 100021	E-AV	184.000	537.630	102.420	2456.883	4.863	6.157	3.97

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100023	C-B	126.346	194.523	47.976	1540.914	7.000	4.488	3.21
100027	D-B	133.660	75.820	33.912	1519.528	2.966	3.286	3.51
100030	R-B	140.584	50.007	153.135	1219.549	3.237	5.200	1.89
100035	K-B	124.645	37.383	60.934	1214.060	3.582	3.620	2.41
5 100036	G-B	152.864	87.596	369.681	1305.790	3.163	4.391	2.71
100041	M-B	168.422	42.376	143.372	1502.562	6.667	3.692	3.10
100042	M-B	138.754	211.387	266.568	2899.870	8.233	8.559	3.98
100043	C-B	111.743	30.883	17.447	1168.327	5.488	2.343	2.76
100045	V-B	164.940	241.063	221.249	2010.088	9.097	4.512	3.76
10 100051	K-B	154.743	222.409	131.122	1600.554	4.863	4.075	3.43
100055	D-B	146.817	110.018	84.447	1827.668	3.188	4.439	5.72
100065	G-B	134.349	72.429	112.148	1455.905	3.353	3.002	4.02
100071	R-C	135.011	74.407	30.352	1485.573	2.820	3.120	3.29
100073	G-C	146.785	63.761	43.312	1530.873	3.027	2.040	< 0.75
15 100074	L-C	151.514	80.248	49.917	1647.700	5.036	3.286	3.11

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND $\alpha 2$ -MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha 2$ -MG (mg/ml)
100078	R-C	163.144	213.365	45.839	1399.880	3.393	2.297	2.61
100081	A-C	147.915	45.862	56.686	1346.315	4.056	4.707	2.37
100084	G-C	144.665	43.130	116.238	1736.670	5.337	6.702	< 0.75
100091	P-C	171.782	215.249	33.321	1999.807	5.096	3.835	3.75
5 100093	M-C	133.786	35.499	105.726	1499.707	5.983	5.876	2.85
100099	S-C	174.239	49.159	28.163	1392.574	4.381	5.876	3.18
100100	D-C	181.284	68.095	82.324	1613.489	5.552	5.547	3.71
100103	J-C	151.396	74.849	55.720	1666.282	4.771	3.955	3.57
100104	C-C	128.182	38.890	13.719	1100.784	3.798	5.324	2.21
10 100106	S-C	170.685	41.811	37.449	1280.697	5.096	4.464	2.52
100107	J-C	103.835	27.397	32.334	1416.481	2.910	3.740	2.57
100108	J-C	148.294	145.629	52.822	1884.389	4.484	3.597	4.62
100115	S-DLT	134.784	108.134	96.415	1597.696	8.860	5.225	3.48
100121	R-D	149.335	74.878	34.109	1759.752	5.912	3.716	4.36
15 100124	R-D	134.766	130.367	35.486	1569.219	2.489	4.877	1.21

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND $\alpha 2$ -MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha 2$ -MG (mg/ml)
100125	B-D	170.790	67.078	97.770	2245.776	7.261	5.876	3.63
100126	D-D	134.313	117.116	65.560	1970.476	2.775	3.788	3.04
100129	E-D	159.707	60.388	28.962	1651.995	5.195	4.902	3.41
100131	J-D	155.166	119.774	31.852	1579.186	3.015	3.405	2.73
5 100133	M-D	146.280	24.371	75.729	2098.560	3.225	3.788	1.74
100135	H-E	167.472	41.600	66.767	1369.735	3.200	3.429	4.20
100137	D-E	158.406	25.104	68.740	1346.906	3.828	3.405	2.69
100139	S-E	139.877	38.484	42.708	1388.605	4.215	2.814	4.05
100140	L-E	158.942	30.695	181.056	1585.482	7.476	4.977	3.32
10 100141	W-E	136.761	185.300	179.774	2045.873	9.097	9.872	2.89
100142	R-E	119.383	62.037	16.170	888.744	3.286	2.673	2.07
100143	D-E	131.717	33.779	29.179	1072.170	2.978	3.144	1.91
100147	D-E	132.426	77.159	35.912	1285.138	3.515	2.696	2.12
100150	D-E	120.207	19.056	155.043	1488.729	2.298	4.196	3.66
15 100151	D-F	125.885	35.735	51.625	1243.711	3.447	3.525	3.05

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100155	C-H	146.728	29.137	54.330	1242.340	4.733	3.573	3.15
100159	M-H	136.303	31.336	106.675	1567.716	4.151	3.525	3.51
100161	JF-F	155.052	1710.890	572.598	1966.460	6.226	4.634	4.27
100163	M-F	153.221	785.420	211.173	2167.501	8.269	5.698	3.57
5 100175a	B-G	148.403	55.347	130.093	1282.502	5.296	4.537	2.79
100181	M-G	137.986	69.735	31.119	1384.651	2.813	3.097	1.89
100183	D-G	168.842	181.909	58.358	1499.596	3.101	3.405	3.86
100186	M-G	184.148	2258.120	347.854	5271.196	11.670	7.756	3.82
100200	R-G	148.660	158.906	143.510	1939.499	5.530	6.080	4.02
10 100208	R-G	156.210	94.021	36.624	1379.174	5.339	4.366	< 0.75
100221	ND-H	117.196	38.393	88.913	1375.112	2.610	5.274	3.76
100222	J-H	106.131	34.544	31.603	1054.973	2.580	2.955	3.40
100229a	S-H	125.123	53.139	73.989	1567.731	3.828	3.788	2.78
100237	J-J	140.718	397.625	578.952	1824.407	11.836	6.675	2.46
15 100268	C-L	155.864	76.151	86.977	2060.140	10.963	6.310	3.66

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100270	T-L	176.060	24.738	38.749	1579.990	2.549	2.625	3.35
100278	S-L	153.789	48.840	52.524	1367.051	3.039	3.167	1.91
100279	L-LG	163.352	139.202	34.492	1223.652	2.921	4.099	2.37
100287	R-L	164.414	636.110	232.253	3285.078	14.450	9.448	4.19
5 100291	MS-L	152.863	197.500	42.925	2144.445	2.180	1.993	3.78
100293	D-L	147.479	104.509	27.209	1559.538	3.151	2.696	1.83
100302	A-M	201.715	1021.070	159.330	3317.515	12.498	8.383	4.83
100306	JT-M	125.203	115.289	83.960	1722.987	3.842	2.413	4.16
100307	D-M	128.095	23.612	10.882	1378.118	3.339	3.382	1.42
10 100313a	J-McA	164.201	192.417	76.599	1966.287	3.828	5.523	3.89
100315	K-MF	153.427	113.449	112.430	2118.580	3.515	4.561	4.46
100317	M-McM	165.245	94.693	144.632	1611.495	4.588	3.215	1.84
100320	D-M	159.724	782.150	106.161	1581.345	13.215	5.597	2.33
100322	R-M	120.098	39.914	51.444	1443.789	2.809	3.644	3.21
15 100327	K-R	168.143	194.521	101.231	1738.827	4.295	3.835	4.35

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,

PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100336	ES -M	165.374	135.711	36.329	1556.071	2.932	2.508	3.49
100347	E-M	173.070	69.889	16.945	1710.951	4.808	3.859	3.29
100348	A-M	207.186	75.06	301.583	1334.475	3.299	4.585	3.36
100350	J-M	154.867	4.418	22.250	1388.371	4.087	3.238	1.60
5 100358	A-M	140.022	15.549	88.786	1247.147	4.502	4.682	1.57
100365	A-M	96.324	26.329	43.344	1170.887	5.381	2.040	2.12
100367	B-M	161.274	30.273	46.174	1469.088	3.089	2.790	3.59
100388	A-P	230.782	275.681	938.015	4245.175	20.496	9.669	5.98
100395	D-R	125.908	24.226	15.309	1299.599	2.478	4.415	3.52
10 100397a	R-R	179.186	56.479	100.853	1455.947	9.769	4.172	3.62
100398	C-R	151.391	29.397	11.833	1100.156	2.652	2.932	2.28
100403	L-P	179.146	321.607	350.713	2061.218	8.938	4.977	2.28
100404	ML-P	179.163	1060.240	141.902	2248.495	5.959	4.123	4.07
100414	S-S	184.451	70.941	40.126	1048.761	2.549	3.026	3.18
15 100424	A-P	158.538	62.439	167.519	1320.841	5.509	5.448	3.88

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1, PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100443	J-R	112.348	40.703	21.470	1045.054	2.663	5.647	4.34
100450	M-R	186.892	200.744	203.399	1287.107	2.586	4.040	3.77
100472	T-S	127.877	119.759	21.867	797.753	2.787	3.962	1.97
100474	J-S	118.319	55.427	19.699	939.914	2.287	3.524	3.15
5 100482	J-T	125.011	33.428	39.749	1099.832	5.159	4.386	2.56
100483	J-T	136.978	18.006	33.467	1003.831	2.586	4.485	2.06
100488	M-T	178.106	51.908	123.723	1345.743	4.464	6.053	3.84
100495	J-V	180.283	219.322	93.680	1734.925	8.387	4.064	3.82
100502	D-V	114.380	153.296	22.649	1136.728	2.774	4.757	4.47
10 100503	J-W	149.928	453.095	92.918	1422.332	5.581	4.114	3.15
100513	M-W	117.649	53.027	100.853	1335.730	5.092	4.757	3.74
100528	D-W	189.040	37.248	26.158	1103.495	4.642	3.215	2.90
100530	M-W	106.100	40.234	20.490	924.291	3.297	3.988	2.29
100534	D-A	135.702	37.357	91.410	1421.072	3.374	3.135	3.26
15 100539	M-DB	167.910	474.960	104.470	2158.238	4.957	5.566	3.99

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1, PIIINP, COLLAGEN TYPE IV AND α 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	α 2-MG (mg/ml)
100540	A-B	135.980	54.815	183.589	1881.935	3.068	4.485	3.34
100546	D-F	113.363	75.942	47.682	1207.833	2.633	3.055	1.57
100557	T-L	121.746	58.286	27.704	1297.060	3.842	2.060	1.83
100560	C-N	194.265	142.696	91.041	2338.303	6.332	3.267	4.25
5 100564	D-R	169.241	250.713	65.865	2407.901	5.502	4.534	3.52
100569	J-DC	160.145	64.634	43.744	1349.485	5.962	3.421	3.78
100572	K-K	162.517	260.632	209.581	1729.746	9.292	6.191	2.70
100585	K-Z	171.277	162.336	126.433	2030.404	8.907	4.881	4.71
100594	R-M	114.193	216.295	42.261	1678.540	3.545	3.602	3.26
10 100603	M-S	114.071	61.460	34.332	1693.071	4.464	4.188	2.92
100611	P-F	178.856	269.956	92.721	896.077	3.960	3.524	3.96
100614	J-McA	204.794	245.159	322.970	3470.966	11.393	6.814	5.32
100617	C-W	159.292	51.343	38.850	1504.544	5.859	4.974	3.60
100630	E-AV	140.072	34.778	59.454	1091.420	1.969	3.161	3.34
15 100637	R-B	179.987	59.477	137.723	2077.095	4.726	3.002	4.24

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,
PIIINP, COLLAGEN TYPE IV AND $\alpha 2$ -MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha 2$ -MG (mg/ml)
101013	T-H	177.189	507.415	237.499	1556.336	9.381	8.910	3.36
101118	G-S	163.553	282.057	150.713	2348.845	6.231	7.161	3.93
101137	S-S	175.291	2049.970	715.601	3318.137	11.450	8.458	3.71
101257	M-F	155.324	40.730	49.441	1082.686	3.240	3.421	1.71
5 101275	J-C	121.598	143.292	45.227	1454.509	1.915	2.788	4.69
101284	R-F	123.312	16.825	42.048	1072.891	4.386	3.679	1.64
101321	H-P	180.159	180.091	367.681	2235.931	11.450	10.045	3.37
101322	A-P	133.640	269.262	244.520	2015.508	9.119	8.458	4.01
101335	L-S	164.947	64.238	208.719	2545.531	11.018	6.930	3.90
10 101336	L-S	156.847	741.300	72.664	1387.281	6.097	4.683	3.73
101351	P-F	128.701	12.461	99.050	1557.519	5.486	4.161	4.01
101441	R-H	104.993	176.909	354.512	1338.637	5.271	6.445	3.50
101478	G-S	142.642	63.543	149.266	1296.972	3.859	4.411	2.64
101565	D-A	132.117	74.355	187.250	1206.734	5.962	4.782	2.92

Clinical performance parameters were analyzed for the combinations of markers best able to differentiate the presence of significant fibrosis (F2-F4) from no/mild fibrosis (F0-F1) using various statistical algorithms. The statistical algorithms analyzed included univariate analysis, receiver operating characteristic curves (ROC), logistic regression, discriminant function analysis, and factorial design optimization.

The results of ROC analysis are shown in Table 5. The area under the curve (AUC) values represent relative diagnostic value of a single marker at the indicated cut-off. As can be seen by the decreasing AUC values, HA was shown to have the best diagnostic value when used alone at the indicated cut-off, followed by PIIINP, TIMP-1, α 2-MG and collagen type IV.

TABLE 5				
ROC ANALYSIS				
	AUC	Sensitivity	Specificity	Cutoff
HA	0.821	90.0%	62.0%	35.5 ng/ml
PIIINP	0.777	90.8%	39.2%	3.0 ng/ml
TIMP-1	0.773	90.8%	43.0%	1190.6 ng/ml
α 2-macroglobulin	0.722	90.5%	34.6%	2.4 mg/ml
Collagen Type IV-7S	0.726	90.8%	24.1%	2.79 ng/ml
YKL-40	0.696	90.8%	19.0%	34.5 ng/ml
Laminin	0.524	90.7%	16.5%	125.2 ng/ml

Clinical performance parameters for various combinations of fibrosis markers are shown in Table 6. The best subsets, including single markers as well as combinations of two to four markers and algorithms for discriminating F0-F1 from F2-F4, were generated by logistic regression. The markers included PIIINP, α 2-MG,

laminin and type IV collagen. As shown in Table 5, the diagnostic performance parameters (sensitivity, specificity, PPV and NPV) were similar for the two, three and four-marker combinations identified by logistic regression in the study population, which had a fibrosis prevalence of about 60% (see lines 2-4 and 6-9).

As shown in Table 6, line 5, stepwise discriminant function analysis (SAS) resulted in identification of the 3-marker subset (PIIINP, α 2-MG and laminin). The clinical performance of this combination was similar to the marker combinations identified using logistic regression.

Design of experiments software (DOE KISS, Build 8, Air Academy Associates) was used to simultaneously optimize the cutoffs of multiple variables to obtain the best performance of the panel of tests in predicting fibrosis. Using DOE KISS, a computer-aided central composite design for a combination of markers was generated; this design matrix consisted of a series of combinations of cutoffs for each of the markers in the combination. The results from these experiments (sensitivity, specificity and accuracy) for differentiating F0-F1 from F2-F4 fibrosis were recorded in the design sheet in DOE. Regression analysis was performed for each of the parameters (sensitivity, specificity and accuracy) to give cutoff values for each of the variables in the combination to achieve maximum performance for that parameter.

The five markers with best diagnostic performance in an ROC analysis (highest AUC) were HA,

PIIINP, TIMP-1, α 2-MG and type IV collagen (see Table 5). Cutoffs for each of the markers in this 5-marker panel were optimized for maximum accuracy. The results shown in Table 6, line 10, indicate that, at the optimum

5 accuracy (69.6%), the specificity was too low to be useful (32.9%) while the sensitivity was high (94.8%). Similar results were obtained when the markers were optimized for sensitivity or specificity. Regression analysis showed that TIMP-1 did not have a significant

10 effect on the accuracy, sensitivity or specificity of this 5-marker panel.

20087488 020802

TABLE 6
PERFORMANCE OF VARIOUS MARKER PANELS

Markers	Method/ Model	Sens.	Spec.	Prevalence 59.3%		Prevalence 20%	
				PPV	NPV	PPV	NPV
1	HA						
	N/A	64.96%	82.05%	84.44%	60.95%	47.50%	90.35%
	cut-off 60ug/ml						78.63%
2	PIIINP						
	Logistic	74.36%	58.97%	73.11%	60.53%	31.18%	90.20%
	Best subset of 1						62.05%
3	PIIINP, AMG						
	Logistic	80.53%	63.63%	76.86%	68.54%	31.18%	90.20%
	Best subset of 2						62.05%
4	PIIINP, AMG, Laminin						
	Logistic	78.76%	67.53%	78.44%	67.95%	37.75%	92.71%
	Best subset of 3						69.78%
5	PIIINP, AMG, Laminin						
	Discriminant	78.76%	67.53%	78.44%	67.95%	37.75%	92.71%
	Stepwise selection						69.78%
6	PIIINP, AMG, Laminin, Coll IV						
	Logistic	78.76%	64.94%	77.11%	67.09%	35.96%	92.44%
	Best subset of 4						67.70%
7	PIIINP, AMG, lminin, YKL-40						
	Logistic	77.87%	67.53%	78.25%	67.05%	37.48%	92.43%
	Second best						69.60%
	subset of 4						
8	PIIINP, AMG, Coll IV, YKL-40						
	Logistic	78.76%	70.13%	79.82%	68.76%	39.73%	92.96%
	Third best						71.86%
	subset of 4						
9	PIIINP, AMG, TIMP-1						
	Logistic	78.76%	64.94%	77.11%	67.09%	35.96%	92.44%
	"Forced" selection						67.70%

TABLE 6
PERFORMANCE OF VARIOUS MARKER PANELS

Markers	Method/ Model	Sens.	Spec.	Prevalence 59.3%			Prevalence 20%		
				PPV	NPV	Acc.	PPV	NPV	Acc.
10 HA, PIIINP, AMG Coll IV, TIMP-1	DOE (for acc.) N/A	94.78%	32.91%	67.28%	81.25%	69.59%	26.1%	96.2%	45.3%
11 HA, PIIINP, Coll IV, AMG	DOE (for acc.) N/A	79.13%	75.95%	82.73%	71.43%	77.84%	45.13%	93.57%	76.59%
12 HA, CollIV, AMG	DOE (for acc.) N/A	95.65%	29.11%	66.27%	82.14%	68.56%	25.2%	96.4%	42.4%
13 HA, PIIINP, AMG	DOE (for acc.) N/A	78.26%	75.95%	82.57%	70.59%	77.32%	44.86%	93.32%	76.41%
14 HA, AMG (B)	DOE (for acc.) N/A	84.35%	73.42%	82.20%	76.32%	79.90%	44.24%	94.94%	75.60%
15 TIMP-1, HA, AMG	DOE (for acc.)	83.48%	75.95%	83.48%	75.95%	80.41%	46.46%	94.84%	77.46%
16 YKL-40, HA, AMG	DOE (for acc.)	82.61%	75.95%	83.33%	75.00%	79.90%	46.20%	94.59%	77.28%
17 HA, reflex with Coll IV and AMG	Logistic Positives reflexed	86.32%	70.51%	81.45%	77.46%	80.00%	42.26%	95.38%	73.68%

A similar 4-marker panel was analyzed by DOE as shown in Table 6, line 11. With TIMP-1 excluded, the four-marker panel was optimized for accuracy (77.8%) to give a sensitivity and specificity of 79.1% and 79.5%,
5 respectively. These results demonstrate that the four-marker panel of HA, PIIINP, α 2-MG and collagen IV has more value in differentiating F0-F1 fibrosis from F2-F4 fibrosis than a five-marker panel made up of HA, PIIINP, α 2-MG, collagen IV and TIMP-1.

10 Several three-marker subsets of the four-marker panel also were analyzed by DOE. Line 12 shows the results obtained for the combination of HA, collagen and α 2-MG with the results optimized for accuracy. This three-marker panel gave a very low specificity of less
15 than 30%. In contrast, when a three-marker panel made up of HA, PIIINP and α 2-MG was optimized for accuracy, performance was similar to the four-marker panel (compare lines 13 and 11 of Table 6).

A similar analysis of the two-marker panel of
20 HA and α 2-MG gave the results shown in line 14 of Table 6. This combination gave an improvement in specificity over the three-marker panel of HA, PIIINP and α 2-MG (84.4% compared to 78.3%).

TIMP-1, which was observed to be a good
25 discriminator of fibrosis in the univariate analysis, was added to the two-marker panel. As shown in line 15, the performance of the HA, α 2-MG and TIMP-1 three-marker panel was similar to that obtained with the two-marker panel, and the sensitivity was improved as compared to
30 the three-marker HA/PIIINP/ α 2-MG panel (83.5% sensitivity

20007188, 020002

compared to 78.3%). Furthermore, in preliminary regression analysis, TIMP-1 contributed significantly to discrimination of fibrosis in a study population with a high prevalence of severe fibrosis.

5 Another three-marker panel, made up of HA, α 2-MG and YKL-40, also was optimized for accuracy in differentiating F0-F1 from F2-F4 fibrosis. As shown in Table 6, line 16, this three marker panel had a performance similar to the α 2-MG/HA/TIMP-1 panel.

10 In sum, these results indicate that a α 2-MG/HA/TIMP-1 or α 2-MG/HA/YKL-40 panel can be useful in differentiating F0-F1 from F2-F4 fibrosis.

EXAMPLE II

DUAL OPTIMIZATION STRATEGY FOR ANALYSIS OF THE

15 α 2-MG/HA/TIMP THREE-MARKER PANEL

This example describes the use of multiple cut-offs for α 2-MG-, HA- and TIMP-1 to achieve a relatively high degree of accuracy in a subset of a total patient population assayed.

20 Using the three-marker panel α 2-MG/HA/TIMP-1 with cutoffs for α 2-MG, HA- and TIMP-1 set at 35 ng/ml, 2 mg/ml, and 1000 ng/ml, respectively, samples were determined to be positive when all three of their variables were above the cut-off values, and were
25 therefore negative when one or more of the α 2-MG, HA or TIMP-1 levels were below the assigned cut-off value. As shown in Table 7, in the 194 patient population, there were a total of 72 negative results, 15 of which were

10087188 032302

false negatives, giving a negative predictive value (NPV) of 79% at the study prevalence of about 60% fibrosis%.

- At a prevalence of 30%, which is typical of the prevalence in a hepatology clinic, the negative
- 5 predictive value is over 92%, which is useful in ruling out the presence of F2-F4 fibrosis (likelihood ratio 0.22).

20250303 10:00:00

TABLE 7

PERFORMANCE OF α 2-MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS WITH VARIOUS DISEASE PREVALENCES

Optimization for specificity to rule in fibrosis

X1 = 2.0 mg/ml for α 2-MG Y1 = 60.00 ng/ml for HA Z1 = 1575.00 for TIMP-1 ng/ml

Prevalence 0.598		Prevalence 0.300		Prevalence 0.200		Prevalence 0.100	
Fib + Fib -		Fib + Fib -		Fib + Fib -		Fib + Fib -	
Test +	53 1 54	137 9 146	91 10 102	46 12 57			
Test Equiv	48 20 68	124 179 304	83 205 288	41 231 272			
Sens.	101 21 122	261 188 450	174 215 390	87 242 329			
Spec.	52.48% LR + 11.02	52.48% LR + 11.02	52.48% LR + 11.02	52.48% LR + 11.02			
PPV	95.24% LR - 0.50	95.24% LR - 0.50	95.24% LR - 0.50	95.24% LR - 0.50			
NPV	98.15%	93.86%	89.91%	79.84%			
Accuracy	29.41%	59.11%	71.25%	84.80%			
	59.84%	70.40%	76.12%	83.93%			

TABLE 7

PERFORMANCE OF α 2-MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS
WITH VARIOUS DISEASE PREVALENCES

Final performance after dual optimization

	Prevalence 0.598		Prevalence 0.300		Prevalence 0.200		Prevalence 0.100	
	Fib +	Fib -	Fib +	Fib -	Fib +	Fib -	Fib +	Fib -
	53	1	137	9	91	10	46	12
	63	77	163	691	109	790	54	888
	116	78	300	700	200	800	100	900
	194		1000		1000		1000	
Sens.	45.69%	LR + 35.64	45.69%	LR + 35.64	45.69%	LR + 35.64	45.69%	LR + 35.64
Spec.	98.72%	LR - 0.55	98.72%	LR - 0.55	98.72%	LR - 0.55	98.72%	LR - 0.55
PPV	98.15%		93.86%		89.91%		79.84%	
NPV	55.00%		80.92%		87.91%		94.24%	
Accuracy	67.01%		82.81%		88.11%		93.42%	

TABLE 7

PERFORMANCE OF α 2-MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS WITH VARIOUS DISEASE PREVALENCES

	Population 194 Prevalence 0.598	Population 1000 Prevalence 0.300	Population 1000 Prevalence 0.200	Population 1000 Prevalence 0.100
5	False Negative 15	False Negative 39	False Negative 26	False Negative 13
	False Positive 1	False Positive 9	False Positive 10	False Positive 12
10	Total Incorrect 16	Total Incorrect 48	Total Incorrect 36	Total Incorrect 24
	Total Correct 110	Total Correct 649	Total Correct 676	Total Correct 703
	2nd Round			
15	Fib + Equivocal 48	Fib + Equivocal 124	Fib + Equivocal 83	Fib + Equivocal 41
	Fib - Equivocal 20	Fib - Equivocal 179	Fib - Equivocal 205	Fib - Equivocal 231
20	Total Equivocal 68	Total Equivocal 304	Total Equivocal 288	Total Equivocal 272

TABLE 7

PERFORMANCE OF α 2-MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS WITH VARIOUS DISEASE PREVALENCES

Final performance after dual optimization without equivocal

Prevalence 0.598		Prevalence 0.300		Prevalence 0.200		Prevalence 0.100	
Fib + Fib -		Fib + Fib -		Fib + Fib -		Fib + Fib -	
Test +	53 1 54	137 9 146	91 10 102	46 12 57			
Test -	15 57 72	39 512 550	26 585 610	13 658 671			
	68 58 126	176 521 696	117 595 712	59 669 728			
Sens.	77.94% LR + 45.21	77.94% LR + 45.21	77.94% LR + 45.21	77.94% LR + 45.21			
Spec.	98.28% LR - 0.22	98.28% LR - 0.22	98.28% LR - 0.22	98.28% LR - 0.22			
PPV	98.15%	93.86%	89.91%	79.84%			
NPV	79.17%	92.95%	95.76%	98.07%			
Accuracy	87.30%	93.14%	94.93%	96.64%			
% of pop.		% of pop.		% of pop.			
Test +	27.8%	Test + 14.6%	Test + 10.2%	Test + 5.7%			
Test -	37.1%	Test - 55.0%	Test - 61.0%	Test - 67.1%			
Equiv	35.1%	Equiv 30.4%	Equiv 28.8%	Equiv 27.2%			

10

15

Furthermore, of the 122 test positives using the 35 ng/ml, 2 mg/ml, and 1000 ng/ml cut-offs, 21 of the test positives were false, giving a positive predictive value (PPV) of 82.8%. However, at a more typical
5 prevalence of 30% fibrosis, the positive predictive value falls to about 58% (see Table 7). Thus, in a population with a typical prevalence, a positive result would not have sufficient predictive value to be useful as a diagnostic.

10 In order to increase the positive predictive value for at least a subset of the total patient population, samples positive by the primary analysis were further evaluated for positivity for the three markers using a second set of cut-off values which were higher
15 than the first set. By evaluating those samples positive after a primary analysis at higher cutoffs, the severe fibrosis samples within this group can be determined to be positive with a relatively high predictive value. Those samples that test negative by the secondary
20 evaluation are considered "indeterminate" in that their fibrosis status cannot be determined with good predictive value.

Table 7 shows performance of the α 2-MG/HA/TIMP-1 panel assay with the dual optimization
25 strategy. The primary cut-offs were set at 2.0 mg/ml, 35 ng/ml and 1000 ng/ml to achieve a relatively high sensitivity in the primary analysis. Any samples having all three of their α 2-MG, HA and TIMP-1 levels above the assigned cut-off values were indicated to be positive.
30 The 122 test positives obtained by the primary analysis were re-evaluated using 2.0 mg/ml, 60 ng/ml and 1575

ng/ml as the α 2-MG, HA and TIMP-1 cut-offs and the criteria that the samples must have α 2-MG, HA and TIMP-1 values above the assigned cut-off values to be positive.

Using the second set of cut-off values, 54 of the 122 patients were determined to be positive, only 1 of which was a false positive. The positive predictive value was 98.2% at 59.8% fibrosis prevalence, and was 93.9% at the more typical 30% fibrosis prevalence. In sum, of the 194 patients, 72 were classified as negative and 54 were classified as positive, while 68 samples had indeterminate results and could not be definitively classified. Furthermore, when the indeterminate samples are excluded, the three-marker assay has a positive predictive value of more than 93% and a negative predictive value of close to 93% in a typical population having a 30% fibrosis prevalence.

Table 8 shows a comparison of the performance of the α 2-MG/HA/TIMP-1 three-marker panel with the six marker panel described in Poynard et al., Lancet 357:1069 (2001).

TABLE 8

COMPARISON OF PERFORMANCE OF $\alpha 2$ -MG/HA/TIMP-1 PANEL WITH 6 MARKER PANEL OF POYNARD ET AL.

Prometheus		Biopsy		Poynard et al.	
	Fib +	Fib -		Fib +	Fib -
Test +	53	1	Test + (>.08)	45	5
Test -	15	57	Test - (<0.20)	13	106
Equivocal	48	20	Equivocal	80	90
Total Pop.	116	78	Total Pop.	138	201
					339
Prevalence	.05979		Prevalence	0.4071	
Sensitivity	0.7794		Sensitivity	0.7759	
Specificity	0.9828		Specificity	0.9550	
PPV	0.9815		PPV	0.900	
NPV	0.7917		NPV	0.8908	
Accuracy	0.8730		Accuracy	0.8935	
% Equivoc	.03505	68/194	% Equivoc	.0515	170/339
False Pos	1	of 54 test +	False Pos	5	of 50 test +
False Neg	15	of 72 test -	False Neg	13	of 119 test -
					10.00%
					10.92%
Prometheus		Biopsy		Poynard et al.	
	Fib +	Fib -		Fib +	Fib -
Test +	186	8	Test + (>.08)	133	15
Test -	53	433	Test - (<0.20)	38	313
Equivocal	168	152	Equivocal	236	265
	407	593	Total Pop.	407	593
					1000
Prevalence	0.4071		Prevalence	0.4071	
Sensitivity	0.7794		Sensitivity	0.7759	
Specificity	0.9828		Specificity	0.9550	
PPV	.9607		PPV	0.9000	
NPV	.8917		NPV	0.8907	
Accuracy	.09113		Accuracy	0.8935	
% test pos	19.4%		% test pos	14.7%	
% test neg	48.6%		% test neg	35.1%	
% Equivoc	32.0		% Equivoc	50.1%	
False Pos	8	of 194 test +	False Pos	15	of 147 test +
False Neg	53	of 486 test -	False Neg	38	of 35 test -
					10.00%
					10.93%
* alpha2-macroglobulin, alpha2-globulin, total bilirubin, gama-globulin, apo A1 and GGT					

These results indicate that the $\alpha 2$ -MG/HA/TIMP-1 three-marker panel can be useful in differentiating F0-F1 fibrosis from F2-F4 fibrosis with very good accuracy.

These results further indicate that a combination

- 5 fibrosis marker assay can be useful in determining the fibrosis status of a portion of the patients tested with very good accuracy, while the remaining patients are candidates for biopsy.

EXAMPLE III

10 ASSAYS FOR $\alpha 2$ -MACROGLOBULIN, HYALURONIC ACID AND TISSUE INHIBITOR OF METALLOPROTEINASES-1

A. Quantitation of human $\alpha 2$ -macroglobulin ($\alpha 2$ -MG)

- Serum levels of human $\alpha 2$ -macroglobulin were quantitated using the Beckman Array® 360 System as
15 follows to determine $\alpha 2$ -MG levels in the range of 0.75-270 mg/ml.

- The Beckman Array® 360 system was used for determination of $\alpha 2$ -MG concentrations. This system utilizes a nephelometer which measures the rate of
20 light-scatter formation resulting from an immunoprecipitation reaction between $\alpha 2$ -MG antigen in a sample with antibody to human $\alpha 2$ -MG. After passing a beam of light through the solution in a flow cell, the intensity of light scattered by the formed macromolecular
25 particles of insoluble complexes suspended in solution is detected and measured by the nephelometer. The increase in light scatter resulting from the antigen-antibody reaction is converted to a peak rate signal proportional to the $\alpha 2$ -MG concentration in the sample. The resulting

10087 188 022802
208220 887 28007

formation of complexes and the consequent change in the intensity of scattered light occurs at a rate that increases gradually at first, then rapidly, and finally proceeds through a peak rate of change for the component
5 being analyzed.

Serum samples were drawn from fasting individuals and generally physically separated from cells within 2 hours from the time of collection as set forth in NCCLS publication H 18-A. Samples not assayed within
10 72 hours were stored frozen at -15°C to -20°C . Frozen samples were at most thawed one time. Grossly hemolyzed, highly lipemic or turbid specimens were rejected for further analysis.

Reagents were removed from storage at 4°C and
15 used immediately. Buffers and Diluents were mixed thoroughly by inversion prior to being added to the instrument. Set-up, priming and calibration were performed according to the manufacturer's instructions with samples diluted 1:36. Relatively concentrated
20 samples such as undiluted samples or 1:6 dilutions were generally avoided. Grossly lipemic sample were diluted 1:2 with diluent before assaying. Dust particles or other particulate matter, which can result in extraneous light-scattering signals, in the reaction solution were
25 avoided. Prior to assaying samples, any air bubbles or foam in the sample cups and reagent bottle were removed by using a disposable transfer pipette or pipette tip to aspirate the bubbles. DTT was avoided in the work area.

Samples were analyzed for $\alpha 2$ -MG concentration
30 as follows. The Reagent Wheel (left wheel) on the

10087188, 033302

instrument was loaded with AMG antiserum in space #2. Dilution segments were loaded with 150 μ L of control or sample in the wells on the larger side of the fan shaped segments. Segments and initial dilution control/sample
5 cups were marked for identification. Bubbles were avoided while controls and serum samples were loaded.

Vigil™ Protein Control Levels 1 and 3 (3 drops) was placed in cups 1 and 3, respectively. Biorad Liquichek™ Immunology Control Level 2 (150 μ L) was placed
10 in cup 2. Patient samples (150 μ L) were added to sequential cups. Segments were placed on right wheel beginning at position #1. Evaporation covers were placed over Reagent and Sample Wheels.

On the Master Screen menu, the RESULTS RECALL
15 (F3) was selected before (F4) CLR CUR RUN. After returning to the MASTER SCREEN, the SAMPLE PROGRAM (F1) was selected. ENTER was selected when Reagent wheel #1 appeared and at each cup number. The control ID or sample Acc.# was entered. Test "2" was selected, and
20 SAVE CUP (F1) was selected for each cup. START was selected to begin the analysis. At the end of the run, (Y) was selected in response to CLEAR CURRENT RUN & START NEXT RUN.

Results were reported by the Beckman Array® 360
25 in mg/dl using whole numbers in the Pros System. Samples were diluted routinely by the instrument 1:36. Samples greater than 750 mg/dl were assayed at a 1:216 dilution by the instrument. Samples having a concentration less than 75 mg/dl at a 1:36 dilution are reported as
30 <75 mg/dl. At initial dilutions the Beckman analytical

20087183.022002

range was 75-750 mg/dl, while the extended range was 75-27,000 mg/dl. The range for normal individuals as verified at Prometheus Laboratories was 103-274 mg/dl.

Quality control was performed as follows.

- 5 Three levels of controls were used: low, medium and high. Controls were within 2 standard deviations, except that runs were accepted with two controls within 2 standard deviations and the third control between 2 and 3 standard deviations. The controls used were Beckman Vigil I and
- 10 III and Biorad Level II. Controls were assayed with each sample run.

- 15 The assay is calibrated every 14 days, and also when changes in reagent lots occur or when a major change has occurred in the instrument. Linearity is confirmed every 6 months with appropriate linearity material. This is done to ensure consistent performance over time and to comply with State and National standards.

- 20 Assay calibration verification is performed every 6 months to ensure consistency over time. A minimum of five verification samples including minimum, mid-point, and maximum concentrations are evaluated every 6 months. The coefficient of variation (%CV) of the verification sample results must be less than 15% in order to report out patient sample results.

25 B. Quantitation of Hyaluronic Acid (HA)

Serum levels of HA were determined using the Hyaluronic Acid (HA) Quantitative test kit (Catalog #029001) from Corgenix essentially as follows.

20250718 09:03:02

Serum samples were stored at -70°C . Multiple freeze/thaw cycles were avoided, with a maximum of 4 freeze/thaw cycles per sample. The kits were stored at $2-8^{\circ}\text{C}$.

5 Prior to use, the kit and patient samples were equilibrated to room temperature ($18-28^{\circ}\text{C}$). The pouch of coated strips also was equilibrated to room temperature before opening. Wash solution (0.01 M PBS , $\text{pH } 7.35 \pm 0.1$) was prepared by diluting the 33X PBS wash
10 concentrate with distilled water and adjusting the pH of the final solution to $\text{pH } 7.35 \pm 0.1$.

All blanks, standards, controls and samples were assayed in duplicate. A water blank for calibration of the spectrophotometer was included with each plate and
15 remained empty until addition of $200\text{ }\mu\text{l}$ water immediately prior to reading. Reaction buffer without serum sample was used for the reagent blank, which represented the 0 ng/ml HA reference solution, and was treated the same as patient samples and reference solutions in subsequent
20 assay steps. Three known patient samples (low, middle and high) were run with each assay. In addition, 50 ng/ml HA , 100 ng/ml HA , 200 ng/ml HA , 500 ng/ml HA and 800 ng/ml HA reference solutions supplied with each kit were assayed as described further below.

25 HA reference solutions and patient samples were diluted 1:11 by addition of $25\text{ }\mu\text{l}$ reference solution or sample to $250\text{ }\mu\text{l}$ of reaction buffer and mixed by gentle vortexing. The diluted reference, samples and controls were added ($100\text{ }\mu\text{l}$) to each well. The water blank
30 remained empty. The plate was covered and incubated for

60 minutes at room temperature. After the incubation was complete, the contents of the wells were removed by aspiration. Plates were washed four times with 1X wash solution while avoiding the plates drying out between washes. The plate was blotted vigorously on paper towels to remove residual buffer after the last wash.

HRP-conjugated HA binding protein solution (100 μ l) was added to all wells except the water blank before covering the plate and incubating for 30 minutes at room temperature. After the incubation was complete, the plate was washed four times as described above. Substrate solution (100 μ l 3,3',5,5'-tetramethylbenzidine and hydrogen peroxide, stabilized) was then added to each well except for the water blank. The covered plate was then incubated for 30 minutes at room temperature. The plate was kept in the dark.

The OD₆₅₀ of the 800 ng/ml HA standard was determined. For an OD less than 0.500, the substrate incubation was continued and the OD monitored to determine if the OD had reached 0.500. For an OD greater than 0.500 or after one hour of substrate incubation even if the OD had not reached 0.500, the reactions were terminated by addition of 100 μ l of Stopping Solution (0.36 N sulfuric acid) to each well except the water blank. The stop solution was added in the same order and at approximately the same rate as addition of the substrate solution. Before reading the optical densities, 200 μ l distilled water was added to the water blank. The OD of each well was read at 450 nm (650 nm reference) within one hour after "zeroing" the plate reader against the water blank.

20087188.022802

The following criteria were used to determine if the assay was reliable. The mean OD value of the reagent blank (zero standard) was less than 0.10. Readings greater than 0.10 were considered indicative of possible substrate or reagent contamination, and results were not reported under these conditions. The mean OD value of the 500 ng/ml HA reference was 0.800 or greater. Controls for the three known patient samples were within the following ranges: Low control: 78.6 to 117.2 ng/ml. Mid control: 148.5 to 214.1 ng/ml. High control: 297.8 to 460.7 ng/ml. Samples with HA concentrations greater than 800 ng/ml were further diluted and assayed a second time to obtain a more accurate result.

The known patient controls and samples were determined from a standard 4-parameter curve generated using Softmax and reported in ng/ml. The patient values were not reported if the concentration exceeded the concentration of the highest standard. For patient values greater than the concentration of the highest standard at a 1:11 dilution, samples were assayed at a 1:55 dilution and, if necessary, at higher dilution.

The HA ELISA assay is evaluated every six months to ensure consistent performance over time. A minimum of five samples with previously known HA values are evaluated in a blinded fashion to the operator. For the assay performance to be acceptable, results for negative samples must be negative, and results for positive samples must be positive and yield results within 15% of the previously obtained values. If greater than 20% of the validation samples fail the performance criteria, troubleshooting is implemented, and the assay

is not used to report patient data until acceptable assay performance are reestablished.

C. Quantitation of Tissue Inhibitor of
Metalloproteinases-1 (TIMP-1)

5 Serum levels of TIMP-1 were determined using the Biotrak™ test kit (Catalog# RPN2611) from Amersham Pharmacia Biotech (Piscataway, NJ) essentially as follows.

10 Kit contents were thawed and equilibrated to 20-25°C. Serum samples were stored frozen at -70°C. Repeated freeze-thaw cycles of the samples were minimized, with a maximum of six freeze-thaw cycles.

15 Assay reagents were prepared as follows and stored at 2-8°C for at most 7 days. Assay buffer 1 (0.1 M phosphate buffer, pH 7.5, with 0.9% (w/v) sodium chloride, 0.1% (w/v) BSA and 0.1% Tween-20) was prepared by adding distilled water to the assay buffer concentrate and adjusting the final volume to 100 ml.

20 Anti-TIMP-1 horseradish peroxidase conjugate was prepared in assay buffer 1 essentially as follows. To the stock bottle containing lyophilized conjugate, 11 ml diluted assay buffer 1 was added; the contents were mixed gently until completely dissolved while avoiding vigorous agitation and foaming. Wash buffer (0.1 M phosphate buffer, pH 7.5, containing 0.05% Tween-20) was prepared by adding distilled water to the wash buffer concentrate and bringing the final volume to 500 ml, followed by thorough mixing.

10087188 023302

The 100 ng/ml TIMP-1 stock solution was prepared as follows and stored at 2-8°C. The lyophilized TIMP-1 standard was reconstituted in 0.1 M phosphate buffer, pH 7.5, containing 0.9% (w/v) sodium chloride, 5 0.1% (w/v) bovine serum albumin and 0.1% Tween-20 to make a standard TIMP-1 stock solution of 100 ng/ml. The contents were mixed gently until completely dissolved without vigorous agitation or foaming. Additional standards (1.565, 3.13, 6.25, 12.5, 25 and 50 ng/ml) for 10 a standard curve were prepared fresh before each assay by two-fold serial dilution of the 100 ng/ml stock solution into assay buffer 1 in 1.2 ml dilution tubes. A zero standard (blank) was also prepared.

The pouch containing the microtiter plate was 15 opened after equilibration to room temperature. All samples and standards were assayed in duplicate, and standards for a standard curve were present on each plate. On each plate, seven standards, two controls and a maximum of different 39 samples were present in 20 duplicate.

Samples were diluted 1:120 in tubes by mixing 595 µl assay buffer 1 with 5 µl serum. The dilutions were mixed by vortexing. Using a multichannel pipettor, 100 µl of blank, standards and diluted samples 25 were added to individual wells on a microtiter plate. The plate was covered with the lid provided and incubated at room temperature for exactly two hours. Following the two hour incubation, the contents of the wells were aspirated, and each well was washed four times with wash 30 buffer, with complete filling and aspiration of the wells

10087488 022802

after each wash. After the final wash, the plates were blotted on paper towels to remove residual wash buffer.

Peroxidase conjugate (100 μ l) was added to each well using a multichannel pipettor, and the covered plate
5 incubated at room temperature for exactly two hours. After the incubation, the wells were aspirated and washed as before. Immediately upon conclusion of the incubation, 100 μ l of room temperature equilibrated TMB substrate (3,3',5,5'-tetramethylbenzidine/ hydrogen
10 peroxide in 20% (v/v) dimethylformamide) was added to each well. The plates were covered and incubated for exactly 30 minutes at room temperature. In some cases, the reactions were monitored at 630 nm. The reactions were stopped by addition of 100 μ l 1 M sulfuric acid to
15 all wells. Absorbance was determined at 450 nm within 30 minutes.

Control and patient samples values were determined using a standard curve (4-parameter curve fit) generated using Softmax. Concentration values from the
20 standard curve were multiplied by the dilution factor (120) to obtain actual concentrations, reported in ng/ml. Quality of the assay was confirmed using known serum samples. The low control was in the range of 668.1 to 979.9 ng/ml. The high control was in the range of 2677.9
25 to 3300.2 ng/ml. Patient values generally did not exceed the concentration in ng/ml of the highest standard. Where the patient value was greater than the concentration of the highest standard at a 1:120 dilution, the result was reported as greater than 120
30 times the concentration of the highest standard.

2025071808.023002

The TIMP-1 ELISA assay is validated every six months to ensure consistent performance over time. A minimum of five samples with previously known values are evaluated in a blinded fashion to the operator. Results
5 for negative samples must be negative. Results for positive samples must be positive and must yield results within 15% of the previously obtained values. Where greater than 20% of the validation samples fail the performance criteria, troubleshooting is implemented.
10 Further patient data are not reported until acceptable assay performance is reestablished.

All journal article, reference and patent citations provided above, in parentheses or otherwise, whether previously stated or not, are incorporated herein
15 by reference in their entirety.

Although the invention has been described with reference to the examples provided above, it should be understood that various modifications can be made without departing from the spirit of the invention. Accordingly,
20 the invention is limited only by the claims.

20087188.032802