



1

METHODS OF DIAGNOSING LIVER FIBROSIS

BACKGROUND OF THE INVENTION

FIELD OF THE INVENTION

The present invention relates generally to the  
5 fields of hepatology and fibrosis and, more specifically,  
to a panel of serological markers which together are  
diagnostic of liver fibrosis.

BACKGROUND INFORMATION

Progressive fibrosis of the liver, kidney,  
10 lungs and other organs frequently results in organ  
failure that leads to organ transplantation or death,  
affecting millions in the United States and worldwide.  
Hepatic fibrosis, for example, is the leading  
non-malignant gastrointestinal cause of death in the  
15 United States, and the progression of fibrosis is the  
single most important determinant of morbidity and  
mortality in patients with chronic liver disease.  
Furthermore, the process of fibrosis is common to liver  
diseases of many etiologies, including chronic viral  
20 hepatitis B and C, autoimmune liver disease such as  
autoimmune hepatitis, alcoholic liver disease, fatty  
liver disease; primary biliary cirrhosis; and  
drug-induced liver disease. The fibrosis seen in these  
disorders results from chronic insults to the liver such  
25 as viral infection, alcohol or drugs.

Hepatitis C, for example, is one of the leading  
causes of chronic liver disease in the United States,

where an estimated 3.9 million people are chronically infected with hepatitis C virus (HCV) and approximately 30,000 new cases of acute HCV occur each year (Alter, Semin. Liver Dis. 15:5-14 (1995)). The prevalence of hepatitis C is estimated to be 1.8% in the United States, with up to 10,000 deaths per year likely resulting from chronic hepatitis C infection (Alter, *supra*, 1995).

While hepatic fibrosis is a reversible process resulting in the accumulation of extracellular matrix, liver cirrhosis is an irreversible process characterized by thick bands of matrix which completely encircle the parenchyma to form nodules. Untreated, fibrosis of the liver leads to cirrhosis and eventually end-stage liver disease or cancer. Cirrhosis of the liver is a common condition that frequently goes undetected. For example, in a large sample of the general Danish population, the prevalence of liver cirrhosis was 4.5%, of which one-third were undiagnosed at the time of death (Graudal, J. Intern. Med. 230:165-171 (1991)).

Timely and accurate diagnosis of liver fibrosis is important to effective medical treatment. As an example, patients with hepatitis C and cirrhosis are less likely to respond to treatment with  $\alpha$ -interferon compared to patients with less advanced disease (Davis, Hepatology 26 (Supp. 1):122-127S). Similarly, treatments for chronic HCV infection can be contra-indicated in patients with histologically advanced and decompensated disease (NIH Consensus Development Conference Panel Statement, Hepatology 26 (Suppl. 1):25-105S (1997)). The importance of early diagnosis is further emphasized by the serious early complications such as variceal rupture that are

associated with cirrhosis; these complications can be prevented by early detection of cirrhosis (Calés and Pasqual, Gastroenterol. Clin. Biol. 12:245-254 (1988)).

Diagnosis of the presence or severity of fibrotic liver disease is difficult, with liver biopsy currently the most reliable method available. Unfortunately, liver biopsy has several limitations: pain in about 30% of patients; the risk of severe complications such as hemorrhage or infection; a death rate of 3 in 10,000; and the cost of hospitalization (Nord, Gastrointest. Endosc. 28:102-104 (1982); Cadranel et al., Hepatology 32:47-481 (2000); and Poynard et al., Can. J. Gastroenterol. 14:543-548 (2000)). Furthermore, slowly progressive diseases such as hepatitis C require repeated biopsies for continual assessment of disease progression, thus compounding the risks and costs of the procedure. Finally, biopsy can fail to detect disease because of the heterogeneous distribution of pathological changes in the liver; it is not surprising, then, that false negatives are seen in a significant percentage of cases biopsied (Nord, *supra*, 1982).

For years there has been a search for biochemical or serological markers which reflect fibrotic processes in liver disease and which can serve as a surrogate for liver biopsy. However, the performance of any single marker has not been good enough to substitute for the biopsy procedure in detecting or staging fibrosis. Thus, there is a need for a non-invasive method of diagnosing the presence or severity of liver fibrosis. The present invention satisfies this need by providing a convenient and reliable method for detection

of liver fibrosis that is suitable for serial testing. Related advantages are provided as well.

SUMMARY OF THE INVENTION

The present invention provides a method of  
5 diagnosing the presence or severity of liver fibrosis in  
an individual by detecting  $\alpha$ 2-macroglobulin ( $\alpha$ 2-MG) in a  
sample from the individual; detecting hyaluronic acid  
(HA) in a sample from the individual; detecting tissue  
inhibitor of metalloproteinases-1 (TIMP-1) in a sample  
10 from the individual; and diagnosing the presence or  
severity of liver fibrosis in the individual based on the  
presence or level of  $\alpha$ 2-MG, HA and TIMP-1. A method of  
the invention can be useful, for example, for  
differentiating no or mild (F0-F1) liver fibrosis from  
15 moderate to severe (F2-F4) liver fibrosis.

The methods of the invention for diagnosing the  
presence or severity of liver fibrosis can be useful in a  
variety of patient populations including, but not limited  
to, those with viral hepatitis, autoimmune liver disease  
20 such as autoimmune hepatitis, alcoholic liver disease,  
fatty liver disease and drug-induced liver disease. In  
one embodiment, a method of the invention is used to  
diagnose the presence or severity of liver fibrosis in an  
individual infected with hepatitis C virus.

25 A variety of means can be useful for detecting  
 $\alpha$ 2-MG, HA and TIMP-1 in a sample. In one embodiment, the  
invention is practiced by determining the level of  $\alpha$ 2-MG  
protein in a sample from the individual to be diagnosed  
using, for example, one or more  $\alpha$ 2-MG-specific binding

agents such as anti- $\alpha$ 2-MG antibodies. In another embodiment, a method of the invention is practiced by determining the level of  $\alpha$ 2-MG activity in a sample from the individual.

5 A variety of means also can be used in a method of the invention to detect hyaluronic acid in a sample. In one embodiment, the invention is practiced by determining the level of HA in a sample, for example, using one or more HA-specific binding agents such as  
10 HA-binding proteins (HABPs) or anti-HA antibodies.

Similarly, a variety of means can be used in a method of the invention to detect TIMP-1 in a sample. In one embodiment, the invention is practiced by determining the level of TIMP-1 protein in a sample from the  
15 individual to be diagnosed. The level of TIMP-1 protein can be determined, for example, using one or more TIMP-1-specific binding agents such as anti-TIMP-1 antibodies. In another embodiment, the invention is practiced by assaying for TIMP-1 activity in a sample  
20 from the individual to be diagnosed.

The invention provides, for example, a method of diagnosing the presence or severity of liver fibrosis in an individual by determining the level of  $\alpha$ 2-MG protein in a sample from the individual; determining the  
25 level of HA in a sample from the individual; and determining the level of TIMP-1 protein in a sample from the individual; and diagnosing the presence or severity of liver fibrosis in the individual based on the levels of  $\alpha$ 2-MG protein, HA and TIMP-1 protein. If desired, the

level of  $\alpha$ 2-MG protein, HA and TIMP-1 protein each can be determined using an enzyme-linked assay.

A variety of samples can be useful in practicing the methods of the invention including, for 5 example, blood, serum, plasma, urine, saliva and liver tissue. In one embodiment, a single sample is obtained from the individual to be diagnosed. Such a sample can be, for example, a serum sample. Such a sample also can be, for example, a tissue sample, for example, a liver 10 biopsy sample.

The present invention further provides a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual. The method includes the steps of contacting an appropriate 15 dilution of a sample from the individual with anti- $\alpha$ 2-MG antibody under conditions suitable to form a first complex of  $\alpha$ 2-MG and anti- $\alpha$ 2-MG antibody; washing the first complex to remove unbound molecules; determining the amount of  $\alpha$ 2-MG-containing first complex; contacting 20 an appropriate dilution of a sample from the individual with a HA-binding protein under conditions suitable to form a second complex of HA and HA-binding protein; washing the second complex to remove unbound molecules; determining the amount of HA-containing second complex; 25 contacting an appropriate dilution of a sample from the individual with anti-TIMP-1 antibody under conditions suitable to form a third complex of TIMP-1 and anti-TIMP-1 antibody; washing the third complex to remove unbound molecules; determining the amount of 30 TIMP-1-containing third complex; and differentiating no or mild liver fibrosis from moderate to severe liver

fibrosis in the individual based on the amounts of  $\alpha$ 2-MG, HA and TIMP-1-containing complexes.

The methods of the invention can be practiced by detecting the three markers  $\alpha$ 2-MG, HA and TIMP-1,

5 without detecting additional serological markers, or can be combined with a detection method for one or more additional markers. Thus, in one embodiment, the invention is practiced by detecting  $\alpha$ 2-MG, HA and TIMP-1 and also detecting at least one of the following markers

10 of fibrosis: N-terminal procollagen III propeptide (PIIINP), laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- $\beta$ 1, IL-10, apoA1, apoA2 or apoB. In a further embodiment, the presence or severity of liver

15 fibrosis is diagnosed by detecting  $\alpha$ 2-MG, HA, TIMP-1 and YKL-40 in a sample from an individual.

The present invention also provides a method of monitoring the efficacy of anti-fibrotic therapy in a patient by detecting  $\alpha$ 2-macroglobulin in a sample from a

20 patient administered an anti-fibrotic therapy; detecting hyaluronic acid (HA) in a sample from the patient; detecting tissue inhibitor of metalloproteinases-1 (TIMP-1) in a sample from the patient; and determining the presence or severity of liver fibrosis in the patient

25 based on the presence or level of  $\alpha$ 2-MG, HA and TIMP-1, thereby monitoring the efficacy of anti-fibrotic therapy. Such a method can further include, if desired, comparing the presence or severity of liver fibrosis determined in step (d) to the presence or severity of liver fibrosis in

30 the patient at an earlier time. The methods of the invention can be used to monitor, for example, the

progression or regression of fibrosis over time in a patient treated with one or more anti-fibrotic therapies, or to compare, for example, the efficacies of two or more anti-fibrotic therapies.

- 5           In one embodiment, at most three markers of fibrosis are detected. In another embodiment, the method includes the step of detecting in a sample from the patient at least one marker selected from the group consisting of: PIIINP, laminin, tenascin, collagen  
10 type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- $\beta$ 1, IL-10, apoA1, apoA2, and apoB.

A variety of means can be useful for detecting  $\alpha$ 2-MG, HA and TIMP-1 in a method of the invention. Step  
15 (a) can be practiced, for example, by determining the level of  $\alpha$ 2-MG protein in the sample. In one embodiment, the level of  $\alpha$ 2-MG protein is determined using one or more anti- $\alpha$ 2-MG antibodies. Step (b) can be practiced, for example, by determining the level of HA in the  
20 sample. In one embodiment, the level of HA is determined using one or more HA-binding proteins. Step (c) can be practiced, for example, by determining the level of TIMP-1 protein in said sample. In one embodiment, the level of TIMP-1 protein is determined using one or more  
25 anti-TIMP-1 antibodies.

Further provided herein is a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual by determining an  $\alpha$ 2-MG level in a sample from the individual;  
30 determining a HA level in a sample from the individual;

determining a TIMP-1 level in a sample from the individual; and diagnosing the individual as having no or mild liver fibrosis when the  $\alpha$ 2-MG level is below an  $\alpha$ 2-MG cut-off value X1, the HA level is below a HA 5 cut-off value Y1 or the TIMP-1 level is below a TIMP-1 cut-off value Z1; diagnosing the individual as having moderate to severe liver fibrosis when the  $\alpha$ 2-MG level is above an  $\alpha$ 2-MG cut-off value X2, the HA level is above a HA cut-off value Y2 and the TIMP-1 level is above a 10 TIMP-1 cut-off value Z2; and diagnosing remaining individuals as having an indeterminate status.

The methods of the invention based on dual cut-off values for the levels of the  $\alpha$ 2-MG, HA and TIMP-1 markers can be useful in differentiating no or mild liver 15 fibrosis from moderate to severe liver fibrosis in a variety of patient populations. The methods of the invention can be useful, for example, in diagnosing an individual having a liver disease such as viral hepatitis, autoimmune liver disease such as autoimmune 20 hepatitis, alcoholic liver disease, fatty liver disease or drug-induced liver disease. In one embodiment, the methods of the invention are used to differentiate no or mild liver fibrosis from moderate to severe liver fibrosis in an individual infected with hepatitis C 25 virus. Samples useful in the methods of the invention include, but are not limited to, blood, serum, plasma, urine, saliva and liver tissue. In one embodiment, a method of the invention is practiced by determining the  $\alpha$ 2-MG level, HA level and TIMP-1 level in one or more 30 serum samples from the individual to be diagnosed.

Thus, the present invention provides, for example, a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual in which the differentiation is based on an X1 5 cut-off value between 1.8 and 2.2 mg/ml; a Y1 cut-off value between 31 and 39 ng/ml; a Z1 cut-off value between 900 and 1100 ng/ml; an X2 cut-off value between 1.8 and 2.2 mg/ml; a Y2 cut-off value between 54 and 66 ng/ml; and a Z2 cut-off value between 1415 and 1735 ng/ml. In a 10 particular embodiment, the differentiation is based on an X1 cut-off value of 2.0 mg/ml; a Y1 cut-off value of 35 ng/ml; a Z1 cut-off value of 1000 ng/ml; an X2 cut-off value of 2.0 mg/ml; a Y2 cut-off value of 60 ng/ml; and a Z2 cut-off value of 1575 ng/ml. In another embodiment, 15 the differentiation is based on an X1 cut-off value of 2.0 mg/ml; a Y1 cut-off value of 37 ng/ml; a Z1 cut-off value of 1100 ng/ml; an X2 cut-off value of 2.0 mg/ml; a Y2 cut-off value of 60 ng/ml; and a Z2 cut-off value of 1575 ng/ml. In a further embodiment, X1, Y1, 20 Z1, X2, Y2 and Z2 are selected such that, in a population having up to 30% liver fibrosis prevalence, at least 65% of individuals in the population are diagnosed as having no/mild fibrosis or moderate/severe fibrosis with an accuracy of at least 90%. In another embodiment, X1, Y1, 25 Z1, X2, Y2 and Z2 are selected such that, in a population having up to 30% liver fibrosis prevalence, at least 65% of individuals in said population are diagnosed as having no/mild fibrosis or moderate/severe fibrosis with a positive predictive value of at least 90% and a negative 30 predictive value of at least 90%. In yet a further embodiment, X1, Y1, Z1, X2, Y2 and Z2 are selected such that, in a population having up to 10% fibrosis prevalence, at least 70% of individuals in the population

are diagnosed as having no/mild fibrosis or moderate/severe fibrosis with an accuracy of at least 90%.

The present invention also provides a method of  
5 diagnosing the presence or severity of liver fibrosis in  
an individual by comparing a level of a first fibrotic  
marker X in the individual to a cut-off value X<sub>1</sub> to  
determine whether the individual is positive for the  
first fibrotic marker X; comparing a level of a second  
10 fibrotic marker Y in the individual to a cut-off value Y<sub>1</sub>  
to determine whether the individual is positive for the  
second fibrotic marker Y; and diagnosing the presence or  
severity of liver fibrosis in the individual based on  
positivity or negativity for X and Y, where, in a  
15 population with up to 40% fibrosis prevalence, at least  
65% of individuals in the population are diagnosed with  
an accuracy of at least 90%.

A method of the invention can include, if  
desired, comparing a level of a third fibrotic marker Z  
20 in the individual to a cut-off value Z<sub>1</sub> to determine  
whether the individual is positive for the third fibrotic  
marker Z and diagnosing the presence or severity of liver  
fibrosis in the individual based on positivity or  
negativity for X, Y and Z. In one embodiment, the first  
25 fibrotic marker is  $\alpha$ 2-MG, the second fibrotic marker is  
HA, and the third fibrotic marker is TIMP-1.

In another embodiment, the levels of at least  
three fibrotic markers are compared, and, in a further  
embodiment, the levels of exactly three fibrotic markers  
30 are compared. In additional embodiments, the levels of

at least four or at least five fibrotic markers are compared. A method of the invention can be useful, for example, to differentiate no or mild liver fibrosis from moderate to severe liver fibrosis.

5           In a specific embodiment, a method of the invention serves to diagnose at least 65% of individuals in a population with up to 30% fibrosis prevalence with an accuracy of at least 93%. In a further embodiment, a method of the invention serves to diagnose at least 70%  
10 of individuals in a population with up to 20% fibrosis prevalence with an accuracy of at least 94%. In yet a further embodiment, a method of the invention serves to diagnose at least 70% of individuals in a population with up to 10% fibrosis prevalence with an accuracy of at  
15 least 96%.

The present invention further provides a method of diagnosing the presence or severity of liver fibrosis in an individual by comparing a level of a first fibrotic marker X in the individual to a cut-off value X<sub>1</sub> to  
20 determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second fibrotic marker Y in the individual to a cut-off value Y<sub>1</sub> to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or  
25 severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where the cut-off values X<sub>1</sub> and Y<sub>1</sub> are optimized individually to give a desired performance characteristic.

If desired, a method of the invention can  
30 include the steps of comparing a level of a third

fibrotic marker Z in the individual to a cut-off value Z1 to determine whether the individual is positive for the third fibrotic marker Z and diagnosing the presence or severity of liver fibrosis in the individual based on 5 positivity or negativity for X, Y and Z, where the cut-off values X1, Y1 and Z1 are optimized individually to give a desired performance characteristic. In one embodiment, levels of  $\alpha$ 2-MG, HA and TIMP-1 are compared. In another embodiment, the cut-off values are optimized 10 using design of experiments (DOE) analysis. In further embodiments, the levels of exactly three, at least three, at least four, or at least five fibrotic markers are compared. A method of the invention can be useful, for example, in differentiating no or mild liver fibrosis 15 from moderate to severe liver fibrosis.

Further provided by the invention is a method of diagnosing the presence or severity of liver fibrosis in an individual by comparing a level of a first fibrotic marker X in the individual to two cut-off values X1 and 20 X2 to determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second fibrotic marker Y in the individual to two cut-off values Y1 and Y2 to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing 25 the presence or severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where the cut-off values X1, Y1, X2 and Y2 are optimized individually to give a desired performance characteristic. A method of the invention can further 30 include the steps of comparing a level of a third fibrotic marker Z in the individual to two cut-off values Z1 and Z1 to determine whether the individual is positive

for the third fibrotic marker Z; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z, where the cut-off values X<sub>1</sub>, Y<sub>1</sub>, Z<sub>1</sub>, X<sub>2</sub>, Y<sub>2</sub> and Z<sub>2</sub> are  
5 optimized individually to give a desired performance characteristic. Cut-off values can be conveniently optimized, for example, using DOE analysis.

#### BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1 shows the nucleic acid sequence (SEQ ID NO: 1) and corresponding amino acid sequence (SEQ ID NO: 2) for mature human  $\alpha$ 2-macroglobulin available from Genbank accession M36501.  
10

Figure 2 shows the nucleic acid sequence (SEQ ID NO: 3) and corresponding amino acid sequence (SEQ ID NO: 4) for human tissue inhibitor of metalloproteinases-1 (TIMP-1) available from Genbank accession NM\_003254.  
15

#### DETAILED DESCRIPTION OF THE INVENTION

As disclosed herein, the serum levels of a number of biochemical markers were analyzed in a patient population with confirmed hepatitis C and having a known Metavir stage (fibrosis score) of F0 to F4, where F0 represents very low or no fibrosis; F1, F2 and F3 represent intermediate fibrosis stages; and F4 represents severe fibrosis (Knodell et al., *Hepatology* 1:431-435 (1981)). See Tables 2 and 3. Using Design of Experiments (DOE) analysis for simultaneous variation of multiple cut-off values, a four-marker panel made up of hyaluronic acid (HA), PIIINP, collagen type IV and

$\alpha_2$ -macroglobulin ( $\alpha_2$ -MG) was identified which was capable of differentiating F0-F1 (no or mild) fibrosis from F2-F4 (moderate to severe) fibrosis with an accuracy of about 77% in a patient population with a fibrosis prevalence of 60%.

As further disclosed herein in Example I, two three-marker panels,  $\alpha_2$ -MG/HA/TIMP-1 and  $\alpha_2$ -MG/HA/YKL-40, also performed well in differentiating F0-F1 fibrosis from F2-F4 fibrosis when cut-offs were optimized using 10 DOE analysis. In particular, the  $\alpha_2$ -MG/HA/TIMP-1 and  $\alpha_2$ -MG/HA/YKL-40 panels each performed better than the four-marker panel and were capable of differentiating F0-F1 from F2-F4 fibrosis with about 80% accuracy in the study population. As can be seen in Table 6, line 15, 15 for example, the  $\alpha_2$ -MG/HA/TIMP-1 panel performed with a sensitivity of 83.48% and a specificity of 75.95% in the study population having 60% fibrosis prevalence. These results demonstrate that the  $\alpha_2$ -MG/HA/TIMP-1 three-marker panel can be useful for differentiating no or mild 20 fibrosis from moderate to severe fibrosis.

Based on these findings, the present invention provides a method of diagnosing the presence or severity of liver fibrosis in an individual by detecting  $\alpha_2$ -MG in a sample from an individual; detecting HA in a sample 25 from the individual; detecting TIMP-1 in a sample from the individual; and diagnosing the presence or severity of liver fibrosis in the individual based on the presence or levels of  $\alpha_2$ -MG, HA and TIMP-1. A method of the invention can be useful, for example, for differentiating 30 no or mild (F0-F1) liver fibrosis from moderate to severe (F2-F4) liver fibrosis.

Liver and other fibrotic disorders

The methods of the invention can be useful for diagnosing the presence or severity of liver fibrosis in a variety of individuals including those at risk for, or

5 having one or more symptoms of, a liver disorder characterized by fibrosis. The methods of the invention can be used to diagnose liver fibrosis in an individual having, for example, viral hepatitis such as hepatitis A, B or C virus or a human immunodeficiency

10 virus (HIV) such as HIV-1; chronic persistent hepatitis or chronic active hepatitis; autoimmune liver disease such as autoimmune hepatitis; alcoholic liver disease; fatty liver disease; primary biliary cirrhosis; primary sclerosing cholangitis, biliary atresia; liver disease

15 resulting from medical treatment (drug-induced liver disease); or a congenital liver disease. The methods of the invention can be extremely useful, for example, in alleviating concerns of potential liver damage due to methotrexate treatment. Periodic monitoring of liver

20 fibrosis in individuals treated with methotrexate or other drugs associated with risk of liver damage can be conveniently performed using the non-invasive methods of the invention, without the risks associated with liver biopsy.

25 In one embodiment, the methods of the invention are useful for differentiating individuals having a Metavir score of F0 or F1 from individuals having a Metavir score of F2, F3 or F4. Metavir scoring is a well accepted system for grading liver biopsy specimens and is

30 described in Knodell, *supra*, 1981. F0 is equivalent to the absence of fibrosis; F1 signifies portal fibrosis

without septa. F2 signifies portal fibrosis with a few septa. F3 signifies numerous septa without cirrhosis. F4 signifies cirrhosis.

It is understood that the methods of the invention are useful for diagnosing the presence or severity of fibrosis associated with a variety of fibrotic disorders, including but not limited to liver fibrosis, pulmonary fibrosis, kidney fibrosis, prostate fibrosis and breast fibrosis. The methods of the invention can be applied, without limitation, to diagnosing the presence or severity of pulmonary fibrosis, for example, idiopathic pulmonary fibrosis or emphysema; kidney fibrosis; bladder fibrosis; periureteric fibrosis or retroperitoneal fibrosis; endomyocardial fibrosis, aortic aneurysm disease; rheumatoid diseases such as rheumatoid arthritis or systemic lupus erythematosus; or another fibrotic disorder such as Alzheimer's disease. It is understood that a  $\alpha$ 2-MG/HA/TIMP-1,  $\alpha$ 2-MG/HA/YKL-40 or  $\alpha$ 2-MG/HA/TIMP-1/YKL-40 panel or other combination of markers disclosed herein as useful for diagnosing the presence or severity of liver fibrosis also can be used to diagnose the presence or severity of fibrosis in another disorder.

It is understood that the diagnostic methods of the invention are applicable to a variety of individuals including individuals with chronic or active disease, individuals with one or more symptoms of fibrotic disease, asymptomatic or healthy individuals and individuals at risk for one or more fibrotic diseases. It further is clear to the skilled person that the

methods of the invention can be useful, for example, to corroborate an initial diagnosis of disease or to gauge the progression of fibrosis in an individual with a previous definitive diagnosis of fibrotic disease. The 5 methods of the invention can be used to monitor the status of fibrotic disease over a period of time and further can be used, if desired, to monitor the efficacy of therapeutic treatment. If desired, the results obtained from a sample from an individual undergoing 10 therapy can be compared, for example, to the individual's baseline results prior to treatment, to results earlier during treatment, or to a historic or reference value.

Samples

A variety of samples can be useful in 15 practicing the methods of the invention including, for example, blood, serum, plasma, urine, saliva and liver tissue. In one embodiment, a single sample is obtained from the individual to be diagnosed. Such a sample can be, for example, a serum sample.

20 As used herein, the term "sample" means a biological specimen that contains one or more fibrotic markers such as  $\alpha$ 2-MG, HA or TIMP-1. A sample can be, for example, a fluid sample such as whole blood, plasma, saliva, urine, synovial fluid or other bodily fluid, or a 25 tissue sample such as a lung, liver, kidney, prostate or breast tissue sample. One skilled in the art understands that fluid samples can be diluted, if desired, prior to analysis.

One skilled in the art understands that a single sample can be obtained from the individual to be diagnosed and can be subdivided prior to detecting  $\alpha_2$ -MG-, HA- and TIMP-1. One skilled in the art also 5 understands that, if desired, two or more samples can be obtained from the individual to be diagnosed and that the samples can be of the same or a different type. In one embodiment,  $\alpha_2$ -MG-, HA- and TIMP-1 each are detected in serum samples. In another embodiment, a single serum 10 sample is obtained from an individual and subdivided prior to detecting  $\alpha_2$ -MG-, HA- and TIMP-1.

#### $\alpha_2$ -macroglobulin

The methods of the invention rely, in part, on detecting  $\alpha_2$ -macroglobulin in a sample.  $\alpha_2$ -MG is a 15 conserved, highly abundant component of plasma that functions as a broad spectrum protease-binding protein to clear active proteases from tissue fluids. Unlike active site protease inhibitors, members of the  $\alpha_2$ -macroglobulin family do not inactivate the catalytic activity of their 20 protease substrates but act by physical entrapment of the target protease within the folds of the  $\alpha_2$ -MG family member.  $\alpha_2$ -MG is itself cleaved by target proteases; reorganization of the  $\alpha_2$ -MG molecule results in sequestering of the target protease within an internal 25 pocket of the  $\alpha_2$ -MG molecule (Starkey et al., Biochem. J. 131:823-831 (1973)). While an  $\alpha_2$ -MG entrapped protease is sterically prevented from interacting with macromolecular substrates such as 30 proteins, it remains active against low molecular mass substrates, such as amide and ester compounds, able to diffuse into the  $\alpha_2$ -MG cage to access the enzymatic site.

Thus,  $\alpha$ 2-MG activity is characterized, in part, by the ability to inhibit proteolytic activity but not amidolytic activity of a protease substrate.  $\alpha$ 2-MG also is characterized by the ability to shield entrapped 5 proteases from antibodies and high molecular mass active site inhibitors. For example, trypsin bound by  $\alpha$ 2-MG is protected from inhibition by soybean trypsin inhibitor (STI).

In contrast to the restricted specificity of 10 active-site protease inhibitors,  $\alpha$ 2-MG acts on a broad spectrum of proteases with diverse substrate specificity and catalytic activity. Such target proteases include trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin and papain. Substrate diversity is 15 determined, in part, by the  $\alpha$ 2-MG "bait" region, a highly flexible and solvent-exposed sequence of 30-40 residues that contains at least one site sensitive to cleavage by each of the major classes of proteolytic enzyme.

As used herein, the term " $\alpha$ 2-macroglobulin" is 20 synonymous with " $\alpha$ 2-MG" and means a protein with significant structural homology to human  $\alpha$ 2-MG (SEQ ID NO: 2) and having broad spectrum protease inhibitory activity.  $\alpha$ 2-MG contains a unique thiol ester bond that is inactivated by small primary amines such as 25 methylamine. Thus,  $\alpha$ 2-MG activity can be characterized, in part, by methylamine-sensitive protease inhibitory activity.  $\alpha$ 2-MG can be distinguished, if desired, from other members of the  $\alpha$ 2-macroglobulin family such as related protease-binding proteins and C3, C4 and C5 of 30 the complement system (Sottrup-Jensen, " $\alpha$ 2-Macroglobulin and Related Thiol Ester Plasma Proteins," in Putnam

(Ed.), The Plasma Proteins: Structure, Function and Genetic Control Second edition, Orlando: Academic Press (1987), pages 191-291. It is understood that an assay for detecting  $\alpha_2$ -MG can be specific for  $\alpha_2$ -MG or can 5 additionally detect one or more other members of the  $\alpha_2$ -macroglobulin family.

The methods of the invention rely, in part, on detecting  $\alpha_2$ -macroglobulin in a sample. As used herein, the phrase "detecting  $\alpha_2$ -MG" means any quantitative or 10 qualitative assay for determining the presence of  $\alpha_2$ -MG. As used herein, the phrase "determining the level of  $\alpha_2$ -MG" means any direct or indirect quantitative assay for  $\alpha_2$ -MG.

Similarly, detecting any specified fibrotic 15 marker in a sample means determining whether the marker is present in the sample, said fibrotic marker having a positive or negative correlation with liver fibrosis or with another fibrotic disorder such as are described herein above. It is understood that detection can refer 20 to non-quantitative analysis, for example, the presence or absence of a particular trait, variable or biochemical or serological substance.

Diagnosis is based on analyzing the sample for the presence or level of the fibrotic marker or other 25 characteristic and comparing it to a reference value, where the reference value serves to assist in differentiating those with a fibrotic disorder from other individuals. Where the fibrotic marker is a biochemical or serological marker, determining a "level" in a sample 30 means quantifying the fibrotic marker by determining, for

example, the relative or absolute amount of RNA, protein or activity of the fibrotic marker. Thus, determining a level in a sample encompasses, without limitation, analysis of relative and absolute RNA, protein and 5 activity levels as well as other direct and indirect measurements of the fibrotic marker as discussed further below. It is understood that any assay useful for determining a "level" of a fibrotic marker also is useful for "detecting" the marker.

10 A variety of assays for detecting  $\alpha$ 2-MG are known in the art and include direct and indirect assays for  $\alpha$ 2-MG RNA,  $\alpha$ 2-MG protein and  $\alpha$ 2-MG activity.  $\alpha$ 2-MG can be detected, or an  $\alpha$ 2-MG level can be determined, for example, by analysis of  $\alpha$ 2-MG mRNA levels using routine 15 techniques such as Northern analysis or RT-PCR, or other methods based on hybridization to a nucleic acid sequence that is complementary to a portion of the  $\alpha$ 2-MG coding sequence. For example, conditions and probes for Northern analysis and RNA slot blot hybridization of 20  $\alpha$ 2-MG RNA in human samples are described in Ortego et al., Exp. Eye Res. 65:289-299 (1997), and Simon et al., Cancer Res. 56:3112-3117 (1996), respectively.

$\alpha$ 2-MG also can be detected, or an  $\alpha$ 2-MG level can be determined, by assaying for  $\alpha$ 2-MG protein by a 25 variety of methods. Immunoassays, including radioimmunoassays, enzyme-linked immunoassays and two-antibody sandwich assays as described further below, are useful in the methods of the invention. For example, in nephelometry assays, complexes of  $\alpha$ 2-MG and anti- $\alpha$ 2-MG 30 antibody result in increased light scatter that is converted to a peak rate signal, which is a function of

- the sample  $\alpha_2$ -MG concentration.  $\alpha_2$ -MG also can be detected, for example, by laser immunonephelometry using a Behring Nephelometer Analyzer (Fink et al., J. Clin. Chem. Clin. Biol. Chem. 27:261-276 (1989)) and rabbit
- 5 anti-human  $\alpha_2$ -MG antiserum as described in Naveau et al., Dig. Diseases Sci. 39:2426-2432 (1994), or using the nephelometry assay commercially available from Beckman Coulter (Brea, CA; kit #449430). Furthermore, monoclonal and polyclonal anti- $\alpha_2$ -MG antibodies useful in
- 10 immunoassays can be readily obtained from a variety of sources. As examples, affinity purified goat anti-human  $\alpha_2$ -MG and peroxidase-labeled goat anti-human  $\alpha_2$ -MG antibodies suitable for immunoassays such as ELISA assays and western blotting are available from Cedarlane
- 15 Laboratories Limited (Ontario, Canada; CL20010AP and CL20010APH) and Affinity Biologicals Incorporated (Ontario, Canada; GAA2M-AP and GAA2M-APH). Levels of  $\alpha_2$ -MG protein also can be determined by quantifying the amount of purified  $\alpha_2$ -MG protein.
- 20 Purification of  $\alpha_2$ -macroglobulin can be achieved, for example, by HPLC, alone or in combination with mass spectrophotometry, or as described, for example, in Hall and Roberts, Biochem. J. 171:27-38 (1978) or Imber and Pezzo, J. Biol. Chem. 256:8134-8139 (1981)).
- 25 Quantitation can be determined by well known methods including Bradford assays, Coomassie blue staining and assays for radiolabeled protein.

A variety of assays for  $\alpha_2$ -MG activity also can be useful for detecting  $\alpha_2$ -MG or determining a level of

30  $\alpha_2$ -MG in a sample according to a method of the invention.  $\alpha_2$ -MG can be detected or a level of  $\alpha_2$ -MG can be determined indirectly, for example, as a function of

inhibition of target protease activity, without a corresponding inhibition of amidolytic activity. As discussed above,  $\alpha$ 2-MG-bound proteases retain the ability to hydrolyze amide and ester bonds of small substrates,

5 even while high molecular mass substrates such as proteins cannot be hydrolyzed (see, for example, Armstrong et al., Develop. Compar. Immunol. 23:375-390 (1999)). As an example,  $\alpha$ 2-MG can be detected or the level of  $\alpha$ 2-MG can be determined by assaying for

10 inhibition of trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin, or papain activity without inhibition of amidolytic activity. Convenient substrates to be analyzed include  $^{14}$ C-labeled casein and  $^{125}$ I-fibrin.

The characteristic of broad protease substrate specificity distinguishes  $\alpha$ 2-MG from inhibitors of protease active sites. Based on this characteristic,  $\alpha$ 2-MG can be detected or the level of  $\alpha$ 2-MG can be determined by assaying for inhibition of the activity of two or more proteases with different active site specificities.  $\alpha$ 2-MG can be detected or the level of  $\alpha$ 2-MG in a sample can be determined, for example, by analyzing the reduction in protease activity of two or more target proteases such as two or more of the following proteases: trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin and papain. Labeled protease substrates such as  $^{14}$ C-casein or  $^{125}$ I-fibrin can be useful in such methods (Armstrong et al., *supra*, 1999).

$\alpha$ 2-MG also can be detected or the level of  $\alpha$ 2-MG determined based on the ability of  $\alpha$ 2-MG to shield a bound protease from an antibody or a high molecular

weight inhibitor. A target protease such as trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin, or papain can be added to a plasma sample. Following removal of unbound protease, for example, by

- 5 immunoprecipitation with anti-protease antibody, the amount of protease bound by  $\alpha$ 2-MG can be determined using a low molecular mass amide or ester substrate. The amount of hydrolyzed low molecular mass substrate is an indicator of the amount of protected,  $\alpha$ 2-MG-bound,
- 10 protease and, therefore, of the concentration of  $\alpha$ 2-MG. Similarly, a sample can be reacted first with a protease such as trypsin and subsequently with excess protease inhibitor such as soybean trypsin inhibitor before assaying residual trypsin activity with a low molecular
- 15 mass substrate, such as the amide BApNA ( $N^{\alpha}$ -benzoyl-DL-arginine *p*-nitroanilide (Ganrot, Clin. Chem. Acta 14:493-501 (1966); Armstrong et al., J. Exp. Zool. 236:1-9 (1985)). Trypsin not sequestered by  $\alpha$ 2-MG is inactivated by the trypsin inhibitor, with only  $\alpha$ 2-MG-protected
- 20 trypsin remaining capable of substrate hydrolysis. Thus, a positive reaction in a soybean trypsin inhibitor assay detects  $\alpha$ 2-MG and is a quantitative measure of the amount of  $\alpha$ 2-MG (Armstrong et al., *supra*, 1999). One skilled in the art understands that the presence of low molecular
- 25 mass protease inhibitors capable of inactivating  $\alpha$ 2-MG-bound enzyme can affect the results obtained with such an assay. It is further understood that these and other routine assays for  $\alpha$ 2-MG activity, as well as  $\alpha$ 2-MG RNA or protein levels, can be useful for detecting  $\alpha$ 2-MG
- 30 or determining a level of  $\alpha$ 2-MG in a method of the invention.

Hyaluronic acid

The methods of the invention further rely, in part, on detecting hyaluronic acid or determining a level of hyaluronic acid in a sample. Hyaluronic acid, also known as hyaluronate or hyaluronan, is a high molecular weight polysaccharide with an unbranched backbone made up of alternating glucuronic acid and  $\beta(1,3)$ -N-acetylglucosamine moieties linked by  $\beta$ -1,4 linkages. Hyaluronic acid can have a length of a few to more than 1,000 dimeric units, with each dimeric unit having a molecular weight of about 450 D. Hyaluronic acid, which is produced principally by fibroblasts and other specialized connective tissue cells, plays a structural role in the connective tissue matrix. Furthermore, hyaluronic acid is widely distributed throughout the body and can be found as a free molecule in, for example, plasma, synovial fluid and urine. In plasma, hyaluronic acid has a relatively short half-life.

Serum HA levels can be elevated in liver diseases including cirrhosis (Bramley et al., J. Hepatol. 20 13:8-13 (1991); Ueno et al., Gastroenterol. 105:475-481 (1993); Oberti et al., Gastroenterol. 113:1609-1616 (1997); and McHutchison et al., J. Gastroenterol. Hepatol. 15:945-951 (2000)). Serum HA levels also can be elevated during synovial inflammation and cartilage destruction seen in rheumatoid arthritis; these levels have been found to correlate with disease activity and degree of synovial involvement (Konttinen et al., Clin. Chimica Acta 193:39-48 (1990); Poole et al., Arthritis Rheum. 37:1030-1038 (1994); Goldberg et al., Arthritis Rheum. 34: 799-807 (1991); and Emile et al., J. 25 30

Rheum. 23:974-978 (1996)). Elevated serum levels of HA also can be present, for example, in patients with osteoarthritis (OA), progressive systemic sclerosis (PSS) and systemic lupus erythematosus (SLE).

5 As used herein, the term "hyaluronic acid" is synonymous with "HA" and means a polymer of two or more dimeric units of alternating glucuronic acid and  
β(1,3)-N-acetylglucosamine moieties linked by β-1,4 linkages. As used herein, the phrase "detecting HA"  
10 means any quantitative or qualitative assay for determining the presence of HA, and the phrase "determining the level of HA" means any direct or indirect quantitative assay for HA. In view of the above, it is understood that the phrase "detecting HA"  
15 encompasses "determining the level of HA."

HA can be detected or a level of HA can be determined using one of a variety of well known assays based on HA-binding proteins or anti-HA antibodies, or by quantitation of purified HA. HA-binding proteins, for  
20 example, can be useful in detecting HA; a radiometric assay for HA based on <sup>125</sup>I-labelled HA-binding protein is available from Pharmacia (Guechot et al., Clin. Chem. 42:558-563 (1996). Other commercial assays based on HA-binding proteins are available, for example, from  
25 Corgenix (Westminster, CO; kit 029001). In addition, HA can be detected or a level of HA can be determined using hyaluronectin as described in Maingonnat and Delpech, Ann. Clin. Biochem. 28:305-306 (1991), or using the kit available from Nalgenunc International (Rochester, NY;  
30 Delpech and Bertrand, Anal. Biochem. 149:555-565 (1985)). Assays for detecting HA or determining a level of HA

include a variety of competitive and non-competitive binding assays, for example, competitive binding assays using  $^{125}\text{I}$ -labeled HA binding protein; competitive binding assays based on alkaline phosphatase

5 labeled-hyaluronectin (HN); and non-competitive binding assays based on peroxidase-labeled proteoglycan or peroxidase-labeled HA-binding protein, among others (Lindquist et al., Clin. Chem. 38:127-132 (1992)). See, also, Delpech and Bertrand, *supra*, 1985; Engstrom-Laurent

10 et al., Scand. J. Clin. Lab. Invest. 45:497-504 (1985); Brandt et al., Acta Otolaryn. 442 (Suppl.):31-35 (1987); Goldberg, Anal. Biochem. 174:448-458 (1988); Chichibu et al., Clin. Chim. Acta 181:317-324 (1989); Li et al., Conn. Tissue Res. 19:243-254 (1989); Poole et al., Arth.

15 Rheum. 33:790-799 (1990); Poole et al., J. Biol. Chem. 260:6020-6025 (1985); and Laurent and Tengblad, Anal. Biochem. 109:386-394 (1980)). Assays for detecting HA or determining a level of HA in a sample can be performed using a variety of immunoassay formats,

20 including radioimmunoassays and enzyme-linked immunoassays. Anti-HA antiserum useful in immunoassays can be, for example, affinity purified sheep anti-HA antiserum available from Biotrend (Cologne, Germany; #5029-9990).

25 A level of HA also can be determined by purifying HA from a sample, and quantifying the amount of purified polysaccharide. High performance liquid chromatography can be used alone or in conjunction with mass spectrophotometry. As an example, HPLC can be used

30 to determine HA levels after digestion of samples containing an internal standard with hyaluronidase, separation by a reversed phase octadecylsilyl column and

elution with 0.01 M tetrabutylammonium phosphate-acetonitrile (83:17, v/v) at pH 7.35 (Payan et al., J. Chromatogr. 566:9-18 (1991)).

HA levels have been shown to correlate with 5 hyaluronidase levels (Bray et al., Am. Rev. Respir. Dis. 3:284-288 (1991)). Thus, HA can be detected or a level of HA can be determined indirectly by assaying for hyaluronidase activity. Assays for hyaluronidase activity are known in the art, as described in Bray et 10 al., *supra*, 1991. One skilled in the art understands that these and other routine assays for determining hyaluronidase or HA levels are encompassed by the phrases "detecting HA" and "determining the level of HA" and can be useful in diagnosing the presence or severity of liver 15 fibrosis according to a method of the invention.

#### TIMP-1

The methods of the invention also are based on detecting TIMP-1 in a sample and, in particular embodiments, on determining a level of TIMP-1 in a 20 sample. Tissue inhibitors of metalloproteinases (TIMPs) regulate the activity of the matrix metalloproteinases (MMPs), which are an important group of ECM-degradative enzymes that include gelatinase A (MMP-2) and gelatinase B (MMP-9). In normal liver, matrix components such as 25 collagens, fibronectin, laminin, tenascin, undulin and entactin are constantly remodeled by matrix degrading enzymes to control deposition of extracellular matrix. Elevation of TIMP levels results in inhibition of MMP activity and favors the accumulation of extracellular 30 matrix. The TIMPs, which include TIMP-1, TIMP-2, TIMP-3

and TIMP-4, interact with the matrix metalloproteinases with a 1:1 stoichiometry and inhibit metalloprotease activity through reversible non-covalent binding. TIMP-1, TIMP-2 and TIMP-3 have similar MMP-inhibitory activities, inhibiting the proteolytic activity of collagenase, gelatinase, stromelysin, proteoglycanase and metalloelastases although their localization and regulation differ (Cawston et al., "Protein Inhibitors of Metalloproteinases" in Barrett and Salvesen (Eds), 10 Proteinase Inhibitors Amsterdam Elsevier pages 589-610 (1986)).

Human TIMP-1 is a 184 amino acid sialoglycoprotein with a molecular weight of 28.5 kDa (Murphy et al., Biochem. J. 195:167-170 (1981); Dockerty et al., Nature 318:66-69 (1985); and Bodden et al., J. Biol. Chem. 269:18943-18952 (1994)). TIMP-1 inhibits all active metalloproteinases, for example, interstitial collagenase MMP-1 as well as stromelysin and gelatinase B (MMP-9). The nucleic acid sequence (SEQ ID NO: 3) and 20 corresponding amino acid sequence (SEQ ID NO: 4) of human TIMP-1 are shown in Figure 2.

As used herein, the term "tissue inhibitor of metalloproteinase-1" is synonymous with "TIMP-1" and means a protein with significant structural homology to 25 human TIMP-1 (SEQ ID NO: 4) that inhibits the proteolytic activity of metalloproteinases with a specificity similar to human TIMP-1. The presence of human TIMP-1 can be conveniently detected by the presence of epitopes reactive with a known specific anti-TIMP-1 antibody such 30 as 7-6C1 or 7-23G9.

As used herein, the phrase "detecting TIMP-1" means any quantitative or qualitative assay for determining the presence of TIMP-1, and the phrase "determining the level of TIMP-1" means any direct or 5 indirect quantitative assay for TIMP-1. In view of the above, it is understood that the phrase "detecting TIMP-1" encompasses "determining the level of TIMP-1."

Assays for detecting TIMP-1 and for determining a level of TIMP-1 include well known assays for TIMP-1 10 RNA, protein and enzymatic activity. Methods of determining TIMP-1 RNA levels by Northern analysis or RT-PCR are well known in the art (Yoshiji et al., Int. J. Cancer 69:131-134 (1996); Janowska-Wieczorek et al., Exp. Hematol. 28:1274-1285 (2000); and Groft et al., Br. J. Cancer 85:55-63 (2001)) as described further below. 15 TIMP-1 protein can be detected or the level of TIMP-1 protein can be conveniently determined, for example, by radioimmunoassay as described in Brophy et al., Biochem. Biophys. Res. Comm. 167:898-903 (1990) or by two-antibody 20 sandwich assay as described in Murawaki et al., Clinica Chimica Acta 218:47-58 (1993). Plasma concentrations of TIMP-1 protein can be assayed by ELISA with a kit commercially available from Amersham Pharmacia (see, also Example III). Levels of TIMP-1 protein also can be 25 determined by quantifying the amount of purified TIMP-1 protein. Purification of TIMP-1 can be achieved, for example, by HPLC, alone or in combination with mass spectrophotometry, or as described, for example, in Murphy et al., Biochem. J. 195:167-170 (1981), or 30 Stricklin and Welgus, J. Biol. Chem. 258:12252-12258 (1983). TIMP-1 also can be detected or a level of TIMP-1 determined by assaying for inhibition of the activity of

one or more metalloproteases, for example, using reverse gelatin zymography as described in Kossakowska et al., Amer. J. Pathology 153:1895-1902 (1998). Assays for TIMP-1 RNA, protein or activity are described further 5 hereinbelow, and one skilled in the art understands that these and other routine assays for detecting TIMP-1 are encompassed by the methods of the invention.

Rule-in/Rule-out analysis

As disclosed herein, two sets of cut-off values 10 can be used to increase the accuracy of an assay based on the  $\alpha$ 2-MG/HA/TIMP-1 three-marker panel. As set forth in Example II, a first set of cut-off values for  $\alpha$ 2-MG, HA and TIMP-1 were selected based on optimization for sensitivity in order to first rule out fibrosis, followed 15 by analysis of the "positive" population using a second set of cut-off values optimized for specificity to determine the presence of significant fibrosis. Table 7 shows the results of the dual optimization strategy on the 194 HCV patient study population. The primary 20 cut-offs were set at 2.0 mg/ml, 35 ng/ml and 1000 ng/ml for  $\alpha$ 2-MG, HA and TIMP-1, respectively, to achieve a high sensitivity in the primary analysis. Any samples with all three of  $\alpha$ 2-MG, HA and TIMP-1 levels above the primary cut-off values were tentatively indicated to be 25 positive for F2-F4 fibrosis and were further evaluated using a second set of cut-off values of 2.0 mg/ml, 60 ng/ml and 1575 ng/ml for  $\alpha$ 2-MG, HA and TIMP-1, respectively, which were obtained by optimizing for specificity.

Using the second set of cut-off values optimized for high specificity, 54 of the 122 patients initially designated as positive for F2-F4 fibrosis were confirmed positive, only one of which was a false positive. In sum, of the 194 HCV patients in the study population, 72 were classified as negative (having F0-F1 fibrosis) and 54 were classified as positive (having F2-F4 fibrosis), while 68 samples had indeterminate results and were not classified. When the indeterminate samples were excluded, the  $\alpha$ 2-MG/HA/TIMP-1 panel performed with a positive predictive value of about 98% and a negative predictive value of about 79%. Furthermore, in a more typical patient population having 30% fibrosis prevalence, the same panel performs with positive and negative predictive values of close to 93%. These results indicate that the use of primary and secondary cut-off levels, whereby sensitivity is initially optimized followed by optimization for specificity, can increase the overall accuracy of a three-marker test, resulting in a panel test with about 93% accuracy for non-indeterminate samples, which make up about 70% of the samples tested.

Thus, the present invention provides a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual by determining an  $\alpha$ 2-MG level in a sample from the individual; determining a HA level in a sample from the individual; determining a TIMP-1 level in a sample from the individual; and diagnosing the individual as having no or mild liver fibrosis when the  $\alpha$ 2-MG level is below an  $\alpha$ 2-MG cut-off value X<sub>1</sub>, the HA level is below a HA cut-off value Y<sub>1</sub> or the TIMP-1 level is below a TIMP-1

cut-off value Z1; diagnosing the individual as having moderate to severe liver fibrosis when the  $\alpha$ 2-MG level is above an  $\alpha$ 2-MG cut-off value X2, the HA level is above a HA cut-off value Y2 and the TIMP-1 level is above a 5 TIMP-1 cut-off value Z2; and diagnosing remaining individuals as having an indeterminate status.

The methods of the invention based on dual cut-off values for the levels of the  $\alpha$ 2-MG, HA and TIMP-1 markers can be useful in differentiating no or mild liver 10 fibrosis from moderate to severe liver fibrosis in a variety of patient populations. Such methods can be useful, for example, in diagnosing an individual having a liver disease such as viral hepatitis, autoimmune liver disease such as autoimmune hepatitis, alcoholic liver 15 disease, fatty liver disease or drug-induced liver disease. In one embodiment, a method of the invention is used to differentiate no or mild liver fibrosis from moderate to severe liver fibrosis in an individual infected with hepatitis C virus. Samples useful in a 20 method of the invention based on dual cut-off values include, but are not limited to, blood, serum, plasma, urine, saliva and liver tissue. In one embodiment, a method of the invention is practiced by determining the  $\alpha$ 2-MG level, HA level and TIMP-1 level in one or more 25 serum samples.

In a further embodiment, the present invention provides a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual, where the differentiation is based on an X1 30 cut-off value between 1.8 and 2.2 mg/ml; a Y1 cut-off value between 31 and 39 ng/ml; a Z1 cut-off value between

900 and 1100 ng/ml; an X2 cut-off value between 1.8 and 2.2 mg/ml; a Y2 cut-off value between 54 and 66 ng/ml; and a Z2 cut-off value between 1415 and 1735 ng/ml. In another embodiment, the differentiation is based on an X1 5 cut-off value of 2.0 mg/ml; a Y1 cut-off value of 35 ng/ml; a Z1 cut-off value of 1000 ng/ml; an X2 cut-off value of 2.0 mg/ml; a Y2 cut-off value of 60 ng/ml; and a Z2 cut-off value of 1575 ng/ml. In yet another embodiment, the differentiation is based on an X1 cut-off 10 value of 2.0 mg/ml; a Y1 cut-off value of 37 ng/ml; a Z1 cut-off value of 1100 ng/ml; an X2 cut-off value of 2.0 mg/ml; a Y2 cut-off value of 60 ng/ml; and a Z2 cut-off value of 1575 ng/ml. In a further embodiment, X1, Y1, Z1, X2, Y2 and Z2 are selected such that, in a population 15 having up to 30% liver fibrosis prevalence, at least 65% of individuals in the population are diagnosed as having no or mild fibrosis or moderate to severe fibrosis with an accuracy of at least 90%. In yet a further embodiment, X1, Y1, Z1, X2, Y2 and Z2 are selected such 20 that, in a population having up to 10% liver fibrosis prevalence, at least 70% of individuals in the population are diagnosed as having no or mild fibrosis or moderate to severe fibrosis with an accuracy of at least 90%.

As set forth above, the methods of the 25 invention are highly accurate for determining the presence or severity of fibrosis in a subgroup of the entire patient population assayed. For example, as shown in Table 7, the methods of the invention perform with more than 93% accuracy in determining the F0-F1 or F2-F4 30 fibrosis status in about 70% of a patient population having a liver fibrosis prevalence of 30%. The remaining 30% of the patient population are indicated to have an

indeterminate status. As used herein, the term "indeterminate status" means that the individual cannot be confidently diagnosed with sufficient predictive value.

5           As used herein, the term "X1" or "X2" refers to an  $\alpha_2$ -MG cut-off value, against which an experimental  $\alpha_2$ -MG sample level is compared. Similarly, as used herein, the term "Y1" or "Y2" refers to an HA cut-off value, against which an experimental HA level is  
10 compared. The term "Z1" or "Z2," as used herein, refers to a TIMP-1 cut-off value against which an experimental TIMP-1 level is compared. X1, Y1 and Z1 cut-offs are combined to determine the presence or severity of fibrosis in a sample. Similarly, X2, Y2 and Z1 cut-off  
15 values are combined to determine the presence or severity of fibrosis in a sample. A sample having an  $\alpha_2$ -MG level less than X1, an HA level less than Y1, or a TIMP-1 level less than Z1 is classified as having F0-F1 fibrosis. A sample having an  $\alpha_2$ -MG level above X1, an HA level above  
20 Y1, and a TIMP-1 level above Z1 is possibly positive for F2-F4 fibrosis and warrants further analysis.  
Furthermore, a sample having an  $\alpha_2$ -MG level above X2, an HA level above Y2, and a TIMP-1 level above Z2 is  
25 classified as having F2-F4 fibrosis. A sample having an  $\alpha_2$ -MG level above X1, an HA level above Y1, and a TIMP-1 level above Z1 but one or more levels below X2, Y2 or Z2 is classified as having an "indeterminate status." It is understood that X2 generally is equal to or greater than X1; Y2 generally is equal to or greater than Y1; and Z2  
30 generally is equal to or greater than Z1.

One skilled in the art can select  $\alpha$ 2-MG, HA and TIMP-1 cut-offs X<sub>1</sub>, Y<sub>1</sub>, Z<sub>1</sub>, X<sub>2</sub>, Y<sub>2</sub> and Z<sub>2</sub> to achieve one or more clinically useful parameters, such as a desired sensitivity or specificity, or a desired negative predictive value, positive predictive value or accuracy for a patient population having a particular fibrosis prevalence. Factorial Design Optimization, also known as Design of Experiments, methodology can be used, for example, to select the appropriate cut-off values. As disclosed herein in Example II, optimization software (DOE Keep It Simple Statistically from Air Academy Associates (Colorado Springs, CO) was used in a central composite design experiment to simultaneously vary the three cut-offs X<sub>1</sub>, Y<sub>1</sub> and Z<sub>1</sub>, and then to simultaneously vary the three cut-offs X<sub>2</sub>, Y<sub>2</sub> and Z<sub>2</sub>. In particular, the  $\alpha$ 2-MG cut-off was varied from 2.0 to 5.0 mg/ml; the HA cut-off was varied from 25-75 ng/ml; and the TIMP-1 cut-off was varied from 1000-1700 ng/ml. By comparing the test results determined for the 194 patients in the database (see Table 4) with the assigned X<sub>1</sub>, Y<sub>1</sub> and Z<sub>1</sub> cut-offs, each of the 194 samples were determined to be a true positive, true negative, false positive or false negative, and the clinical parameters of sensitivity, specificity, negative predictive value, positive predictive value and accuracy were determined for the study patient population. Although determination of the  $\alpha$ 2-MG, HA and TIMP-1 cut-off values is illustrated herein using the DOE KISS program, one skilled in the art understands that other computer programs for identifying cooperative interactions among multiple variables and for performing simultaneous equation calculations also can be used. For example, ECHIP optimization software, available from ECHIP, Incorporated (Hockessin, DE), or

Statgraphics optimization software, available from STSC, Incorporated (Rockville, MD), also can be useful in determining  $\alpha_2$ -MG, HA and TIMP-1 cut-off values useful in the methods of the invention.

5           The clinical parameters of sensitivity, specificity, negative predictive value, positive predictive value and accuracy are calculated using true positives, false positives, true negatives and false negatives. A "true positive" sample is a sample positive  
10 for the indicated stage of fibrosis according to clinical biopsy, which is also diagnosed positive according to a method of the invention. A "false positive" sample is a sample negative for the indicated stage of fibrosis by biopsy, which is diagnosed positive according to a method  
15 of the invention. Similarly, a "false negative" is a sample positive for the indicated stage of fibrosis by biopsy, which is diagnosed negative according to a method of the invention. A "true negative" is a sample negative for the indicated stage of fibrosis by biopsy, and also  
20 negative for fibrosis according to a method of the invention. See, for example, Motulsky (Ed.), Intuitive Biostatistics New York: Oxford University Press (1995).

As used herein, the term "sensitivity" means the probability that a diagnostic method of the invention  
25 gives a positive result when the sample is positive, for example, fibrotic with a Metavir score of F2-F4. Sensitivity is calculated as the number of true positive results divided by the sum of the true positives and false negatives. Sensitivity essentially is a measure of  
30 how well a method correctly identifies those with fibrotic disease. In a method of the invention, the  $X_1$ ,

Y<sub>1</sub>, Z<sub>1</sub>, X<sub>2</sub>, Y<sub>2</sub> and Z<sub>2</sub> values can be selected such that the sensitivity of diagnosing an individual is at least about 70%, and can be, for example, at least 75%, 80%, 85%, 90% or 95% in at least 60% of the patient population assayed, or in at least 65%, 70%, 75% or 80% of the patient population assayed.

As used herein, the term "specificity" means the probability that a diagnostic method of the invention gives a negative result when the sample is not positive, for example, not of Metavir fibrosis stage F2-F4. Specificity is calculated as the number of true negative results divided by the sum of the true negatives and false positives. Specificity essentially is a measure of how well a method excludes those who do not have fibrosis. In a method of the invention, the cut-off values X<sub>1</sub>, Y<sub>1</sub>, Z<sub>1</sub>, X<sub>2</sub>, Y<sub>2</sub> and Z<sub>2</sub> can be selected such that, when the sensitivity is at least about 70%, the specificity of diagnosing an individual is in the range of 70-100%, for example, at least 75%, 80%, 85%, 90% or 95% in at least 60% of the patient population assayed, or in at least 65%, 70%, 75% or 80% of the patient population assayed. As illustrated in Example II, a specificity of greater than 98% and a sensitivity of about 77% were achieved in the non-indeterminate patient population, which was about 70% of the patient population having a fibrosis prevalence of 30%.

The term "negative predictive value," as used herein, is synonymous with "NPV" and means the probability that an individual diagnosed as not having fibrosis actually does not have the disease. Negative predictive value can be calculated as the number of true

negatives divided by the sum of the true negatives and false negatives. Negative predictive value is determined by the characteristics of the diagnostic method as well as the prevalence of fibrosis in the population analyzed.

- 5 In a method of the invention, the  $\alpha$ 2-MG, HA and TIMP-1 cut-off values can be selected such that the negative predictive value in a population having a liver fibrosis prevalence of up to 10% is in the range of 75-99% and can be, for example, at least 80%, at least 85%, at
- 10 least 90%, or at least 95%, in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed. The  $\alpha$ 2-MG, HA and TIMP-1 cut-off values also can be selected such that the negative predictive value in a population
- 15 having a liver fibrosis prevalence of up to 20% is in the range of 75-99% and can be, for example, at least 80%, at least 85%, at least 90%, or at least 95%, in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population
- 20 assayed. In addition,  $\alpha$ 2-MG, HA and TIMP-1 cut-off values can be selected such that the negative predictive value in a population having a liver fibrosis prevalence of up to 30% is in the range of 75-99% and can be, for example, at least 80%, at least 85%, at least 90%, or at
- 25 least 95%, in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed.

The term "positive predictive value," as used herein, is synonymous with "PPV" and means the probability that an individual diagnosed as having fibrosis actually has the condition. Positive predictive value can be calculated as the number of true positives

divided by the sum of the true positives and false positives. Positive predictive value is determined by the characteristics of the diagnostic method as well as the prevalence of fibrosis in the population analyzed.

- 5 In a method of the invention, the  $\alpha$ 2-MG, HA and TIMP-1 cut-off values can be selected such that, in a patient population having up to 10% liver fibrosis prevalence, the positive predictive value of the method is at least about 75%, and can be at least 80%, at least 85%, at  
10 least 90% or at least 95% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed. The  $\alpha$ 2-MG, HA and TIMP-1 cut-off values also can be selected such that, in a patient population having up to 20% liver  
15 fibrosis prevalence, the positive predictive value of the method is at least about 75%, and can be at least 80%, at least 85%, at least 90% or at least 95% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population  
20 assayed. Similarly, the  $\alpha$ 2-MG, HA and TIMP-1 cut-off values can be selected such that, in a patient population having up to 30% liver fibrosis prevalence, the positive predictive value of the method is at least about 75%, and can be at least 80%, at least 85%, at least 90% or at  
25 least 95% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed.

Predictive values, including negative and positive predictive values, are influenced by the  
30 prevalence of the disease in the population analyzed. In the methods of the invention, the cut-off values X<sub>1</sub>, Y<sub>1</sub>, Z<sub>1</sub>, X<sub>2</sub>, Y<sub>2</sub> and Z<sub>2</sub> can be selected to produce a desired

clinical parameter for a clinical population with a particular liver fibrosis prevalence. For example, cut-off values can be selected for a liver fibrosis prevalence of up to 10%, 12%, 15%, 18%, 20%, 25% or 30% 5 which can be seen, for example, in a hepatologist's office. Cut-off values also can be selected for a liver fibrosis prevalence of up to 1%, 2%, 3%, 4%, 5%, 6%, 7% or 8%, which can be representative of the fibrosis prevalence seen in a general practitioner's office.

10 As used herein, the term "accuracy" means the overall agreement between the diagnostic method and the disease state. Accuracy is calculated as the sum of the true positives and true negatives divided by the total number of sample results and is affected by the 15 prevalence of fibrosis in the population analyzed. The  $\alpha$ 2-MG, HA and TIMP-1 cut-off values can be selected such that the accuracy of a method of the invention in a patient population having a liver fibrosis prevalence of up to 10% is at least about 80% and can be, for example, 20 at least 85%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, or 99% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed. The  $\alpha$ 2-MG, HA and TIMP-1 25 cut-off values also can be selected such that the accuracy of a method of the invention in a patient population having a liver fibrosis prevalence of up to 20% is at least about 80% and can be, for example, at least 85%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, or 99% in at least 60% of the patient population assayed, 30 for example, in at least 65%, 70%, 75% or 80% of the patient population assayed. Similarly, the  $\alpha$ 2-MG, HA and TIMP-1 cut-off values can be selected such that the

accuracy of a method of the invention in a patient population having a liver fibrosis prevalence of up to 30% is at least about 80% and can be, for example, at least 85%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%,  
5 or 99% in at least 60% of the patient population assayed, for example, in at least 65%, 70%, 75% or 80% of the patient population assayed.

Methods not limited to specific markers

The present invention also provides a method of  
10 diagnosing the presence or severity of liver fibrosis in an individual by comparing a level of a first fibrotic marker X in the individual to a cut-off value X<sub>1</sub> to determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second  
15 fibrotic marker Y in the individual to a cut-off value Y<sub>1</sub> to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where, in a  
20 population with up to 40% fibrosis prevalence, at least 65% of individuals in the population are diagnosed with an accuracy of at least 90%.

A method of the invention can include, if desired, comparing a level of a third fibrotic marker Z  
25 in the individual to a cut-off value Z<sub>1</sub> to determine whether the individual is positive for the third fibrotic marker Z and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z. In one embodiment, the first

fibrotic marker is  $\alpha$ 2-MG, the second fibrotic marker is HA, and the third fibrotic marker is TIMP-1.

In another embodiment, the levels of at least three fibrotic markers are compared, and, in a further 5 embodiment, the levels of exactly three fibrotic markers are compared to their respective cut-off values. In additional embodiments, the levels of at least four or at least five fibrotic markers are compared. A method of the invention can be useful, for example, to 10 differentiate no or mild liver fibrosis from moderate to severe liver fibrosis.

In a specific embodiment, a method of the invention serves to diagnose at least 65% of individuals in a population with up to 30% fibrosis prevalence with 15 an accuracy of at least 93%. In a further embodiment, a method of the invention serves to diagnose at least 70% of individuals in a population with up to 20% fibrosis prevalence with an accuracy of at least 94%. In yet a further embodiment, a method of the invention serves to 20 diagnose at least 70% of individuals in a population with up to 10% fibrosis prevalence with an accuracy of at least 96%.

The methods of the invention provide unparalleled performance in diagnosing the presence or 25 severity of liver fibrosis. While not all patients are provided with a diagnosis, the majority are diagnosed with extremely good accuracy. As an example, in a patient population with about 40% fibrosis prevalence, almost 70% of the population are diagnosed with more than 30 91% accuracy and with a positive predictive value of more

than 96% and a negative predictive value of more than 89%. This excellent performance contrasts with alternative methods such as the method of Pynnard et al., Lancet 357:1069 (2001). Using the method of Pynnard et al. based on analysis of the six markers  $\alpha$ 2-MG,  $\alpha$ 2-globulin, total bilirubin,  $\gamma$ -globulin, apoA1 and GGT, only about 50% of a population having about 40% fibrosis prevalence are diagnosed, and only with an accuracy of about 89% (see Table 8). Thus, the methods of the invention provide an improvement, in that a significantly greater percentage of a patient population (about 70% as compared to about 50%) are diagnosed, and with an accuracy of more than 91% as compared to an accuracy of around 89% (see Table 8). Due to the novel performance characteristics of a method of the invention, biopsy is typically unnecessary in at least 65% of a patient population, and the patients diagnosed can have confidence in a diagnosis that is more than 90% accurate.

Like other methods of the invention, a method of the invention based on comparison of at least two fibrotic markers can be used to diagnose the presence or severity of liver fibrosis in an individual having or suspected of having any liver disorder, including viral hepatitis, autoimmune liver disease such as autoimmune hepatitis, alcoholic liver disease, fatty liver disease or drug-induced liver disease, or any of the other liver diseases described herein above. Similarly a method of the invention based on comparison of at least two fibrotic markers can be used to diagnose the presence or severity of fibrotic disorders including pulmonary fibrosis, kidney fibrosis, prostate fibrosis, breast

fibrosis or a rheumatoid disease, or another fibrotic disorder described herein or known in the art.

A method of the invention relies on comparison of the level of a fibrotic marker to a predetermined 5 cut-off value. For markers that positively correlate with fibrosis, positivity is indicated by a level that is greater than the predetermined cut-off value. For markers that negatively correlate with fibrosis, positivity is indicated by a level that is less than the 10 predetermined cut-off value. Cut-off values useful in the methods of the invention can be determined as described herein, for example, using design of experiments (DOE) analysis.

As for the other diagnostic methods of the 15 invention, these methods can be practiced using a variety of fibrotic markers known in the art or described herein. Such fibrotic markers include, without limitation,  $\alpha$ 2-MG, HA, TIMP-1, PIIINP, laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 20 complex, sFas ligand, TGF- $\beta$ 1, IL-10, apoA1, apoA2 or ApoB. Additional serological, biochemical, clinical and echographic fibrotic markers are described herein above or are known in the art and can be included in any combination in a method of the invention. Furthermore, 25 it is understood that comparison of the first and second fibrotic markers and any additional fibrotic markers can be performed simultaneously or in any order and using any combination of assay formats.

As described above, the "level" of a fibrotic 30 marker can be a relative or absolute amount of, for

example, RNA, protein or activity and can be a direct or indirect measurement of the fibrotic marker. In addition, the value of the level can be obtained from a secondary source, such as a physician or diagnostic 5 laboratory or can be determined using any convenient sample and assay, including but not limited to those described herein above. Methods useful in determining the level of a fibrotic marker in order to perform the comparisons included in the methods of the invention 10 encompass, for example, hybridization methods such as RT-PCR and RNA blot analysis, immunoassays including enzyme-linked immunosorbent assays (ELISAs) and radioimmunoassays (RIAs), sandwich immunoassays, quantitative western blotting and other standard assays 15 for determining protein levels, and, where applicable, assays for the activity of the fibrotic marker. Such assays are routine in the art and described herein above.

The present invention further provides a method of diagnosing the presence or severity of liver fibrosis 20 in an individual by comparing a level of a first fibrotic marker X in the individual to a cut-off value X<sub>1</sub> to determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second fibrotic marker Y in the individual to a cut-off value Y<sub>1</sub> 25 to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where the cut-off values X<sub>1</sub> and Y<sub>1</sub> are optimized individually to give a 30 desired performance characteristic.

If desired, a method of the invention can include the steps of comparing a level of a third fibrotic marker Z in the individual to a cut-off value Z<sub>1</sub> to determine whether the individual is positive for the 5 third fibrotic marker Z and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z, where the cut-off values X<sub>1</sub>, Y<sub>1</sub> and Z<sub>1</sub> are optimized individually to give a desired performance characteristic. In one 10 embodiment, levels of  $\alpha$ 2-MG, HA and TIMP-1 are compared. In other embodiments, the levels of exactly three, at least three, at least four, or at least five fibrotic markers are compared. A method of the invention can be useful, for example, in differentiating no or mild liver 15 fibrosis from moderate to severe liver fibrosis. Cut-off values can be optimized as described herein, for example, using DOE analysis.

Further provided by the invention is a method of diagnosing the presence or severity of liver fibrosis 20 in an individual by comparing a level of a first fibrotic marker X in the individual to two cut-off values X<sub>1</sub> and X<sub>2</sub> to determine whether the individual is positive for the first fibrotic marker X; comparing a level of a second fibrotic marker Y in the individual to two cut-off 25 values Y<sub>1</sub> and Y<sub>2</sub> to determine whether the individual is positive for the second fibrotic marker Y; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X and Y, where the cut-off values X<sub>1</sub>, Y<sub>1</sub>, X<sub>2</sub> and Y<sub>2</sub> are optimized 30 individually to give a desired performance characteristic. Such performance characteristics include

particular sensitivities, specificities, PPVs, NPVs and accuracies, as described herein above.

A method of the invention can further include the steps of comparing a level of a third fibrotic marker Z in the individual to two cut-off values Z<sub>1</sub> and Z<sub>2</sub> to determine whether the individual is positive for the third fibrotic marker Z; and diagnosing the presence or severity of liver fibrosis in the individual based on positivity or negativity for X, Y and Z, where the cut-off values X<sub>1</sub>, Y<sub>1</sub>, Z<sub>1</sub>, X<sub>2</sub>, Y<sub>2</sub> and Z<sub>2</sub> are optimized individually to give a desired performance characteristic. In a method of the invention, cut-off values can be conveniently optimized, for example, using DOE analysis.

#### 15 Methodology

A variety of means can be useful for detecting  $\alpha_2$ -MG, HA and TIMP-1 and for determining a level of  $\alpha_2$ -MG, HA and TIMP in a sample. In one embodiment, the invention is practiced by determining the level of  $\alpha_2$ -MG protein in a sample from the individual to be diagnosed using, for example, one or more  $\alpha_2$ -MG-specific binding agents such as anti- $\alpha_2$ -MG antibodies. In another embodiment, a method of the invention is practiced by assaying for  $\alpha_2$ -MG activity in a sample from the individual.

A variety of means also can be used in a method of the invention to detect HA or determine a level of HA in a sample. In one embodiment, the invention is practiced by determining the level of HA in a sample

using one or more HA-specific binding agents such as HA-binding proteins or anti-HA antibodies.

Similarly, a variety of means can be used in a method of the invention to detect TIMP-1 or determine a 5 level of TIMP-1 in a sample. In one embodiment, the invention is practiced by determining the level of TIMP-1 protein in a sample from the individual to be diagnosed. The level of TIMP-1 protein can be determined, for example, using one or more TIMP-1-specific binding agents 10 such as anti-TIMP-1 antibodies. In another embodiment, the invention is practiced by assaying for TIMP-1 activity in a sample from the individual to be diagnosed.

In a particular embodiment, the invention provides a method of diagnosing the presence or severity 15 of liver fibrosis in an individual by determining the level of  $\alpha$ 2-MG protein in a sample from the individual; determining the level of HA in a sample from the individual; and determining the level of TIMP-1 protein in a sample from the individual; and diagnosing the 20 presence or severity of liver fibrosis in the individual based on the levels of  $\alpha$ 2-MG protein, HA and TIMP-1 protein. If desired, the level of  $\alpha$ 2-MG protein, HA and TIMP-1 protein each can be determined using an enzyme-linked assay.

25 In a further embodiment, the present invention provides a method of differentiating no or mild liver fibrosis from moderate to severe liver fibrosis in an individual by contacting an appropriate dilution of a sample from the individual with anti- $\alpha$ 2-MG antibody under 30 conditions suitable to form a first complex of  $\alpha$ 2-MG and

anti- $\alpha$ 2-MG antibody; washing the first complex to remove unbound molecules; determining the amount of  $\alpha$ 2-MG-containing first complex; contacting an appropriate dilution of a sample from the individual with a HA-  
5 binding protein under conditions suitable to form a second complex of HA and HA-binding protein; washing the second complex to remove unbound molecules; determining the amount of HA-containing second complex; contacting an appropriate dilution of a sample from the individual with  
10 anti-TIMP-1 antibody under conditions suitable to form a third complex of TIMP-1 and anti-TIMP-1 antibody; washing the third complex to remove unbound molecules; determining the amount of TIMP-1-containing third complex; and differentiating no or mild liver fibrosis  
15 from moderate to severe liver fibrosis in the individual based on the amounts of  $\alpha$ 2-MG, HA and TIMP-1-containing complexes.

It is understood that detecting  $\alpha$ 2-MG, HA and TIMP-1, or detecting  $\alpha$ 2-MG, HA and YKL-40, as discussed  
20 further below, can be accomplished by assaying for the amount of protein or polysaccharide directly, or, in the case of  $\alpha$ 2-MG and TIMP-1, can be determined by assaying for RNA levels or enzyme activity of a protease regulated by  $\alpha$ 2-MG or TIMP-1. Similarly, where one or more  
25 additional fibrotic markers is detected in a method of the invention, the marker can be assayed directly, or a precursor such as RNA, or a breakdown or proteolytic product, or an activity correlated with levels of the marker can be assayed. It is understood that determining  
30 a level of  $\alpha$ 2-MG, HA, TIMP-1 and YKL-40, or a level of any additional marker of fibrosis, can be performed using absolute values, for example, for RNA or protein levels

or enzyme activity, or can be determined as relative values in comparison to one or more reference values.

It further is understood that each of the three fibrotic marker assays ( $\alpha$ 2-MG/HA/TIMP-1 or 5  $\alpha$ 2-MG/HA/YKL-40), as well as any additional assays, is performed independently of the others, in any order, and that any combination of assay formats is encompassed by the invention. As an example, a level of  $\alpha$ 2-MG and HA can be determined by assaying for the concentration of 10  $\alpha$ 2-MG and HA while a level of TIMP-1 is determined by assaying for TIMP-1 enzyme activity. As another example, a level of  $\alpha$ 2-MG can be determined using a radioimmunoassay, while levels of HA and TIMP-1 are determined using enzyme-linked assays. One skilled in 15 the art understands that detection of the three fibrotic markers ( $\alpha$ 2-MG/HA/TIMP-1 or  $\alpha$ 2-MG/HA/YKL-40) and detection of any additional markers can be performed simultaneously or in any order. Furthermore, a single sample such as a serum sample can be obtained from an 20 individual and subdivided into three portions for detecting  $\alpha$ 2-MG, HA and TIMP-1 or  $\alpha$ 2-MG, HA and TIMP-1, or the markers can be detected using different samples, which can be of the same or a different type and can be undiluted or diluted to the same or different extents. 25 Where two or more samples are used, the samples are usually obtained from the individual within a relatively short time frame, for example, several days to several weeks.

RNA methods

Hybridization methods can be used to detect α2-MG or TIMP-1 mRNA or determine the level of α2-MG or TIMP-1 mRNA or the mRNA of another fibrotic marker useful in the invention such as YKL-40. Numerous methods are well known in the art for determining mRNA levels by specific or selective hybridization with a complementary nucleic acid probe. Such methods include solution hybridization procedures as well as solid-phase hybridization procedures in which the probe or sample is immobilized on a solid support. Specific examples of useful methods include amplification methods such as target and signal amplification methods and include PCR (polymerase chain reaction) and reverse-transcriptase-PCR (RT-PCR); transcription mediated amplification (Gen-Probe Incorporated; San Diego, CA); branched chain DNA (bDNA) amplification (Bayer Diagnostics; Emeryville, CA); strand displacement amplification (SDA; Becton Dickinson; Franklin Lakes, NJ); and ligase chain reaction (LCR) amplification (Abbott Laboratories; Abbott Park, IL). Additional methods useful in the invention include RNase protection; Northern analysis or other RNA blot, dot blot or membrane-based technology; dip stick; pin; and two-dimensional array immobilized onto a chip.

Conditions are well known in the art for quantitative determination of mRNA levels using both solution and solid phase hybridization procedures as described, for example, in Ausubel et al., Current Protocols in Molecular Biology (Supplement 47), John Wiley & Sons, New York (1999).

The polymerase chain reaction (PCR) RT-PCR can be useful in the methods of the invention. PCR or RT-PCR can be performed with isolated RNA or crude or partially fractionated samples, for example, cells pelleted from a 5 whole blood sample. PCR methods are well known in the art as described, for example, in Dieffenbach and Dveksler, PCR Primer: A Laboratory Manual, Cold Spring Harbor Press, Plainview, New York (1995). Multisample formats such as two-dimensional arrays offer the 10 advantage of analyzing numerous different samples in a single assay. Solid-phase dip stick-based methods also can be useful in the invention and offer the advantage of being able to rapidly analyze a fluid sample and obtain an immediate result.

15 Probes for detecting  $\alpha$ 2-MG and TIMP-1 mRNA or for determining  $\alpha$ 2-MG and TIMP-1 mRNA levels are well known in the art. One skilled in the art can use, for example, a probe corresponding to some or all of the human  $\alpha$ 2-MG nucleic acid sequence shown in Figure 1 (SEQ 20 ID NO: 1) or some or all of the human TIMP-1 nucleic acid sequence shown in Figure 3, respectively. Appropriate conditions for various assay formats for detecting  $\alpha$ 2-MG and TIMP-1 mRNA or for determining  $\alpha$ 2-MG and TIMP-1 mRNA levels are well known in the art or can be established 25 using routine methods. As an example, conditions and probes for Northern analysis of  $\alpha$ 2-MG RNA in human samples are described, for example, in Ortego et al., *supra*, 1997. As another example, conditions and probes for RNA slot blot hybridization to determine  $\alpha$ 2-MG RNA 30 expression in human samples are described in Simon et al., *supra*, 1996. Similarly, Northern analysis of TIMP-1 RNA in human samples can be performed as described, for

example, in Yoshiji et al., *supra*, 1996; RT-PCR assays for TIMP-1 in human samples also are well known in the art as described, for example, in Janowska-Wieczorek et al., *supra*, 2000, and Groft et al., *supra*, 2001. The 5 skilled person understands that these and other assays can be useful for detecting  $\alpha$ 2-MG, TIMP-1 or YKL-40 RNA or for determining  $\alpha$ 2-MG, TIMP-1 or YKL-40 RNA levels or the levels of other fibrotic markers useful in the methods of the invention.

10 Immunoassays

A variety of immunoassay formats, including competitive and non-competitive immunoassay formats, antigen capture assays and two-antibody sandwich assays also are useful the methods of the invention (Self and 15 Cook, *Curr. Opin. Biotechnol.* 7:60-65 (1996)). In one embodiment, a method of the invention relies on one or more antigen capture assays. In an antigen capture assay, antibody is bound to a solid phase, and sample is added such that  $\alpha$ 2-MG, HA, TIMP-1, YKL-40 or another 20 fibrotic marker antigen is bound by the antibody. After unbound proteins are removed by washing, the amount of bound antigen can be quantitated, if desired, using, for example, a radioassay (Harlow and Lane, *Antibodies A Laboratory Manual* Cold Spring Harbor Laboratory: New 25 York, 1988)). One skilled in the art understands that immunoassays useful in the invention are performed under conditions of antibody excess, or as antigen competitions, to quantitate the amount of antigen and, thus, determine a level of  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40.

30 Enzyme-linked immunosorbent assays (ELISAs) can be useful in the methods of the invention. An enzyme

such as horseradish peroxidase (HRP), alkaline phosphatase (AP),  $\beta$ -galactosidase or urease can be linked, for example, to an anti- $\alpha$ 2-MG, anti-HA, anti-TIMP-1 or anti-YKL-40 antibody or to a secondary antibody for use in a method of the invention. A horseradish-peroxidase detection system can be used, for example, with the chromogenic substrate tetramethylbenzidine (TMB), which yields a soluble product in the presence of hydrogen peroxide that is detectable at 450 nm. Other convenient enzyme-linked systems include, for example, the alkaline phosphatase detection system, which can be used with the chromogenic substrate *p*-nitrophenyl phosphate to yield a soluble product readily detectable at 405 nm. Similarly, a  $\beta$ -galactosidase detection system can be used with the chromogenic substrate *o*-nitrophenyl- $\beta$ -D-galactopyranoside (ONPG) to yield a soluble product detectable at 410 nm, or a urease detection system can be used with a substrate such as urea-bromocresol purple (Sigma Immunochemicals, St. Louis, MO). Useful enzyme-linked primary and secondary antibodies can be obtained from a number of commercial sources such as Jackson Immuno-Research (West Grove, PA) as described further below.

Chemiluminescent detection also can be useful for detecting  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40 or for determining a level of  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40 or another fibrotic marker according to a method of the invention. Chemiluminescent secondary antibodies can be obtained commercially from various sources such as Amersham.

Fluorescent detection also can be useful for detecting  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40 or for determining a level of  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40 or another fibrotic marker in a method of the invention. Useful 5 fluorochromes include, without limitation, DAPI, fluorescein, Hoechst 33258, R-phycocyanin, B-phycoerythrin, R-phycoerythrin, rhodamine, Texas red and lissamine. Fluorescein or rhodamine labeled  $\alpha$ 2-MG-, HA-, TIMP-1- or YKL-40-specific binding agents such as 10 anti- $\alpha$ 2-MG, anti-HA, anti-TIMP-1, or anti-YKL-40 antibodies, or fluorescein- or rhodamine-labeled secondary antibodies can be useful in the invention. Useful fluorescent antibodies can be obtained 15 commercially, for example, from Tago Immunologicals (Burlingame, CA) as described further below.

Radioimmunoassays (RIAs) also can be useful in the methods of the invention. Such assays are well known in the art. For example, Brophy et al., Biochem. Biophys. Res. Comm. 167:898-903 (1990)), describes a 20 radioimmunoassay for detection of TIMP-1, and Pharmacia makes a radiometric assay for quantitation of HA using an  $^{125}$ I-labelled HA-binding protein (Guechot et al., Clin. Chem. 42:558-563 (1996)). Radioimmunoassays can be performed, for example, with  $^{125}$ I-labeled primary or 25 secondary antibody (Harlow and Lane, *supra*, 1988).

A signal from a detectable reagent can be analyzed, for example, using a spectrophotometer to detect color from a chromogenic substrate; a radiation counter to detect radiation, such as a gamma counter for 30 detection of  $^{125}$ I; or a fluorometer to detect fluorescence in the presence of light of a certain wavelength. Where

an enzyme-linked assay is used, quantitative analysis of the amount of  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40 or another fibrotic marker can be performed using a spectrophotometer such as an EMAX Microplate Reader 5 (Molecular Devices; Menlo Park, CA) in accordance with the manufacturer's instructions. It is understood that the assays of the invention can be automated or performed robotically, if desired, and that the signal from multiple samples can be detected simultaneously.

10 The methods of the invention also encompass the use of capillary electrophoresis based immunoassays (CEIA), which can be automated, if desired. Immunoassays also can be used in conjunction with laser-induced fluorescence as described, for example, in Schmalzing and 15 Nashabe, Electrophoresis 18:2184-93 (1997), and Bao, J. Chromatogr. B. Biomed. Sci. 699:463-80 (1997). Liposome immunoassays, such as flow-injection liposome immunoassays and liposome immunosensors, also can be used to detect  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40 or to determine a 20 level of  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40 or another fibrotic marker according to a method of the invention (Rongen et al., J. Immunol. Methods 204:105-133 (1997)).

Sandwich enzyme immunoassays also can be useful in the methods of the invention. In a two-antibody 25 sandwich assay, a first antibody is bound to a solid support, and the antigen is allowed to bind to the first antibody. The amount of  $\alpha$ 2-MG, HA, TIMP-1, YKL-40 or another fibrotic marker antigen is quantitated by measuring the amount of a second antibody that binds the 30 fibrotic marker.

As an example, a two-antibody sandwich immunoassay can be useful to determine a level of TIMP-1 as described in Murawaki et al., *supra*, 1993. Briefly, serum (25 µl) is diluted 41-fold with 10 mM sodium phosphate buffer, pH 7.0 (1.0 ml). The diluted sample (20 µl) is mixed with 0.3 ml of 10 mM sodium phosphate buffer, pH 7.0, containing 50 ng/ml monoclonal antibody (Fab of clone 7-6C1) labeled with horseradish peroxidase, 1% bovine serum albumin, 0.1% Tween 20, 0.1 M NaCl and 0.005% thimerosal. A 0.1 ml aliquot of the mixed solution is transferred to each microplate well previously coated with a second monoclonal antibody (clone 7-23G9) having a different epitope specificity, and the plate incubated for 30 minutes at room temperature without shaking. The plate is washed three times with 0.3 ml 10 mM sodium phosphate buffer, pH 7.0, containing 0.1% Tween 20 and 0.1 M NaCl. Peroxidase activity bound to the plate is assayed by a 15 minute incubation at room temperature with 0.1 ml 0.15 M citric acid sodium phosphate buffer, pH 4.9, containing 0.5 mg/ml o-phenylenediamine and 0.02% H<sub>2</sub>O<sub>2</sub>. After stopping the reaction by addition of 0.1 ml 2 N H<sub>2</sub>SO<sub>4</sub>, the absorbance at 492 nm is measured in a microplate reader using a standard of human serum TIMP-1. Linearity between the amount of TIMP-1 and absorbance at 492 nm is demonstrated by graphing with logarithmic scales and yields an assay range of about 1.5 to 300 µg/well.

Quantitative western blotting also can be used to detect α2-MG, HA, TIMP-1 or YKL-40 or to determine a level of α2-MG, HA, TIMP-1 or YKL-40 or a level of another fibrotic marker antigen in a method of the invention. Western blots can be quantitated by well

known methods such as scanning densitometry. As an example, protein samples are electrophoresed on 10% SDS-PAGE Laemmli gels. Primary murine monoclonal antibodies, for example, against human  $\alpha$ 2-MG, HA, TIMP-1 or YKL-40 are reacted with the blot, and antibody binding confirmed to be linear using a preliminary slot blot experiment. Goat anti-mouse horseradish peroxidase-coupled antibodies (BioRad) are used as the secondary antibody, and signal detection performed using chemiluminescence, for example, with the Renaissance chemiluminescence kit (New England Nuclear; Boston, MA) according to the manufacturer's instructions.

Autoradiographs of the blots are analyzed using a scanning densitometer (Molecular Dynamics; Sunnyvale, CA) and normalized to a positive control. Values are reported, for example, as a ratio between the actual value to the positive control (densitometric index). Such methods are well known in the art as described, for example, in Parra et al., J. Vasc. Surg. 28:669-675 (1998).

#### Sources for antibodies

As described herein above, immunoassays including but not limited to enzyme-linked immunosorbent assays, radioimmunoassays and quantitative western analysis, can be useful in the diagnostic methods of the invention. Such assays rely on one or more antibodies, for example, anti- $\alpha$ 2-MG, anti-HA, anti-TIMP-1 or anti-YKL-40 antibodies. As used herein, the term "antibody" is used in its broadest sense to include polyclonal and monoclonal antibodies, as well as polypeptide fragments of antibodies that retain binding activity for  $\alpha$ 2-MG, HA,

TIMP-1, YKL-40 or the relevant fibrotic marker antigen of at least about  $1 \times 10^5 M^{-1}$ . One skilled in the art understands that antibody fragments such as anti- $\alpha 2$ -MG, anti-HA, anti-TIMP-1 and anti-YKL-40 antibody fragments 5 and including Fab, F(ab')<sub>2</sub> and Fv fragments can retain binding activity for the relevant fibrotic marker antigen and, thus, are included within the definition of the term antibody as used herein. Methods of preparing monoclonal and polyclonal antibodies are routine in the art, as 10 described, for example, in Harlow and Lane, *supra*, 1988.

The term antibody, as used herein, also encompasses non-naturally occurring antibodies and fragments containing, at a minimum, one V<sub>H</sub> and one V<sub>L</sub> domain, such as chimeric antibodies, humanized antibodies 15 and single chain Fv fragments (scFv) that specifically bind  $\alpha 2$ -MG, HA, TIMP-1, YKL-40 or the relevant fibrotic marker antigen. Such non-naturally occurring antibodies can be constructed using solid phase peptide synthesis, produced recombinantly or obtained, for example, by 20 screening combinatorial libraries consisting of variable heavy chains and variable light chains as described by Borrebaeck (Ed.), Antibody Engineering (Second edition) New York: Oxford University Press (1995).

A variety of useful anti- $\alpha 2$ -MG, anti-HA, 25 anti-TIMP-1 and anti-YKL-40 monoclonal and polyclonal antibodies are well known in the art and, in many cases, are commercially available. For example, a nephelometry assay for  $\alpha 2$ -macroglobulin is available from Beckman Coulter (kit #449430), and affinity purified goat 30 anti-human  $\alpha 2$ -MG and peroxidase-labeled goat anti-human  $\alpha 2$ -MG antibodies suitable for ELISA and western blotting

are available, for example, from Cedarlane Laboratories Limited (CL20010AP and CL20010APH) and Affinity Biologicals Incorporated (GAA2M-AP and GAA2M-APH). Similarly, affinity purified sheep anti-HA antiserum can 5 be obtained from Biotrend (#5029-9990).

Anti-human TIMP-1 antibodies also are readily available from a variety of commercial sources. For example, the anti-human TIMP-1 monoclonal antibody 147-6D11 is suitable for ELISA or western blotting 10 analysis and can be obtained from Medicorp, Inc. (Montreal, Canada), and the anti-human TIMP-1 monoclonal antibody MAB970 is available from R&D Systems, Inc., for use, for example, in western blotting or sandwich ELISA assays. MAB970 can be combined, for example, with 15 biotinylated anti-human TIMP-1 antibody (BAF970) from R&D Systems, Inc., for detection of TIMP-1 by sandwich ELISA. In addition, rabbit anti-human TIMP-1 polyclonal 20 antiserum and mouse anti-human monoclonal antibodies suitable, for example, for western blotting with enhanced chemiluminescence detection can be obtained from Research Diagnostics Inc. (RDI-TIMP1abr and RDI-TIMP1-C1).

#### Assays for activity

As discussed above, assays based on the activity of a fibrotic marker also can be useful for 25 detecting  $\alpha$ 2-MG, HA or TIMP-1 or for determining a level of  $\alpha$ 2-MG, HA or TIMP-1 or another fibrotic marker and, therefore, are useful in the methods of the invention. As an example, a variety of assays for  $\alpha$ 2-MG activity can be useful for detecting  $\alpha$ 2-MG or determining a level of 30  $\alpha$ 2-MG in a sample in a method of the invention. Because

$\alpha_2$ -MG-bound proteases display inhibited proteolytic activity but retain the ability to hydrolyze amide and ester bonds of small substrates,  $\alpha_2$ -MG can be detected, or a level determined, by assaying for inhibition of trypsin, subtilisin, chymotrypsin, plasmin, elastase, thermolysin, or papain activity or the activity of another target protease without inhibition of amidolytic activity. Substrates such as labeled casein or labeled fibrin can be useful for assaying for inhibition of target protease activity. Furthermore, based on its broad protease substrate specificity, a level of  $\alpha_2$ -MG can be determined by assaying for inhibition of the activity of two or more target proteases using, for example,  $^{14}\text{C}$ -casein and  $^{125}\text{I}$ -fibrin (Armstrong et al., *supra*, 1999).  $\alpha_2$ -MG also can be detected or a level of  $\alpha_2$ -MG determined based on the ability of  $\alpha_2$ -MG to shield a bound protease from an antibody or a high molecular weight inhibitor. Following reaction of a sample with, for example, trypsin and then trypsin inhibitor, residual trypsin activity is assayed with a low molecular mass substrate such as the amide BApNA (Ganrot, *supra*, 1966; Armstrong et al., *supra*, 1985). Trypsin activity following treatment with trypsin inhibitor is indicative of  $\alpha_2$ -MG. These and other well known assays for  $\alpha_2$ -MG activity can be useful in the methods of the invention.

Similarly, assays for TIMP-1 activity are well known in the art. In particular, one assays for the ability to inhibit protease activity of a matrix metalloproteinase, for example, using reverse gelatin zymography. Reverse gelatin zymography is performed by including a gelatinase such as gelatinase A in a gel mix with the gelatin substrate. Conditioned media, such as

conditioned media from baby hamster kidney cells can be used as a convenient source of gelatinase. Plasma samples are electrophoresed, and the resulting pattern analyzed, for example, with scanning digitization using a 5 Hewlett Packard scanner. TIMP-1 activity is observed as a reduction of gelatin degradation. See, for example, Kossakowska et al., *supra*, 1998. The skilled person recognizes that these and other routine assays for TIMP-1 activity can be useful in the methods of the invention.

10 Additional markers

It is clear that the methods of the invention can be practiced, if desired, by detecting the three markers  $\alpha$ 2-MG, HA and TIMP-1 without assaying for any additional markers or evaluating any other clinical or 15 echographic characteristics. In addition, these three assays can be used as a panel in combination with one or more additional fibrotic marker assays or evaluation of one or more clinical or echographic variables. In specific embodiments, the invention provides a method of 20 diagnosing the presence or severity of liver fibrosis in an individual by detecting  $\alpha$ 2-MG, HA and TIMP-1 in a sample and also detecting at least one of the following markers: PIIINP, laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 25 complex, sFas ligand, TGF- $\beta$ 1, IL-10, apoA1, apoA2 or ApoB. In one embodiment, a method of the invention for diagnosing the presence or severity of liver fibrosis includes the steps of detecting  $\alpha$ 2-MG, HA, TIMP-1 and YKL-40 in a sample. In a further embodiment, a method of

the invention is limited to detecting  $\alpha$ 2-MG, HA, TIMP-1 and YKL-40, and no additional fibrotic markers are detected.

In view of the above, it is clear that assays  
5 for one or more additional biochemical or serological  
markers of fibrosis or evaluation of one or more clinical  
or echographic variables associated with fibrosis can be  
combined with detection of  $\alpha$ 2-MG, HA, and TIMP-1 to  
diagnose the presence or severity of liver fibrosis.  
10 Examples of additional biochemical and serological  
markers include, yet are not limited to, PIIINP, laminin,  
tenascin, collagen type IV, collagen type VI, YKL-40,  
MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- $\beta$ 1,  
IL-10, apoA1, apoA2 and apoB. Additional biochemical and  
15 serological markers useful in the invention include,  
without limitation, fibronectin, pseudocholinesterase,  
manganese superoxide dismutase, *N*-acetyl-  
 $\beta$ -glucosaminidase ( $\beta$ -NAG), glutathione peroxidase,  
connective tissue growth factor (CTGF); platelet derived  
20 growth factor (PDGF), PDGF receptor, inducible nitric  
oxide synthetase, nitrotyrosine, bilirubin, ferritin and  
 $\alpha$ -fetoprotein,  $\gamma$ -glutamyl transpeptidase (GGT), aspartate  
aminotransferase (AST), alanine aminotransferase (ALT),  
AST/ALT ratio, albumin,  $\gamma$ -globulins,  $\beta\gamma$ -block,  
25 prothrombin index, Child-Pugh score, PGA index  
(prothrombin time, GGT concentration and apoA1  
concentration), PGAA index (PGA score with  
 $\alpha$ 2-macroglobulin level), hemoglobin, mean corpuscular  
volume, lymphocyte count, cholesterol, urea, creatinine,  
30 sodium and platelet count.

A clinical or echographic variable also can be a fibrotic "marker" useful in the methods of the invention. Thus, analysis of one or more clinical or echographic variables can be combined with detection of

5  $\alpha$ 2-MG, HA and TIMP-1 to diagnose the presence or severity of liver fibrosis, or another fibrotic disorder as described hereinabove. As examples, such a clinical variable can be patient age or gender or the presence of palmar erythema, Dupuytren's contracture, finger

10 clubbing, spider nevi, firm liver, splenomegaly or collateral circulation. Echographic variables useful in a method of the invention include, for example, liver length (right kidney), irregular liver surface, liver heterogeneity, spleen length, ascites or collateral

15 circulation. See, for example, Oberti et al., Gastroenterol. 113:1609-1616 (1997). It is understood that the analysis of these and other well known clinical or echographic variables can be useful in a method of the invention. Furthermore, a method of the invention

20 encompasses determination of the clinical or echographic variable, for example, liver palpation, or can rely on one or more historic, or previously determined clinical or echographic variables.

Assays for detection of biochemical or

25 serological markers useful in the invention are well known in the art and in many cases commercially available. Such assays include, but are not limited to, amplification based methods such as RT-PCR and other methods for quantitative analysis of RNA levels;

30 immunoassays such as radioimmunoassays, enzyme-linked assays, two-antibody sandwich assays and quantitative western analysis; and assays for biological activity such

as enzyme activity. Assays for PIIINP, laminin, tenascin, collagen type IV, collagen type VI, YKL-40, MMP-3, MMP-2, MMP-9/TIMP-1 complex, sFas ligand, TGF- $\beta$ 1, IL-10, apoA1, apoA2 and apoB are commercially available  
5 from various sources as summarized in Table 1.

TABLE 1 COMMERCIAL SOURCES FOR FIBROTIC MARKER ASSAYS				
	Marker	Company	Assay	Catalog number
5	PIIINP	Orion Diagnostica (Espoo, Finland)	RIA	05903
	laminin	Chemicon Intl. (Temecula, CA)	ELISA	ECM310
	tenascin	Chemicon Intl. (Temecula, CA)	ELISA	ECM320
	collagen IV	Iatron Laboratories (Tokyo, Japan)	RIA	KCAD1
	YKL-40	Metra Biosystems (Mountain View, CA)	ELISA	8020
	MMP-3	Amersham Pharmacia (Piscataway, NJ)	ELISA	RPN 2613
	MMP-2	Amersham Pharmacia (Piscataway, NJ)	ELISA	RPN 2617
	MMP-9/TIMP-1 complex	SBA Sciences (Turku, Finland)	ELISA	MP2215
	sFas ligand	Bender MedSystems Diagnostics (Vienna, Austria)	ELISA	BMS260/2
	TGF- $\beta$ 1	R&D Systems (Minneapolis, MN)	ELISA	DB100
15	IL-10	R&D Systems (Minneapolis, MN)	ELISA	HS100B
	apoA1	AlerChek, Inc. (Portland, ME)	ELISA	A70101
	apoA2	AlerChek, Inc. (Portland, ME)	ELISA	A70102
	apoB	Sigma Diagnostics (St. Louis, MO)	IT*	357-A
* designates immunoturbidimetric				

Assays for additional biochemical or serological markers that can be combined with detection of  $\alpha_2$ -MG, HA and TIMP-1 in a method of the invention also are well known in the art. Fibronectin, for example, can 5 be conveniently assayed by turbidimetric assay available from Roche Diagnostics (Mannheim, Germany). Pseudocholinesterase (PCHE) can be assayed using standard methodology available from Boehringer. Levels of 10  $N$ -acetyl-  $\beta$ -glucosaminidase ( $\beta$ -NAG) can be determined by assaying for enzymatic activity using a kit available 15 from Cortecs diagnostics. Manganese superoxide dismutase (Mn-SOD) levels can be conveniently determined by ELISA using a kit available, for example, from Bender MedSystem. Glutathione peroxidase levels can be determined by assaying for enzymatic activity using, for 20 example, a kit available from Randox Laboratories Ltd (Oceanside, CA).

Total or direct bilirubin, GGT, AST and ALT levels can be determined using an autoanalyser such as 25 Hitachi 917 Automate (Mannheim, Germany) with Roche Diagnostics reagents. Albumin levels can be determined, for example, by the bromocresol green method as described in Doumas et al., Clin. Chim Acta 31:87-96 (1971), and ferritin and  $\alpha$ -fetoprotein levels can be conveniently 30 determined using, for example, an immunoassay available from Boehringer. In addition, levels of  $\alpha_1$  globulin,  $\alpha_2$  globulin,  $\beta$  globulin and  $\gamma$ -globulin can be determined, for example, by serum protein electrophoresis in an automatic system (Hydrasys and Hyrys, Sebia; Issy-Les-Moulineaux, France). Methods of determining prothrombin activity also are well known in the art and include the clotting method available from Organon

Technika (West Orange, NJ). PGA index can be determined as described in Poynard et al., Gastroenterol. 100:1397-1402 (1991), and PGAA index also can be determined by well known methods as described in Naveau et al., Dig. 5 Dis. Sci. 39:2426-2432 (1994)).

Platelet counts, lymphocyte counts, mean corpuscular volume and related variables can be determined by a variety of methodologies using, for example, a Bayer-Technicon H2 analyser (Bayer-Technicon 10 Instruments; Tarrytown, NY). Cholesterol levels can be determined by standard methodologies available, for example, from Boehringer. Thus, it is clear to the skilled person that a variety of methodologies, including but not limited to the above, are well known in the art 15 and can be useful in the diagnostic methods of the invention.

$\alpha$ 2-MG/HA/YKL-40 panel

The present invention also provides a method of diagnosing the presence or severity of liver fibrosis in 20 an individual by detecting  $\alpha$ 2-MG in a sample; detecting HA in a sample; detecting YKL-40 in a sample; and diagnosing the presence or severity of liver fibrosis in the individual based on the presence or level of  $\alpha$ 2-MG, HA and YKL-40. A method of the invention can be useful, 25 for example, or differentiating no or mild (F0-F1) liver fibrosis from moderate to severe (F2-F4) liver fibrosis.

In one embodiment, the invention provides a method of diagnosing the presence or severity of liver fibrosis in an individual by determining the level of

α2-MG protein in a sample from the individual; determining the level of HA in a sample from the individual; determining the level of YKL-40 protein in a sample from the individual; and diagnosing the presence 5 or severity of liver fibrosis in the individual based on the levels of α2-MG protein, HA and YKL-40 protein. If desired, the levels of α2-MG protein, HA and YKL-40 protein each can be determined using an enzyme-linked assay.

10 Thus, in particular embodiments, the present invention provides diagnostic methods which rely, in part, on determining a level of the fibrotic marker YKL-40 in a sample. YKL-40, also known as human cartilage glycoprotein 39 (HC gp-39), is named for a 15 molecular weight of 40 kDa and the amino-terminal sequence of the protein, tyrosine-lysine-leucine (YKL). This glycoprotein, a mammalian member of the chitinase family (18-glycosylhydrolases), is a lectin that binds heparin and chitin and is produced by chondrocytes, 20 synovial cells, activated macrophages, neutrophils and MG-63 osteosarcoma cells (Hakala et al., J. Biol. Chem. 268:25803-15810 (1993); Nyirkos and Golds, Biochem. J. 268:265-268 (1990); Renkema et al., Eur. J. Biochem. 251:504-509 (1998); Volck et al., Proc. Assoc. Am. 25 Physicians 110:351-360 (1998); and Johansen et al., J. Bone Miner. Res. 7:501-511 (1992)). The pattern of YKL-40 expression in normal and diseased tissue indicates that this glycoprotein can function in extracellular matrix remodelling or tissue inflammation (Nyirkos and 30 Golds, *supra*, 1990; Renkema et al., *supra*, 1998; and Verheijden et al., Arthritis Rheum. 40:1115-1125 (1997)). Furthermore, YKL-40 mRNA is expressed in liver, and

initial studies have shown that YKL-40 expression is elevated in patients with chronic liver disease and that increased serum YKL-40 can be associated with fibrosis and fibrogenesis (Johansen et al., Scand. J. 5 Gastroenterol. 32:582-590 (1997); and Johansen, J. Hepatol. 32:911-920 (2000)).

Methods of determining a level of YKL-40 in samples such as serum and synovial fluid are well known in the art. For example, a radioimmunoassay for YKL-40 10 based on a rabbit antibody raised against YKL-40 is described in Johansen et al., Br. J. Rheumatology 32:949-955 (1993). In addition, a sandwich enzyme immunoassay in a microliter stripwell format is commercially available from Metra Biosystems. In the Metra Biosystems 15 assay, the Fab fragment of a biotin-conjugated monoclonal anti-YKL-40 antibody binds to streptavidin on the strip and captures YKL-40 in a sample. Alkaline phosphatase-conjugated polyclonal anti-YKL-40 antiserum binds the captured YKL-40 antigen, and alkaline phosphatase 20 activity is detected with *p*-nitrophenyl phosphate substrate as an indication of YKL-40 concentration. It is understood that the methods of the invention can be practiced with these or other routine assays for detecting or determining a level of YKL-40 RNA or 25 protein.

The following examples are intended to illustrate but not limit the present invention.

**EXAMPLE I****MARKER PANELS FOR NON-INVASIVE DIAGNOSIS**  
**OF LIVER FIBROSIS**

This example demonstrates that several

- 5 serological markers can be combined together as a panel  
which is useful in differentiating fibrosis stages F2, F3  
and F4 from stages F0 and F1 in HCV-infected patients.

Serum samples from 194 HCV patients positive  
for hepatitis C virus by RNA and immuno-analysis and  
10 having elevated alanine aminotransferase (ALT) levels  
were chosen at random from an existing serum library.  
Each of the patients had had a liver biopsy as part of  
their care. Patient samples were chosen to allow  
comparison of other routine blood markers and physical  
15 examination results attendant to routine medical care,  
including HCV viral load.

The inclusion criteria for the study were that  
the patient 1) have a confirmed hepatitis C infection at  
the time of liver biopsy and serum draw; 2) have  
20 undergone a liver biopsy as part of their medical care  
independent of the study; and 3) have previously given  
informed consent. Patients who did not give informed  
consent or who were incarcerated were excluded from the  
study.

25 Fibrosis scores (Metavir stage) for the 194  
patients were established by histopathologic examination  
of a needle biopsy specimen prior to therapy according to  
the criteria set forth in The French Metavir Cooperative  
Study Group, Hepatol. 20:15-20 (1994). All Metavir

fibrosis scores were established by the same pathologist. For all analyses, Metavir scores of F0 and F1 were grouped together as "no/mild" fibrosis, while scores of F2, F3 and F4 were grouped together as "moderate/severe" 5 fibrosis. The fibrosis prevalence in the 194 patient group was determined to be 60% was based on the proportion of F2-F4 scores in the group as shown in Table 2.

10

TABLE 2  
COMPOSITION OF THE 194 HCV PATIENT STUDY  
POPULATION BY FIBROSIS STAGE

15

Fibrosis stage	Number	Total F0-F1 or F2-F4
F0	38	F0-F1 = 78
F1	40	
F2	40	F2-F4 = 116
F3	39	
F4	37	
Total	194	Prevalence = 59.8%

20

As shown in the table above, the panel of HCV patient samples included 37 samples with very high fibrosis stage (F4); 39 samples from patients with very low or zero fibrosis stage (F0); and 158 samples from patients with fibrosis stage F1, F2 or F3.

25

Serum samples were assayed for the presence of several putative fibrosis markers, including laminin, YKL-40, HA, TIMP-1, PIIINP, type IV collagen and  $\alpha$ 2-MG. Assays were performed using commercial kits according to manufacturers' instructions (see Table 3). The results obtained for the 194 samples analyzed for laminin,

YKL-40, HA, TIMP-1, PIIINP, collagen type IV and  $\alpha$ 2-MG are shown in Table 4.

5

TABLE 3  
COMMERCIALLY AVAILABLE KITS FOR DETECTION OF  
FIBROSIS MARKERS

10

Marker	Manufacturer	Assay type	Catalogue number
Laminin	Chemicon Intl. (Temecula, CA)	ELISA	ECM310
YKL-40	Metra Biosystems (Mountain View, CA)	ELISA	8020
HA	Corgenix (Westminster, CO)	ELISA	029001
TIMP-1	Amersham Pharmacia (Piscataway, NJ)	ELISA	RPN 2611
PIIINP	Orion Diagnostica (Espoo, Finland)	RIA	05903
collagen IV	Iatron Laboratories (Tokyo, Japan)	RIA	KCAD1
$\alpha$ 2-MG	Beckman Coulter	Nephelometry	449430

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
5	ID							
100010	B-A	175.244	81.608	15.730	1308.802	2.288	1.737	3.03
100038	P-B	151.888	67.220	9.288	917.104	2.049	2.617	2.01
100044	C-B	187.811	60.757	44.127	1610.690	3.883	3.408	5.03
100059	T-B	232.082	51.002	22.583	1077.343	2.297	1.901	2.29
10 100069	N-C	285.269	131.726	73.851	2381.222	8.034	2.954	4.05
100077	H-C	268.685	47.709	18.066	1122.818	2.260	3.159	1.75
100090	B-C	263.426	26.370	47.339	1380.182	3.526	2.561	3.30
100127	D-D	279.580	166.113	105.505	1180.879	3.343	2.804	3.54
100167	G-F	274.533	482.708	341.132	2523.637	9.745	5.110	3.98
15 100175	B-G	266.903	95.808	27.721	1178.105	4.345	3.911	2.67
100178	M-G	211.613	159.040	25.669	1176.718	2.357	3.795	1.93
100182	T-G	246.686	55.391	8.889	1308.815	2.924	2.468	2.86
100198	A-G	226.372	48.441	13.901	1126.962	2.595	0.819	3.03
100209	J-G	288.524	83.925	5.051	1081.470	5.173	0.801	2.73
20 100229	S-H	253.561	110.020	46.568	1391.433	4.905	4.410	3.12

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
100238	T-J	229.781	38.076	29.516	1190.567	2.626	3.141	2.68
100245	D-K	279.768	270.250	171.481	2310.561	5.876	4.713	4.01
100247	C-K	244.559	131.482	10.219	1405.454	2.297	1.089	2.50
100250	J-K	262.136	101.729	54.821	1155.963	4.192	2.655	4.04
5 100252	J-K	260.998	61.366	57.275	1560.856	3.498	7.040	3.24
100253	E-K	292.189	173.917	168.768	1652.033	9.252	7.336	3.59
100254	W-K	288.551	477.560	102.775	1580.756	5.300	4.188	4.01
100271	M-L	278.201	89.900	69.651	1140.761	4.092	4.066	2.94
100276	V-L	257.309	176.112	12.196	1088.369	1.985	4.147	3.57
10 100284	M-L	224.339	130.263	31.822	1104.885	3.653	2.542	3.14
100290	D-L	199.542	541.552	50.429	1550.943	4.399	6.368	3.40
100294	TRL	281.501	217.328	200.436	2340.630	11.006	7.016	4.20
100301	P-M	285.543	430.475	29.772	1884.061	4.453	3.524	4.26
100313	J-M	301.751	188.062	45.539	1852.125	3.610	4.773	4.22
15 100323	M-M	223.002	626.140	144.334	2382.232	8.873	8.042	4.09

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
100334	K-R	173.320	63.317	38.516	1290.231	4.327	4.753	2.46
100339	S-M	184.085	36.125	12.049	1268.153	1.865	2.211	2.39
100340	K-M	206.582	57.496	22.015	974.475	2.315	0.819	2.80
100341	T-M	257.580	76.834	91.748	1492.882	3.976	4.793	3.00
5 100343	D-M	334.202	140.629	49.322	1098.199	4.092	2.412	3.12
100357	P-M	419.291	27.368	14.971	784.932	4.225	0.801	2.01
100374	K-N	300.366	28.231	36.608	1697.678	5.478	6.690	3.16
100379	T-O	233.496	75.711	26.906	1437.939	2.086	2.437	3.70
100382	C-P	206.796	44.461	6.034	989.007	2.214	4.409	2.14
10 100397	R-R	223.006	66.474	13.912	981.736	4.091	3.814	3.53
100410	R-S	224.775	36.605	37.499	1152.258	4.592	6.466	3.43
100438	D-Q	228.008	149.349	75.452	1682.636	7.734	7.039	3.89
100451	H-R	248.528	526.840	226.386	1961.029	10.505	7.113	3.36
100453	A-R	225.862	65.956	33.169	1421.632	2.786	3.465	3.79
15 100454	O-R	220.481	56.892	32.531	1125.960	3.003	7.014	3.68

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll I V (ng/ml)	$\alpha$ 2-MG (ng/ml)
100456	S-S	241.591	46.274	30.745	1337.355	5.182	2.411	2.81
100466	S-S	210.562	48.605	34.443	1482.475	4.286	4.754	3.27
100470	D-S	229.912	162.039	55.053	1684.159	6.806	3.640	3.10
100485	C-T	229.811	113.523	38.389	1247.547	3.901	3.291	4.11
5 100486	M-T	265.326	281.257	706.557	2716.589	15.362	11.975	2.82
100505	L-W	229.363	33.843	23.305	1149.367	3.397	2.889	3.29
100519	R-W	204.646	68.632	10.431	845.571	3.852	6.815	2.52
100547	S-G	223.959	75.711	8.257	1000.623	4.286	3.090	2.92
100638	J-P	265.819	264.250	68.361	2095.698	15.945	7.418	4.67
10 100640	M-V	170.293	43.770	17.728	1200.584	4.755	5.561	3.50
100006	L-A	135.628	75.349	79.430	1354.782	6.612	3.477	3.14
100009	R-A	157.239	72.429	21.947	932.635	2.080	3.050	1.96
100011	A-B	136.197	251.237	149.932	2004.294	7.600	4.853	3.57
100016	E-A	161.133	272.434	186.536	1900.341	9.341	9.071	3.08
15 100021	E-AV	184.000	537.630	102.420	2456.883	4.863	6.157	3.97

10087188 .062102

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll I V (ng/ml)	$\alpha$ 2-MG (ng/ml)
100023	C-B	126.346	194.523	47.976	1540.914	7.000	4.488	3.21
100027	D-B	133.660	75.820	33.912	1519.528	2.966	3.286	3.51
100030	R-B	140.584	50.007	153.135	1219.549	3.237	5.200	1.89
100035	K-B	124.645	37.383	60.934	1214.060	3.582	3.620	2.41
5 100036	G-B	152.864	87.596	369.681	1305.790	3.163	4.391	2.71
100041	M-B	168.422	42.376	143.372	1502.562	6.667	3.692	3.10
100042	M-B	138.754	211.387	266.568	2899.870	8.233	8.559	3.98
100043	C-B	111.743	30.883	17.447	1168.327	5.488	2.343	2.76
100045	V-B	164.940	241.063	221.249	2010.088	9.097	4.512	3.76
10 100051	K-B	154.743	222.409	131.122	1600.554	4.863	4.075	3.43
100055	D-B	146.817	110.018	84.447	1827.668	3.188	4.439	5.72
100065	G-B	134.349	72.429	112.148	1455.905	3.353	3.002	4.02
100071	R-C	135.011	74.407	30.352	1485.573	2.820	3.120	3.29
100073	G-C	146.785	63.761	43.312	1530.873	3.027	2.040	< 0.75
15 100074	L-C	151.514	80.248	49.917	1647.700	5.036	3.286	3.11

10087188 - 062102

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
100078	R-C	163.144	213.365	45.839	1399.880	3.393	2.297	2.61
100081	A-C	147.915	45.862	56.686	1346.315	4.056	4.707	2.37
100084	G-C	144.665	43.130	116.238	1736.670	5.337	6.702	< 0.75
100091	P-C	171.782	215.249	33.321	1999.807	5.096	3.835	3.75
5 100093	M-C	133.786	35.499	105.726	1499.707	5.983	5.876	2.85
100099	S-C	174.239	49.159	28.163	1392.574	4.381	5.876	3.18
100100	D-C	181.284	68.095	82.324	1613.489	5.552	5.547	3.71
100103	J-C	151.396	74.849	55.720	1666.282	4.771	3.955	3.57
100104	C-C	128.182	38.890	13.719	1100.784	3.798	5.324	2.21
10 100106	S-C	170.685	41.811	37.449	1280.697	5.096	4.464	2.52
100107	J-C	103.835	27.397	32.334	1416.481	2.910	3.740	2.57
100108	J-C	148.294	145.629	52.822	1884.389	4.484	3.597	4.62
100115	S-DLT	134.784	108.134	96.415	1597.696	8.860	5.225	3.48
100121	R-D	149.335	74.878	34.109	1759.752	5.912	3.716	4.36
15 100124	R-D	134.766	130.367	35.486	1569.219	2.489	4.877	1.21

10007188 - 062102

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Collagen IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
100125	B-D	170.790	67.078	97.770	2245.776	7.261	5.876	3.63
100126	D-D	134.313	117.116	65.560	1970.476	2.775	3.788	3.04
100129	E-D	159.707	60.388	28.962	1651.995	5.195	4.902	3.41
100131	J-D	155.166	119.774	31.852	1579.186	3.015	3.405	2.73
5 100133	M-D	146.280	24.371	75.729	2098.560	3.225	3.788	1.74
100135	H-E	167.472	41.600	66.767	1369.735	3.200	3.429	4.20
100137	D-E	158.406	25.104	68.740	1346.906	3.828	3.405	2.69
100139	S-E	139.877	38.484	42.708	1388.605	4.215	2.814	4.05
100140	L-E	158.942	30.695	181.056	1585.482	7.476	4.977	3.32
10 100141	W-E	136.761	185.300	179.774	2045.873	9.097	9.872	2.89
100142	R-E	119.383	62.037	16.170	888.744	3.286	2.673	2.07
100143	D-E	131.717	33.779	29.179	1072.170	2.978	3.144	1.91
100147	D-E	132.426	77.159	35.912	1285.138	3.515	2.696	2.12
100150	D-E	120.207	19.056	155.043	1488.729	2.298	4.196	3.66
15 100151	D-F	125.885	35.735	51.625	1243.711	3.447	3.525	3.05

10087188 - 062102

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/mL)	YKL-40 (ng/mL)	HA (ng/mL)	TIMP-1 (ng/mL)	PIIINP (ng/mL)	Coll IV (ng/mL)	$\alpha$ 2-MG (ng/mL)
100155	C-H	146.728	29.137	54.330	1242.340	4.733	3.573	3.15
100159	M-H	136.303	31.336	106.675	1567.716	4.151	3.525	3.51
100161	JF-F	155.052	1710.890	572.598	1966.460	6.226	4.634	4.27
100163	M-F	153.221	785.420	211.173	2167.501	8.269	5.698	3.57
5 100175a	B-G	148.403	55.347	130.093	1282.502	5.296	4.537	2.79
100181	M-G	137.986	69.735	31.119	1384.651	2.813	3.097	1.89
100183	D-G	168.842	181.909	58.358	1499.596	3.101	3.405	3.86
100186	M-G	184.148	2258.120	347.854	5271.196	11.670	7.756	3.82
100200	R-G	148.660	158.906	143.510	1939.499	5.530	6.080	4.02
10 100208	R-G	156.210	94.021	36.624	1379.174	5.339	4.366	< 0.75
100221	ND-H	117.196	38.393	88.913	1375.112	2.610	5.274	3.76
100222	J-H	106.131	34.544	31.603	1054.973	2.580	2.955	3.40
100229a	S-H	125.123	53.139	73.989	1567.731	3.828	3.788	2.78
100237	J-J	140.718	397.625	578.952	1824.407	11.836	6.675	2.46
15 100268	C-L	155.864	76.151	86.977	2060.140	10.963	6.310	3.66

100087198 • 0621102

TABLE 4

RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
100270	T-L	176.060	24.738	38.749	1579.990	2.549	2.625	3.35
100278	S-L	153.789	48.840	52.524	1367.051	3.039	3.167	1.91
100279	L-LG	163.352	139.202	34.492	1223.652	2.921	4.099	2.37
100287	R-L	164.414	636.110	232.253	3285.078	14.450	9.448	4.19
5 100291	MS-L	152.863	197.500	42.925	2144.445	2.180	1.993	3.78
100293	D-L	147.479	104.509	27.209	1559.538	3.151	2.696	1.83
100302	A-M	201.715	1021.070	159.330	3317.515	12.498	8.383	4.83
100306	JT-M	125.203	115.289	83.960	1722.987	3.842	2.413	4.16
100307	D-M	128.095	23.612	10.882	1378.118	3.339	3.382	1.42
10 100313a	J-MCA	164.201	192.417	76.599	1966.287	3.828	5.523	3.89
100315	K-MF	153.427	113.449	112.430	2118.580	3.515	4.561	4.46
100317	M-MCM	165.245	94.693	144.632	1611.495	4.588	3.215	1.84
100320	D-M	159.724	782.150	106.161	1581.345	13.215	5.597	2.33
100322	R-M	120.098	39.914	51.444	1443.789	2.809	3.644	3.21
15 100327	K-R	168.143	194.521	101.231	1738.827	4.295	3.835	4.35

100027168 - 062102

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
100336	E-S -M	165.374	135.711	36.329	1556.071	2.932	2.508	3.49
100347	E-M	173.070	69.889	16.945	1710.951	4.808	3.859	3.29
100348	A-M	207.186	75.06	301.583	1334.475	3.299	4.585	3.36
100350	J-M	154.867	4.418	22.250	1388.371	4.087	3.238	1.60
5 100358	A-M	140.022	15.549	88.786	1247.147	4.502	4.682	1.57
100365	A-M	96.324	26.329	43.344	1170.887	5.381	2.040	2.12
100367	B-M	161.274	30.273	46.174	1469.088	3.089	2.790	3.59
100388	A-P	230.782	275.681	938.015	4245.175	20.496	9.669	5.98
100395	D-R	125.908	24.226	15.309	1299.599	2.478	4.415	3.52
10 100397a	R-R	179.186	56.479	100.853	1455.947	9.769	4.172	3.62
100398	C-R	151.391	29.397	11.833	1100.156	2.652	2.932	2.28
100403	L-P	179.146	321.607	350.713	2061.218	8.938	4.977	2.28
100404	ML-P	179.163	1060.240	141.902	2248.495	5.959	4.123	4.07
100414	S-S	184.451	70.941	40.126	1048.761	2.549	3.026	3.18
15 100424	A-P	158.538	62.439	167.519	1320.841	5.509	5.448	3.88

10047168 , 062102

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha$ 2-MG (mg/ml)
100443	J-R	112.348	40.703	21.470	1045.054	2.663	5.647	4.34
100450	M-R	186.892	200.744	203.399	1287.107	2.586	4.040	3.77
100472	T-S	127.877	119.759	21.867	797.753	2.787	3.962	1.97
100474	J-S	118.319	55.427	19.699	939.914	2.287	3.524	3.15
5 100482	J-T	125.011	33.428	39.749	1099.832	5.159	4.386	2.56
100483	J-T	136.978	18.006	33.467	1003.831	2.586	4.485	2.06
100488	M-T	178.106	51.908	123.723	1345.743	4.464	6.053	3.84
100495	J-V	180.283	219.322	93.680	1734.925	8.387	4.064	3.82
100502	D-V	114.380	153.296	22.649	1136.728	2.774	4.757	4.47
10 100503	J-W	149.928	453.095	92.918	1422.332	5.581	4.114	3.15
100513	M-W	117.649	53.027	100.853	1335.730	5.092	4.757	3.74
100528	D-W	189.040	37.248	26.158	1103.495	4.642	3.215	2.90
100530	M-W	106.100	40.234	20.490	924.291	3.297	3.988	2.29
100534	D-A	135.702	37.357	91.410	1421.072	3.374	3.135	3.26
15 100539	M-DB	167.910	474.960	104.470	2158.238	4.957	5.566	3.99

10087168 - 062102

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Coll IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
100540	A-B	135.980	54.815	183.589	1881.935	3.068	4.485	3.34
100546	D-F	113.363	75.942	47.682	1207.833	2.633	3.055	1.57
100557	T-L	121.746	58.286	27.704	1297.060	3.842	2.060	1.83
100560	C-N	194.265	142.696	91.041	2338.303	6.332	3.267	4.25
5 100564	D-R	169.241	250.713	65.865	2407.901	5.502	4.534	3.52
100569	J-DC	160.145	64.634	43.744	1349.485	5.962	3.421	3.78
100572	K-K	162.517	260.632	209.581	1729.746	9.292	6.191	2.70
100585	K-Z	171.277	162.336	126.433	2030.404	8.907	4.881	4.71
100594	R-M	114.193	216.295	42.261	1678.540	3.545	3.602	3.26
10 100603	M-S	114.071	61.460	34.332	1693.071	4.464	4.188	2.92
100611	P-F	178.856	269.956	92.721	896.077	3.960	3.524	3.96
100614	J-MCA	204.794	245.159	322.970	3470.966	11.393	6.814	5.32
100617	C-W	159.292	51.343	38.850	1504.544	5.859	4.974	3.60
100630	E-AV	140.072	34.778	59.454	1091.420	1.969	3.161	3.34
15 100637	R-B	179.987	59.477	137.723	2077.095	4.726	3.002	4.24

10087189 - 062102

TABLE 4  
RAW DATA FROM 194 HCV PATIENTS ANALYZED FOR LEVELS OF LAMININ, YKL-40, HA, TIMP-1,  
PIIINP, COLLAGEN TYPE IV AND  $\alpha$ 2-MG

Sample ID	Patient	Laminin (ng/ml)	YKL-40 (ng/ml)	HA (ng/ml)	TIMP-1 (ng/ml)	PIIINP (ng/ml)	Collagen IV (ng/ml)	$\alpha$ 2-MG (ng/ml)
101013	T-H	177.189	507.415	237.499	1556.336	9.381	8.910	3.36
101118	G-S	163.553	282.057	150.713	2348.845	6.231	7.161	3.93
101137	S-S	175.291	2049.970	715.601	3318.137	11.450	8.458	3.71
101257	M-F	155.324	40.730	49.441	1082.686	3.240	3.421	1.71
5 101275	J-C	121.598	143.292	45.227	1454.509	1.915	2.788	4.69
101284	R-F	123.312	16.825	42.048	1072.891	4.386	3.679	1.64
101321	H-P	180.159	180.091	367.681	2235.931	11.450	10.045	3.37
101322	A-P	133.640	269.262	244.520	2015.508	9.119	8.458	4.01
101335	L-S	164.947	64.238	208.719	2545.531	11.018	6.930	3.90
10 101336	L-S	156.847	741.300	72.664	1387.281	6.097	4.683	3.73
101351	P-F	128.701	12.461	99.050	1557.519	5.486	4.161	4.01
101441	R-H	104.993	176.909	354.512	1338.637	5.271	6.445	3.50
101478	G-S	142.642	63.543	149.266	1296.972	3.859	4.411	2.64
101565	D-A	132.117	74.355	187.250	1206.734	5.962	4.782	2.92

Clinical performance parameters were analyzed for the combinations of markers best able to differentiate the presence of significant fibrosis (F2-F4) from no/mild fibrosis (F0-F1) using various 5 statistical algorithms. The statistical algorithms analyzed included univariate analysis, receiver operating characteristic curves (ROC), logistic regression, discriminant function analysis, and factorial design optimization.

10 The results of ROC analysis are shown in Table 5. The area under the curve (AUC) values represent relative diagnostic value of a single marker at the indicated cut-off. As can be seen by the decreasing AUC values, HA was shown to have the best diagnostic value 15 when used alone at the indicated cut-off, followed by PIIINP, TIMP-1,  $\alpha$ 2-MG and collagen type IV.

TABLE 5				
ROC ANALYSIS				
	AUC	Sensitivity	Specificity	Cutoff
20	HA	0.821	90.0%	62.0% 35.5 ng/ml
	PIIINP	0.777	90.8%	39.2% 3.0 ng/ml
	TIMP-1	0.773	90.8%	43.0% 1190.6 ng/ml
	$\alpha$ 2-macroglobulin	0.722	90.5%	34.6% 2.4 mg/ml
25	Collagen Type IV-7S	0.726	90.8%	24.1% 2.79 ng/ml
	YKL-40	0.696	90.8%	19.0% 34.5 ng/ml
	Laminin	0.524	90.7%	16.5% 125.2 ng/ml

Clinical performance parameters for various combinations of fibrosis markers are shown in Table 6. The best subsets, including single markers as well as 30 combinations of two to four markers and algorithms for discriminating F0-F1 from F2-F4, were generated by logistic regression. The markers included PIIINP,  $\alpha$ 2-MG, laminin and type IV collagen. As shown in Table 5, the

diagnostic performance parameters (sensitivity, specificity, PPV and NPV) were similar for the two, three and four-marker combinations identified by logistic regression in the study population, which had a fibrosis prevalence of about 60% (see lines 2-4 and 6-9).

As shown in Table 6, line 5, stepwise discriminant function analysis (SAS) resulted in identification of the 3-marker subset (PIIINP,  $\alpha$ 2-MG and laminin). The clinical performance of this combination was similar to the marker combinations identified using logistic regression.

Design of experiments software (DOE KISS, Build 8, Air Academy Associates) was used to simultaneously optimize the cutoffs of multiple variables to obtain the best performance of the panel of tests in predicting fibrosis. Using DOE KISS, a computer-aided central composite design for a combination of markers was generated; this design matrix consisted of a series of combinations of cutoffs for each of the markers in the combination. The results from these experiments (sensitivity, specificity and accuracy) for differentiating F0-F1 from F2-F4 fibrosis were recorded in the design sheet in DOE. Regression analysis was performed for each of the parameters (sensitivity, specificity and accuracy) to give cutoff values for each of the variables in the combination to achieve maximum performance for that parameter.

The five markers with best diagnostic performance in an ROC analysis (highest AUC) were HA, PIIINP, TIMP-1,  $\alpha$ 2-MG and type IV collagen (see Table 5).

Cutoffs for each of the markers in this 5-marker panel were optimized for maximum accuracy. The results shown in Table 6, line 10, indicate that, at the optimum accuracy (69.6%), the specificity was too low to be 5 useful (32.9%) while the sensitivity was high (94.8%). Similar results were obtained when the markers were optimized for sensitivity or specificity. Regression analysis showed that TIMP-1 did not have a significant effect on the accuracy, sensitivity or specificity of 10 this 5-marker panel.

TABLE 6  
PERFORMANCE OF VARIOUS MARKER PANELS

Markers	Method/ Model	Sens.	Spec.	Prevalence 59.3%			Prevalence 20%		
				PPV	NPV	Acc.	PPV	NPV	Acc.
1	HA	N/A cut-off 60ug/ml	64.96%	82.05%	84.44%	60.95%	47.50%	90.35%	78.63%
2	PIIINP	Logistic Best subset of 1	74.36%	58.97%	73.11%	60.53%	31.18%	90.20%	62.05%
3	PIIINP, AMG	Logistic Best subset of 2	80.53%	63.63%	76.86%	68.54%	31.18%	90.20%	62.05%
4	PIIINP, AMG, Laminin	Logistic Best subset of 3	78.76%	67.53%	78.44%	67.95%	37.75%	92.71%	69.78%
5	PIIINP, AMG, Laminin	Discriminant Stepwise selection	78.76%	67.53%	78.44%	67.95%	37.75%	92.71%	69.78%
6	PIIINP, AMG, Laminin, Coll IV	Logistic Best subset of 4	78.76%	64.94%	77.11%	67.09%	35.96%	92.44%	67.70%
7	PIIINP, AMG, Laminin, YKL-40	Logistic Second best subset of 4	77.87%	67.53%	78.25%	67.05%	37.48%	92.43%	69.60%
8	PIIINP, AMG, Coll IV, YKL-40	Logistic Third best subset of 4	78.76%	70.13%	79.82%	68.76%	39.73%	92.96%	71.86%
9	PIIINP, AMG, TIMP-1	Logistic "Forced" selection	78.76%	64.94%	77.11%	67.09%	35.96%	92.44%	67.70%

TABLE 6  
PERFORMANCE OF VARIOUS MARKER PANELS

Markers	Method/ Model	Sens.	Spec.	Prevalence 59.3%			Prevalence 20%		
				PPV	NPV	Acc.	PPV	NPV	Acc.
10 HA, PIIINP, AMG Coll IV, TIMP-1	DOE (for acc.) N/A	94.78%	32.91%	67.28%	81.25%	69.59%	26.1%	96.2%	45.3%
11 HA, PIIINP, Coll IV, AMG	DOE (for acc.) N/A	79.13%	75.95%	82.73%	71.43%	77.84%	45.13%	93.57%	76.59%
12 HA, CollIV, AMG	DOE (for acc.) N/A	95.65%	29.11%	66.27%	82.14%	68.56%	25.2%	96.4%	42.4%
13 HA, PIIINP, AMG	DOE (for acc.) N/A	78.26%	75.95%	82.57%	70.59%	77.32%	44.86%	93.32%	76.41%
14 HA, AMG (B)	DOE (for acc.) N/A	84.35%	73.42%	82.20%	76.32%	79.90%	44.24%	94.94%	75.60%
15 TIMP-1, HA, AMG	DOE (for acc.)	83.48%	75.95%	83.48%	75.95%	80.41%	46.46%	94.84%	77.46%
16 YKL-40, HA, AMG	DOE (for acc.)	82.61%	75.95%	83.33%	75.00%	79.90%	46.20%	94.59%	77.28%
17 HA, reflex with Coll IV and AMG	Logistic Positives reflexed	86.32%	70.51%	81.45%	77.46%	80.00%	42.26%	95.38%	73.68%

A similar 4-marker panel was analyzed by DOE as shown in Table 6, line 11. With TIMP-1 excluded, the four-marker panel was optimized for accuracy (77.8%) to give a sensitivity and specificity of 79.1% and 79.5%,  
5 respectively. These results demonstrate that the four-marker panel of HA, PIIINP,  $\alpha$ 2-MG and collagen IV has more value in differentiating F0-F1 fibrosis from F2-F4 fibrosis than a five-marker panel made up of HA, PIIINP,  $\alpha$ 2-MG, collagen IV and TIMP-1.

10 Several three-marker subsets of the four-marker panel also were analyzed by DOE. Line 12 shows the results obtained for the combination of HA, collagen and  $\alpha$ 2-MG with the results optimized for accuracy. This three-marker panel gave a very low specificity of less  
15 than 30%. In contrast, when a three-marker panel made up of HA, PIIINP and  $\alpha$ 2-MG was optimized for accuracy, performance was similar to the four-marker panel (compare lines 13 and 11 of Table 6).

A similar analysis of the two-marker panel of  
20 HA and  $\alpha$ 2-MG gave the results shown in line 14 of Table 6. This combination gave an improvement in specificity over the three-marker panel of HA, PIIINP and  $\alpha$ 2-MG (84.4% compared to 78.3%).

TIMP-1, which was observed to be a good  
25 discriminator of fibrosis in the univariate analysis, was added to the two-marker panel. As shown in line 15, the performance of the HA,  $\alpha$ 2-MG and TIMP-1 three-marker panel was similar to that obtained with the two-marker panel, and the sensitivity was improved as compared to  
30 the three-marker HA/PIIINP/ $\alpha$ 2-MG panel (83.5% sensitivity

compared to 78.3%). Furthermore, in preliminary regression analysis, TIMP-1 contributed significantly to discrimination of fibrosis in a study population with a high prevalence of severe fibrosis.

5 Another three-marker panel, made up of HA,  $\alpha$ 2-MG and YKL-40, also was optimized for accuracy in differentiating F0-F1 from F2-F4 fibrosis. As shown in Table 6, line 16, this three marker panel had a performance similar to the  $\alpha$ 2-MG/HA/TIMP-1 panel.

10 In sum, these results indicate that a  $\alpha$ 2-MG/HA/TIMP-1 or  $\alpha$ 2-MG/HA/YKL-40 panel can be useful in differentiating F0-F1 from F2-F4 fibrosis.

#### EXAMPLE II

##### DUAL OPTIMIZATION STRATEGY FOR ANALYSIS OF THE 15 $\alpha$ 2-MG/HA/TIMP THREE-MARKER PANEL

This example describes the use of multiple cut-offs for  $\alpha$ 2-MG-, HA- and TIMP-1 to achieve a relatively high degree of accuracy in a subset of a total patient population assayed.

20 Using the three-marker panel  $\alpha$ 2-MG/HA/TIMP-1 with cutoffs for  $\alpha$ 2-MG, HA- and TIMP-1 set at 35 ng/ml, 2 mg/ml, and 1000 ng/ml, respectively, samples were determined to be positive when all three of their variables were above the cut-off values, and were 25 therefore negative when one or more of the  $\alpha$ 2-MG, HA or TIMP-1 levels were below the assigned cut-off value. As shown in Table 7, in the 194 patient population, there were a total of 72 negative results, 15 of which were

false negatives, giving a negative predictive value (NPV) of 79% at the study prevalence of about 60% fibrosis%. At a prevalence of 30%, which is typical of the prevalence in a hepatology clinic, the negative 5 predictive value is over 92%, which is useful in ruling out the presence of F2-F4 fibrosis (likelihood ratio 0.22).

TABLE 7

PERFORMANCE OF  $\alpha_2$ -MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS  
WITH VARIOUS DISEASE PREVALENCES

Optimization for sensitivity to rule out fibrosis  
 $X_1 = 2.0 \text{ mg/ml}$  for  $\alpha_2$ -MG       $Y_1 = 35.00 \text{ ng/ml}$  for HA       $Z_1 = 1000.00 \text{ ng/ml}$  for TIMP-1

		Prevalence 0.598		Prevalence 0.300		Prevalence 0.200		Prevalence 0.100	
		Fib +	Fib -	Fib +	Fib -	Fib +	Fib -	Fib +	Fib -
Test +	101	21	122	261	188	450	174	215	390
Test -	15	57	72	39	512	550	26	585	610
Sens.	87.07%	LR +	3.23	87.07%	LR +	3.23	87.07%	LR +	3.23
Spec.	73.08%	LR -	0.18	73.08%	LR -	0.18	73.08%	LR -	0.18
PPV	82.79%			58.09%		44.71%		26.43%	
NPV	79.17%			92.95%		95.76%		98.07%	
Accuracy	81.44%			77.27%		75.88%		74.48%	

**TABLE 7**

PERFORMANCE OF  $\alpha$ 2-MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS  
WITH VARIOUS DISEASE PREVALENCES

		Optimization for specificity to rule in fibrosis				Z1 = 1575.00 for TIMP-1 ng/ml					
		X1 = 2.0 mg/ml for $\alpha$ 2-MG Y1 = 60.00 ng/ml for HA		Prevalence 0.598		Prevalence 0.300		Prevalence 0.200		Prevalence 0.100	
		Fib +	Fib -	Fib +	Fib -	Fib +	Fib -	Fib +	Fib -	Fib +	Fib -
Test +	53	1	54	137	9	146	91	10	102	46	12
Test -	48	20	68	124	179	304	83	205	288	41	231
Equiv											272
	101	21	122	261	188	450	174	215	390	87	242
Sens.	52.48%	LR +	11.02	52.48%	LR +	11.02	52.48%	LR +	11.02	52.48%	LR +
Spec.	95.24%	LR -	0.50	95.24%	LR -	0.50	95.24%	LR -	0.50	95.24%	LR -
PPV	98.15%			93.86%			89.91%			79.84%	
NPV	29.41%			59.11%			71.25%			84.80%	
Accuracy	59.84%			70.40%			76.12%			83.93%	

TABLE 7									
PERFORMANCE OF $\alpha$ 2-MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS WITH VARIOUS DISEASE PREVALENCES									
Final performance after dual optimization									
	Prevalence	0.598	Prevalence	0.300	Prevalence	0.200	Prevalence	0.100	
	Fib +	Fib -	Fib +	Fib -	Fib +	Fib -	Fib +	Fib -	
5	53	1	54	137	9	146	91	10	102
	63	77	140	163	691	854	109	790	898
	116	78	194	300	700	1000	200	800	1000
Sens.	45.69%	LR +	35.64	45.69%	LR +	35.64	45.69%	LR +	35.64
Spec.	98.72%	LR -	0.55	98.72%	LR -	0.55	98.72%	LR -	0.55
PPV	98.15%			93.86%			89.91%		79.84%
NPV	55.00%			80.92%			87.91%		94.24%
Accuracy	67.01%			82.81%			88.11%		93.42%

TABLE 7

PERFORMANCE OF  $\alpha_2$ -MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS WITH VARIOUS DISEASE PREVALENCES

Population Prevalence	194 0.598	Population		Population Prevalence	0.300	Population		Population Prevalence	0.200	Population		Population Prevalence	0.100
		1000 False	39 Negative			1000 False	26 Negative			1000 False	13 Negative		
False Negative	15	False Negative	39	False Negative	39	False Negative	26	False Negative	26	False Negative	13	False Negative	13
False Positive	1	False Positive	9	False Positive	9	False Positive	10	False Positive	10	False Positive	12	False Positive	12
Total Incorrect	16	Total Incorrect	48	Total Incorrect	48	Total Incorrect	36	Total Incorrect	36	Total Incorrect	24	Total Incorrect	24
Total Correct	110	87.3%	649	Total Correct	649	Total Correct	676	Total Correct	676	Total Correct	703	Total Correct	703
2nd Round													
15 Equivocal	48	Fib + Equivocal	124	Fib + Equivocal	124	Fib + Equivocal	83	Fib + Equivocal	83	Fib + Equivocal	41	Fib + Equivocal	41
Fib - Equivocal	20	Fib - Equivocal	179	Fib - Equivocal	179	Fib - Equivocal	205	Fib - Equivocal	205	Fib - Equivocal	231	Fib - Equivocal	231
Total Equivocal	68	35.1%	304	Total Equivocal	304	Total Equivocal	288	Total Equivocal	288	Total Equivocal	272	Total Equivocal	272

TABLE 7

PERFORMANCE OF  $\alpha$ 2-MG/HA/TIMP-1 PANEL WITH DUAL OPTIMIZATION STRATEGY IN POPULATIONS  
WITH VARIOUS DISEASE PREVALENCES

Final performance after dual optimization without equivocals									
		Prevalence	0.598	Prevalence	0.300	Prevalence	0.200	Prevalence	0.100
	Fib +	Fib -	Fib +						
5	Test +	53	1	54	137	9	146	10	102
	Test -	15	57	72	39	512	550	26	585
	Sens.	77.94%	LR +	45.21	77.94%	LR +	45.21	77.94%	LR +
	Spec.	98.28%	LR -	0.22	98.28%	LR -	0.22	98.28%	LR -
	PPV	98.15%			93.86%			89.91%	
10	NPV	79.17%			92.95%			95.76%	
	Accuracy	87.30%			93.14%			94.93%	
	% of pop.		% of pop.		% of pop.		% of pop.		% of pop.
	Test +	27.8%		Test +	14.6%		Test +	10.2%	Test +
	Test -	37.1%		Test -	55.0%		Test -	61.0%	Test -
15	Equiv	35.1%		Equiv	30.4%		Equiv	28.8%	Equiv

Furthermore, of the 122 test positives using the 35 ng/ml, 2 mg/ml, and 1000 ng/ml cut-offs, 21 of the test positives were false, giving a positive predictive value (PPV) of 82.8%. However, at a more typical 5 prevalence of 30% fibrosis, the positive predictive value falls to about 58% (see Table 7). Thus, in a population with a typical prevalence, a positive result would not have sufficient predictive value to be useful as a diagnostic.

10 In order to increase the positive predictive value for at least a subset of the total patient population, samples positive by the primary analysis were further evaluated for positivity for the three markers using a second set of cut-off values which were higher 15 than the first set. By evaluating those samples positive after a primary analysis at higher cutoffs, the severe fibrosis samples within this group can be determined to be positive with a relatively high predictive value. Those samples that test negative by the secondary 20 evaluation are considered "indeterminate" in that their fibrosis status cannot be determined with good predictive value.

Table 7 shows performance of the α2-MG/HA/TIMP-1 panel assay with the dual optimization 25 strategy. The primary cut-offs were set at 2.0 mg/ml, 35 ng/ml and 1000 ng/ml to achieve a relatively high sensitivity in the primary analysis. Any samples having all three of their α2-MG, HA and TIMP-1 levels above the assigned cut-off values were indicated to be positive. 30 The 122 test positives obtained by the primary analysis were re-evaluated using 2.0 mg/ml, 60 ng/ml and 1575

ng/ml as the  $\alpha$ 2-MG, HA and TIMP-1 cut-offs and the criteria that the samples must have  $\alpha$ 2-MG, HA and TIMP-1 values above the assigned cut-off values to be positive.

Using the second set of cut-off values, 54 of 5 the 122 patients were determined to be positive, only 1 of which was a false positive. The positive predictive value was 98.2% at 59.8% fibrosis prevalence, and was 93.9% at the more typical 30% fibrosis prevalence. In sum, of the 194 patients, 72 were classified as negative 10 and 54 were classified as positive, while 68 samples had indeterminate results and could not be definitively classified. Furthermore, when the indeterminate samples are excluded, the three-marker assay has a positive predictive value of more than 93% and a negative 15 predictive value of close to 93% in a typical population having a 30% fibrosis prevalence.

Table 8 shows a comparison of the performance of the  $\alpha$ 2-MG/HA/TIMP-1 three-marker panel with the six marker panel described in Poynard et al., Lancet 357:1069 20 (2001).

TABLE 8  
COMPARISON OF PERFORMANCE OF  $\alpha_2$ -MG/HA/TIMP-1 PANEL WITH 6 MARKER PANEL OF POYNARD ET AL.  
Prometheus

		Biopsy				Biopsy	
		Fib +	Fib -	Fib +	( $> 0.08$ )	Fib +	Fib -
Test +	53	1	54	Test + ( $> 0.08$ )	45	5	50
Test -	15	57	72	Test - ( $< 0.20$ )	13	106	119
Equivocal	48	20	68	Equivocal	80	90	170
Total Pop.	116	78	194	Total Pop.	138	201	339
<b>Prevalence</b>	<b>.05979</b>			<b>Prevalence</b>	<b>0.4071</b>		
Sensitivity	0.7794			Sensitivity	0.7759		
Specificity	0.9828			Specificity	0.9550		
PPV	0.9815			PPV	0.900		
NPV	0.7917			NPV	0.8908		
Accuracy	0.8730			Accuracy	0.8935		
% Equivoc	.03505	68/194		% Equivoc	.0515	170/339	
False Pos	1	of 54 test +	1.85%	False Pos	5	of 50 test +	10.00%
False Neg	15	of 72 test -	20.83%	False Neg	13	of 119 test -	10.92%
Prometheus							
		Biopsy				Biopsy	
		Fib +	Fib -	Fib +	( $> 0.08$ )	Fib +	Fib -
Test +	186	8	194	Test + ( $> 0.08$ )	133	15	147
Test -	53	433	486	Test - ( $< 0.20$ )	38	313	351
Equivocal	168	152	320	Equivocal	236	265	501
Total	407	593	1000	Total Pop.	407	593	1000
<b>Prevalence</b>	<b>0.4071</b>			<b>Prevalence</b>	<b>0.4071</b>		
Sensitivity	0.7794			Sensitivity	0.7759		
Specificity	0.9828			Specificity	0.9550		
PPV	.9607			PPV	0.900		
NPV	.8917			NPV	0.8907		
Accuracy	.09113			Accuracy	0.8935		
% test pos	19.4%			% test pos	14.7%		
% test neg	48.6%			% test neg	35.1%		
% Equivoc	32.0			% Equivoc	50.1%		
False Pos	8	of 194 test +	3.93%	False Pos	15	of 147 test +	10.00%
False Neg	53	of 486 test -	10.83%	False Neg	38	of 35 test -	10.93%

\* alpha2-macroglobulin, alpha2-globulin, total bilirubin, gama-globulin, apo A1 and GGT

These results indicate that the  $\alpha$ 2-MG/HA/TIMP-1 three-marker panel can be useful in differentiating F0-F1 fibrosis from F2-F4 fibrosis with very good accuracy. These results further indicate that a combination 5 fibrosis marker assay can be useful in determining the fibrosis status of a portion of the patients tested with very good accuracy, while the remaining patients are candidates for biopsy.

### EXAMPLE III

10 **ASSAYS FOR  $\alpha$ 2-MACROGLOBULIN, HYALURONIC ACID AND TISSUE INHIBITOR OF METALLOPROTEINASES-1**

A. Quantitation of human  $\alpha$ 2-macroglobulin ( $\alpha$ 2-MG)

Serum levels of human  $\alpha$ 2-macroglobulin were quantitated using the Beckman Array® 360 System as 15 follows to determine  $\alpha$ 2-MG levels in the range of 0.75-270 mg/ml.

The Beckman Array® 360 system was used for determination of  $\alpha$ 2-MG concentrations. This system utilizes a nephelometer which measures the rate of 20 light-scatter formation resulting from an immunoprecipitation reaction between  $\alpha$ 2-MG antigen in a sample with antibody to human  $\alpha$ 2-MG. After passing a beam of light through the solution in a flow cell, the intensity of light scattered by the formed macromolecular 25 particles of insoluble complexes suspended in solution is detected and measured by the nephelometer. The increase in light scatter resulting from the antigen-antibody reaction is converted to a peak rate signal proportional to the  $\alpha$ 2-MG concentration in the sample. The resulting

formation of complexes and the consequent change in the intensity of scattered light occurs at a rate that increases gradually at first, then rapidly, and finally proceeds through a peak rate of change for the component 5 being analyzed.

Serum samples were drawn from fasting individuals and generally physically separated from cells within 2 hours from the time of collection as set forth in NCCLS publication H 18-A. Samples not assayed within 10 72 hours were stored frozen at -15°C to -20°C. Frozen samples were at most thawed one time. Grossly hemolyzed, highly lipemic or turbid specimens were rejected for further analysis.

Reagents were removed from storage at 4°C and 15 used immediately. Buffers and Diluents were mixed thoroughly by inversion prior to being added to the instrument. Set-up, priming and calibration were performed according to the manufacturer's instructions with samples diluted 1:36. Relatively concentrated 20 samples such as undiluted samples or 1:6 dilutions were generally avoided. Grossly lipemic sample were diluted 1:2 with diluent before assaying. Dust particles or other particulate matter, which can result in extraneous light-scattering signals, in the reaction solution were 25 avoided. Prior to assaying samples, any air bubbles or foam in the sample cups and reagent bottle were removed by using a disposable transfer pipette or pipette tip to aspirate the bubbles. DTT was avoided in the work area.

Samples were analyzed for α2-MG concentration 30 as follows. The Reagent Wheel (left wheel) on the

instrument was loaded with AMG antiserum in space #2. Dilution segments were loaded with 150  $\mu$ L of control or sample in the wells on the larger side of the fan shaped segments. Segments and initial dilution control/sample 5 cups were marked for identification. Bubbles were avoided while controls and serum samples were loaded.

Vigil™ Protein Control Levels 1 and 3 (3 drops) was placed in cups 1 and 3, respectively. Biorad Liquichek™ Immunology Control Level 2 (150  $\mu$ L) was placed 10 in cup 2. Patient samples (150  $\mu$ L) were added to sequential cups. Segments were placed on right wheel beginning at position #1. Evaporation covers were placed over Reagent and Sample Wheels.

On the Master Screen menu, the RESULTS RECALL 15 (F3) was selected before (F4) CLR CUR RUN. After returning to the MASTER SCREEN, the SAMPLE PROGRAM (F1) was selected. ENTER was selected when Reagent wheel #1 appeared and at each cup number. The control ID or sample Acc.# was entered. Test "2" was selected, and 20 SAVE CUP (F1) was selected for each cup. START was selected to begin the analysis. At the end of the run, (Y) was selected in response to CLEAR CURRENT RUN & START NEXT RUN.

Results were reported by the Beckman Array® 360 25 in mg/dl using whole numbers in the Pros System. Samples were diluted routinely by the instrument 1:36. Samples greater than 750 mg/dl were assayed at a 1:216 dilution by the instrument. Samples having a concentration less than 75 mg/dl at a 1:36 dilution are reported as 30 <75 mg/dl. At initial dilutions the Beckman analytical

range was 75-750 mg/dl, while the extended range was 75-27,000 mg/dl. The range for normal individuals as verified at Prometheus Laboratories was 103-274 mg/dl.

Quality control was performed as follows.

- 5 Three levels of controls were used: low, medium and high. Controls were within 2 standard deviations, except that runs were accepted with two controls within 2 standard deviations and the third control between 2 and 3 standard deviations. The controls used were Beckman Vigil I and  
10 III and Biorad Level II. Controls were assayed with each sample run.

The assay is calibrated every 14 days, and also when changes in reagent lots occur or when a major change has occurred in the instrument. Linearity is confirmed  
15 every 6 months with appropriate linearity material. This is done to ensure consistent performance over time and to comply with State and National standards.

- Assay calibration verification is performed every 6 months to ensure consistency over time. A  
20 minimum of five verification samples including minimum, mid-point, and maximum concentrations are evaluated every 6 months. The coefficient of variation (%CV) of the verification sample results must be less than 15% in order to report out patient sample results.

25 B. Quantitation of Hyaluronic Acid (HA)

Serum levels of HA were determined using the Hyaluronic Acid (HA) Quantitative test kit (Catalog #029001) from Corgenix essentially as follows.

Serum samples were stored at -70°C. Multiple freeze/thaw cycles were avoided, with a maximum of 4 freeze/thaw cycles per sample. The kits were stored at 2-8°C.

5 Prior to use, the kit and patient samples were equilibrated to room temperature (18-28°C). The pouch of coated strips also was equilibrated to room temperature before opening. Wash solution (0.01 M PBS, pH 7.35 +/- 0.1) was prepared by diluting the 33X PBS wash  
10 concentrate with distilled water and adjusting the pH of the final solution to pH 7.35 +/- 0.1.

All blanks, standards, controls and samples were assayed in duplicate. A water blank for calibration of the spectrophotometer was included with each plate and  
15 remained empty until addition of 200 µl water immediately prior to reading. Reaction buffer without serum sample was used for the reagent blank, which represented the 0 ng/ml HA reference solution, and was treated the same as patient samples and reference solutions in subsequent  
20 assay steps. Three known patient samples (low, middle and high) were run with each assay. In addition, 50 ng/ml HA, 100 ng/ml HA, 200 ng/ml HA, 500 ng/ml HA and 800 ng/ml HA reference solutions supplied with each kit were assayed as described further below.

25 HA reference solutions and patient samples were diluted 1:11 by addition of 25 µl reference solution or sample to 250 µl of reaction buffer and mixed by gentle vortexing. The diluted reference, samples and controls were added (100 µl) to each well. The water blank  
30 remained empty. The plate was covered and incubated for

60 minutes at room temperature. After the incubation was complete, the contents of the wells were removed by aspiration. Plates were washed four times with 1X wash solution while avoiding the plates drying out between 5 washes. The plate was blotted vigorously on paper towels to remove residual buffer after the last wash.

HRP-conjugated HA binding protein solution (100  $\mu$ l) was added to all wells except the water blank before covering the plate and incubating for 30 minutes 10 at room temperature. After the incubation was complete, the plate was washed four times as described above. Substrate solution (100  $\mu$ l 3,3',5,5'-tetramethylbenzidine and hydrogen peroxide, stabilized) was then added to each well except for the water blank. The covered plate was 15 then incubated for 30 minutes at room temperature. The plate was kept in the dark.

The OD<sub>650</sub> of the 800 ng/ml HA standard was determined. For an OD less than 0.500, the substrate incubation was continued and the OD monitored to 20 determine if the OD had reached 0.500. For an OD greater than 0.500 or after one hour of substrate incubation even if the OD had not reached 0.500, the reactions were terminated by addition of 100  $\mu$ l of Stopping Solution (0.36 N sulfuric acid) to each well except the water 25 blank. The stop solution was added in the same order and at approximately the same rate as addition of the substrate solution. Before reading the optical densities, 200  $\mu$ l distilled water was added to the water blank. The OD of each well was read at 450 nm (650 nm 30 reference) within one hour after "zeroing" the plate reader against the water blank.

The following criteria were used to determine if the assay was reliable. The mean OD value of the reagent blank (zero standard) was less than 0.10. Readings greater than 0.10 were considered indicative of 5 possible substrate or reagent contamination, and results were not reported under these conditions. The mean OD value of the 500 ng/ml HA reference was 0.800 or greater. Controls for the three known patient samples were within the following ranges: Low control: 78.6 to 117.2 ng/ml. 10 Mid control: 148.5 to 214.1 ng/ml. High control: 297.8 to 460.7 ng/ml. Samples with HA concentrations greater than 800 ng/ml were further diluted and assayed a second time to obtain a more accurate result.

The known patient controls and samples were 15 determined from a standard 4-parameter curve generated using Softmax and reported in ng/ml. The patient values were not reported if the concentration exceeded the concentration of the highest standard. For patient values greater than the concentration of the highest 20 standard at a 1:11 dilution, samples were assayed at a 1:55 dilution and, if necessary, at higher dilution.

The HA ELISA assay is evaluated every six months to ensure consistent performance over time. A minimum of five samples with previously known HA values 25 are evaluated in a blinded fashion to the operator. For the assay performance to be acceptable, results for negative samples must be negative, and results for positive samples must be positive and yield results within 15% of the previously obtained values. If greater 30 than 20% of the validation samples fail the performance criteria, troubleshooting is implemented, and the assay

is not used to report patient data until acceptable assay performance are reestablished.

C. Quantitation of Tissue Inhibitor of Metalloproteinases-1 (TIMP-1)

5            Serum levels of TIMP-1 were determined using the Biotrak™ test kit (Catalog# RPN2611) from Amersham Pharmacia Biotech (Piscataway, NJ) essentially as follows.

Kit contents were thawed and equilibrated  
10 to 20-25°C. Serum samples were stored frozen at -70°C.  
Repeated freeze-thaw cycles of the samples were  
minimized, with a maximum of six freeze-thaw cycles.

Assay reagents were prepared as follows and  
stored at 2-8°C for at most 7 days. Assay buffer 1 (0.1  
15 M phosphate buffer, pH 7.5, with 0.9% (w/v) sodium  
chloride, 0.1% (w/v) BSA and 0.1% Tween-20) was prepared  
by adding distilled water to the assay buffer concentrate  
and adjusting the final volume to 100 ml.

Anti-TIMP-1 horseradish peroxidase conjugate  
20 was prepared in assay buffer 1 essentially as follows.  
To the stock bottle containing lyophilized conjugate, 11  
ml diluted assay buffer 1 was added; the contents were  
mixed gently until completely dissolved while avoiding  
vigorous agitation and foaming. Wash buffer (0.1 M  
25 phosphate buffer, pH 7.5, containing 0.05% Tween-20) was  
prepared by adding distilled water to the wash buffer  
concentrate and bringing the final volume to 500 ml,  
followed by thorough mixing.

The 100 ng/ml TIMP-1 stock solution was prepared as follows and stored at 2-8°C. The lyophilized TIMP-1 standard was reconstituted in 0.1 M phosphate buffer, pH 7.5, containing 0.9% (w/v) sodium chloride, 5 0.1% (w/v) bovine serum albumin and 0.1% Tween-20 to make a standard TIMP-1 stock solution of 100 ng/ml. The contents were mixed gently until completely dissolved without vigorous agitation or foaming. Additional standards (1.565, 3.13, 6.25, 12.5, 25 and 50 ng/ml) for 10 a standard curve were prepared fresh before each assay by two-fold serial dilution of the 100 ng/ml stock solution into assay buffer 1 in 1.2 ml dilution tubes. A zero standard (blank) was also prepared.

The pouch containing the microtiter plate was 15 opened after equilibration to room temperature. All samples and standards were assayed in duplicate, and standards for a standard curve were present on each plate. On each plate, seven standards, two controls and a maximum of different 39 samples were present in 20 duplicate.

Samples were diluted 1:120 in tubes by mixing 595 µl assay buffer 1 with 5 µl serum. The dilutions were mixed by vortexing. Using a multichannel pipettor, 100 µl of blank, standards and diluted samples 25 were added to individual wells on a microtiter plate. The plate was covered with the lid provided and incubated at room temperature for exactly two hours. Following the two hour incubation, the contents of the wells were aspirated, and each well was washed four times with wash 30 buffer, with complete filling and aspiration of the wells

after each wash. After the final wash, the plates were blotted on paper towels to remove residual wash buffer.

Peroxidase conjugate (100  $\mu$ l) was added to each well using a multichannel pipettor, and the covered plate  
5 incubated at room temperature for exactly two hours. After the incubation, the wells were aspirated and washed as before. Immediately upon conclusion of the incubation, 100  $\mu$ l of room temperature equilibrated TMB substrate (3,3',5,5'-tetramethylbenzidine/ hydrogen  
10 peroxide in 20% (v/v) dimethylformamide) was added to each well. The plates were covered and incubated for exactly 30 minutes at room temperature. In some cases, the reactions were monitored at 630 nm. The reactions were stopped by addition of 100  $\mu$ l 1 M sulfuric acid to  
15 all wells. Absorbance was determined at 450 nm within 30 minutes.

Control and patient samples values were determined using a standard curve (4-parameter curve fit) generated using Softmax. Concentration values from the  
20 standard curve were multiplied by the dilution factor (120) to obtain actual concentrations, reported in ng/ml. Quality of the assay was confirmed using known serum samples. The low control was in the range of 668.1 to 979.9 ng/ml. The high control was in the range of 2677.9  
25 to 3300.2 ng/ml. Patient values generally did not exceed the concentration in ng/ml of the highest standard. Where the patient value was greater than the concentration of the highest standard at a 1:120 dilution, the result was reported as greater than 120  
30 times the concentration of the highest standard.

The TIMP-1 ELISA assay is validated every six months to ensure consistent performance over time. A minimum of five samples with previously known values are evaluated in a blinded fashion to the operator. Results 5 for negative samples must be negative. Results for positive samples must be positive and must yield results within 15% of the previously obtained values. Where greater than 20% of the validation samples fail the performance criteria, troubleshooting is implemented. 10 Further patient data are not reported until acceptable assay performance is reestablished.

All journal article, reference and patent citations provided above, in parentheses or otherwise, whether previously stated or not, are incorporated herein 15 by reference in their entirety.

Although the invention has been described with reference to the examples provided above, it should be understood that various modifications can be made without departing from the spirit of the invention. Accordingly, 20 the invention is limited only by the claims.